


# Safety Notice

## Medical Devices

### Interference Of *N*-Acetylcysteine (NAC) with Enzyme-based Assays Used following Paracetamol Overdose



#### Priority 3 – Advisory

HPRA Safety Notice: SN2015(09)

Issue Date: 21 May 2015

| MANUFACTURER / SUPPLIER | HPRA CASE REFERENCE |
|-------------------------|---------------------|
| Several                 | N/A                 |

#### ISSUE

The Health Product Regulatory Authority (HPRA) would like to remind customers of the potential inference of *N*-Acetylcysteine (NAC) with enzyme-based assays used in patients diagnosed with paracetamol overdose, as these can result in biased results.

The HPRA has been notified of a number of field safety corrective actions (FSCAs) undertaken by manufacturers relating to the above issue.

#### ACTION OR RECOMMENDATIONS

The HPRA advises that users:

- 1 Identify these enzyme-based assays used in your laboratories.
- 2 Check your internal protocols for measuring results such as Paracetamol level, Enzyme Creatinine, Lactate, Glucose, and their respective validation methods.

- 3 If in doubt, check with the manufacturer of the device if the enzyme assay has the potential to be affected by biased readings as a result of this inference.
- 4 Consider the need to carry out a review of previously reported assay test results if deemed necessary.

| TARGET GROUPS   |   |
|---|---|
| A&E Departments<br>Hospital Managers / CEOs<br>Laboratory Managers<br>Laboratory staff<br>Chief Medical Scientists<br>Medical Scientists<br>Risk Managers<br>Clinical Directors | Nursing Managers<br>Nursing staff<br>Intensive Care Units<br>Purchasing Managers<br>Supplies Managers<br>Hospital personnel<br>Private Medical Practitioners<br>All wards |

**BACKGROUND**

Paracetamol is a common analgesic and antipyretic drug which, when ingested in excessive concentrations, may lead to acute renal failure and liver damage. Paracetamol overdose is associated with a significant depletion of glutathione. When the intermediate paracetamol metabolite, *N*-acetyl-*p*-benzoquinoneimine (NAPQI) is formed, with the help of glutathione it binds covalently to sulphhydryl groups. But when glutathione eventually is depleted it binds to proteins and enzymes within hepatocytes initiating an inflammatory response and leading to cell death.

*N*-Acetylcysteine (NAC) is the antidote for paracetamol overdose as it was shown to be effective in preventing liver damage when given within 8 hours of the overdose. NAC has the ability to serve as a potent sulphhydryl donor allowing the restoration of depleted glutathione. However, when NAC is administered in recommended therapeutic concentrations it can interfere with commonly used enzymatic-colorimetric assays.

Tests which are commercially available that adopt enzyme mechanisms for analyte detection such as:

- The use of Peroxidases as a catalyst in Trinder reactions, and
- The use of Acyl Amidohyrolase and 8-hydroxyquinoline

have been found to produce potential negative bias in the test results.

The level of interference varies depending on the analyte being tested and the assay in use. NAC has a half-life of approximately 6 hours and the significance of the interference is shown to lessen with time and varies with drug plasma level.

In light of this, the HPRa would like to remind lab users of these interferences and to address any concerns with the respective manufacturers.

## MANUFACTURER / AUTHORISED REPRESENTATIVE/ DISTRIBUTOR CONTACT INFORMATION

Enquiries to the **manufacturer** should be addressed to the contact details found on the device labelling / instructions for use.

## HPRA CONTACT INFORMATION

All **adverse incidents** relating to a medical device should be reported to:

Health Products Regulatory Authority  
Kevin O'Malley House  
Earlsfort Centre  
Earlsfort Terrace  
Dublin 2

Telephone: +353-1-6764971  
Fax: +353-1-6344033  
E-mail: [devicesafety@hpra.ie](mailto:devicesafety@hpra.ie)  
Website: [www.hpra.ie](http://www.hpra.ie)