

# Safety Notice

## Medical Devices

### T34

## Ambulatory Syringe Pumps

Priority 2 – Warning



HPRA Safety Notice: SN2020(03)

Issue Date: 12<sup>th</sup> March 2020

MANUFACTURER / SUPPLIER	HPRA CASE REFERENCE
BD/CME	V42655

ISSUE
<p>BD/CME has identified, through Post Market Surveillance feedback, that the syringe pump motor block mechanism may become affected over time by “wear / tear”. This scenario causes under infusion situations without the pump alarm being generated.</p> <p>Further information is provided in the accompanying follow up field safety notice (FSN).</p>

ACTION OR RECOMMENDATIONS
<p>The HPRA advises that users:</p> <ol style="list-style-type: none"> <li>1 Refer to the accompanying FSN and follow the instructions provided by the manufacturer and the distributor.</li> <li>2 Forward a copy of this Safety Notice and FSN to all relevant personnel within your organisation and to any other organisations/persons to which/whom these devices have been transferred.</li> </ol>

- 3 Acknowledge receipt of the FSN if you have not already done so.
- 4 Report any concerns/incidents involving these devices and/or this remediation activity to the manufacturer and the HPRA.

TARGET GROUPS	
All wards Carers Clinical Directors Community care managers Community nurses Hospices Hospital Managers / CEOs Nursing Homes	Nursing Managers Nursing staff Oncology Nurse Specialists Palliative Care Staff Purchasing Managers Risk Managers Supplies Managers

**BACKGROUND**

The manufacturer has advised that this issue can be observed when the actuator does not move as the lead screw is rotating, resulting in potentially no pressure on the syringe plunger to deliver syringe contents.

BD/CME has advised that the cause of this issue is the actuator nut being worn enough to stop or lose traction on the lead screw due to irregular surface of the lead screw. The "wear / tear" effect can be visually observed as the lead screw is progressively wearing the plastic actuator nut overtime, leaving white plastic debris along the lead screw.

BD/CME have revised the Directions for Use (DFU) to include specific instructions for users to check for plastic white debris on the lead screw during pump set up. The Technical Service Manual (TSM) instructions have also been updated to outline changes to the maintenance and repair processes to check for plastic white debris on the lead screw and to preventively change the lead screw and/or motor block assembly where necessary.

BD/CME advise that the pump has to be checked prior to each use for plastic white debris. Users responsible for pump set-up need to double check that the infusion rate and amount was as prescribed, due to the fact that the pumps will **not alarm** the user in case of under-infusion.

The HPRA is issuing this Safety Notice to raise awareness of this field safety corrective action in Ireland. Please refer to the accompanying FSN for further details.

HPRA advise users of the T34 pump to ensure they are up to date with all recommendations issued by BD/CME in relation to this device. Copies of HPRA Safety Notices 2016(05), 2016(41), 2018(09), 2018(31), 2018(34), 2019(07), 2019(29) and FSNs issued by the manufacturer are available on the HPRA website.

## MANUFACTURER / DISTRIBUTOR CONTACT INFORMATION

Enquiries to the **manufacturer** should be addressed to:

Caesarea Medical Electronics Ltd.  
(now part of BD)  
16 Shacham Street  
Industrial Park Caesarea North  
P.O.Box 3009  
Caesarea 3088900, Israel

Telephone: +44 (0) 800 090 2460  
(UK Number)  
E-mail: [BDUKFieldAction@bd.com](mailto:BDUKFieldAction@bd.com)  
Website: [www.cme-infusion.com](http://www.cme-infusion.com)

Enquiries to the **distributor** should be addressed to:

Rockford Healthcare Limited  
3 The Westway Centre  
Ballymount Avenue  
Dublin 12

Telephone: +353-1-4509050  
E-mail: [fieldsafetyaction@rockford.ie](mailto:fieldsafetyaction@rockford.ie)  
Website: [www.rockford.ie](http://www.rockford.ie)

## HPRRA CONTACT INFORMATION

All **adverse incidents** relating to a medical device should be reported to:

Health Products Regulatory Authority  
Kevin O'Malley House  
Earlsfort Centre  
Earlsfort Terrace  
Dublin 2

Telephone: +353-1-6764971  
E-mail: [devicesafety@hpra.ie](mailto:devicesafety@hpra.ie)  
Website: [www.hpra.ie](http://www.hpra.ie)