Field Safety Notice

Product: Seronorm™ Paediatric
Article no.: 208005
FSN ID: 2016-03
Action: Information

Date issued: 18/4-16

Attention: SERO AS Distributor /SERO AS Customer

Details of affected device:

Product Name: Seronorm™ Paediatric
Article no.: 208005
Lot numbers: 1412555

Description of the problem

The component phenylalanine has been present in the above mentioned product as an added amount since 2007. No analytical data has been presented for this component and the product has not been intended to be used as a trueness control for phenylalanine. SERO has not marketed the product as a control for analysis used in the diagnosis of phenylketonuria and the product has been CE-marked following the conformity assessment procedure for general IVD-products.

SERO has many CE-marked IVD- products since the introduction of the IVD-directive, including also list II B products. During a routine audit of SERO's Quality system from the Notified Body VTT (10.3.2016) concerning products CE-marked with the approval of VTT in accordance with List IIB requirements, a complete review of all SERO products was performed. This included a review of the composition, assigned values and stated intended use. During the audit, it was brought to our attention that customers may potentially use the product Seronorm™ Paediatric as a precision control for analysis used for the diagnosis of phenylketonuria and that the product therefore should have been considered as a candidate for the approval of a Notified Body as stated in list II B of the IVD-directive.

To ensure that the product is not used for the control of analysis intended for the diagnosis of phenylketonuria, we have taken the decision to delete the information about the added amount of this component in the lot-specific product documentation.

The version number of the updated lot-specific product documentation is as follows:

- Seronorm Paediatric lot # 1412555 2016-04

There will be no reference to this component in any other promotion or information material concerning the product.

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Necessary actions to be taken by the distributor / user:

Please inform your customers by distributing this letter to all users of the above mentioned product. New updated lot-specific product documentation is available in the SERO library online (www.sero.no/theserolibrary). A copy of the updated lot-specific product documentation for the concerned lots is attached with this field safety notice.

Any reportable incidents related to the use of this product as a control for analysis intended for the diagnosis of phenylketonuria should be reported to SERO.

Risk or patient impact

Negligible.

Recommended patient follow up:

To this date, SERO has not received any reports of incidents affecting patients related to the use of the product upon diagnosis of phenylketonuria.

Each laboratory must consider the need for patient follow up depending on their use of the product.

Transmission of this Field Safety notice:

This notice needs to be passed on to all those who need to be aware within your organization or to any organization where Seronorm™ Paediatric lot no. 1412555 has been transferred.

Please transfer this field safety notice to other organizations on which this action has an impact (if appropriate).

Contact information

If you have any questions regarding this information, please contact our Customer Service (+47 66 85 89 00).

We kindly ask for your understanding for this information and thank you for your cooperation.

Sincerely

Silje Solheim Johnsen
Head of Product Management

Ase Davanger
Quality Manager

Seronorm™ Paediatric

REF 208005

LOT 1412555

2019-01

6 x 3 mL

IVD

CE

Collaborating laboratories

The analytical data of Seronorm™ Paediatric have been elaborated in collaboration with the following independent laboratories:

Collaborating laboratories:

- The Central Laboratory of Vestfold Hospital, Tønsberg, Norway
- Laboratory of Akershus University Hospital, Lørenskog, Norway
- The Laboratory of SERO AS, Billingstad, Norway

EN

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

After reconstitution, all analytes are stable for 7 days when stored tightly capped at 2-8 °C.

Vials aliquoted and frozen within 30 minutes at ≤ -20 °C in screw cap tubes with O-ring are stable for 1 month.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

Assignment of values

Seronorm™ Paediatric is intended for use as an assayed quality control material/serum to monitor the precision and trueness of laboratory measurement procedures. Analytical values of components which can be used to monitor trueness are given with an uncertainty value (U) and traceability in

accordance with ISO 17511¹⁾ and ISO 18153²⁾. All the analyses have been performed by independent collaborative laboratories using SERO's protocol. The reference values are based on reference methods or transferred values in accordance to certified reference materials. At the time of establishment of this package insert, this reference method or material constitutes the top of the calibration hierarchy for the given components. These values are presented with shaded background.

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor $k = 2$), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. For information, the relevant Rilibak range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

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1) ISO 17511. In Vitro Diagnostic Medical Devices- Measurement of quantities in samples of biological origin- Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.
 2) ISO 18153 In vitro diagnostic medical devices -- Measurement of quantities in biological samples -- Metrological traceability of values for catalytic concentration of enzymes assigned calibrators and control materials.
 3) This is not a Rilibak range, but SERO's assessment