



January 03, 2017

URGENT FIELD SAFETY NOTICE – FSN 29740-B

PRODUCT	REF	SOFTWARE VERSION
UniCel DxH 800 Coulter Cellular Analysis System	629029, B24465, B24802, B68304	All
UniCel DxH 600 Coulter Cellular Analysis System	B23858	All

Dear Beckman Coulter Customer,

Beckman Coulter is sending you this letter regarding a documentation update for the blast Suspect Messages. The Suspect message information has been updated for clarity.

ISSUE:	Beckman Coulter has determined that additional clarification for the Blast Suspect messages is necessary. In rare situations, the UniCel DxH 800 and DxH 600 Coulter Cellular Analysis System may not flag or detect blasts in some blood samples. This is due to limitations in the available technology as well as sample limitations.
IMPACT:	In these situations there could be a delay in the diagnosis and treatment of conditions associated with blasts in the peripheral blood.
ACTION:	<p>Please refer to the following modified information for both the Suspect messages and Limitation sections for the Differential:</p> <p>Blasts are detected, but not enumerated, by internal algorithms using acquired events, histogram and dataplot patterns, and sophisticated statistical methods for all available data for the sample analyzed. A standard trigger value or limit corresponding to enumeration on peripheral smear cannot be established because:</p> <ul style="list-style-type: none"> • Laboratories differ in their desired sensitivity to abnormal flagging and messaging. • Laboratories differ in their definition of blasts. • Mature and immature abnormal cell types may be identified as blasts. • Blasts can be rare events. <p>Blasts can represent a mixed population of cells often associated with specimen abnormalities that alter the white cell population’s pattern distribution in dataplots and histograms away from a normal distribution. The presence of blast cells may trigger other available Suspect messages, and not all blood samples that contain blasts may report a Suspect message.</p>



Beckman Coulter Diagnostics Limited
 Lismeehan
 O’Callaghans Mills
 Co. Clare
 Ireland

Telephone +353 65 6831100
 Fax +353 65 6831122
 e-mail infoie@beckman.com
 Website www.beckman.com

Directors: P Devereux, A Doyle, P Della Puca
 Registered in Ireland no 485929 Registered Office: 70 Sir John Rogerson’s Quay, Dublin 2, Ireland

	<p>A blast Suspect message is not diagnostic. The user should not rely upon instrument results alone to replace the need for manual microscopic review of blood samples, if indicated by other clinical and laboratory features of the patient. Further diagnostic procedures and clinical evaluation must be evaluated for diagnosis.</p> <p>Refer to your Instructions for Use REF B26647, Chapter 6, Data Review, Processing Results, for complete information on all available messaging and flagging options on the system.</p>
RESOLUTION:	<ul style="list-style-type: none"> • The Instructions for Use, B26647, will be updated to Revision AE in March 2017. • The modified information is temporarily available within the online ReadMefile, B44444AE, available on the Beckman Coulter website at https://www.beckmancoulter.com/wsrportal/page/techdocSearch

The national competent authority has been informed of this field safety corrective action.

Share this information with your laboratory staff and retain this notice as part of your laboratory Quality System documentation. If you have forwarded the affected product to another laboratory, please provide them with a copy of this letter.

Complete and return the enclosed response form within 10 days so that we are assured you have received this important communication.

If you have any questions regarding this notice, please contact the Customer Support Hotline at 00353 1407 3082 or techsupportie@beckman.com.

We apologize for any inconvenience to your laboratory.

Yours sincerely,



John Winter
Quality / Regulatory Affairs Manager - Northern Region Europe
Tel +44 (0) 1494 429184/+44 (0)1494 429128
Fax +44 (0) 1494 429182

Enclosed: Vigilance Response Form