Firma XXXXX  
XXX XXXXXX  
XXXXXXXXXXXXXXXX XX

XXXXX XXXXXXXX

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Ansprechpartner

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26. Januar 2015

-

## **Safety Information referring to TRILUX Medical OT Lighting systems with cameras and monitor holders**

Ladies and Gentlemen,

in line with our product monitoring it is our duty to inform you about any possible risks for patients and users, that might arise during handling of the a.m. products.

Trilux Medical received a customer complaint about a drop of an Aurinio LED OT light with the connected spring arm due to improper use. The improper use was to move the light over the stop what lead to a derailing of the pivot bearing.

This risk upon improper use is only applying to lights which are equipped with a camera system with stop as well as with monitor holders with stop. These products are amongst others equipped with stops at the connection point between the spring arm and the boom of the central axis to avoid a rotation of more than 360°.

Aurinio OT Lights without camera systems are not affected!

The stop protects the inside cabling to avoid a damage of same when rotating more than 360°.

An excessive use of physical strength which may result in the shearing-off of the stopper pin may result in a minor possibility for the OT light with the spring arm (monitor holder) to fall down.

To ensure that the light with the spring arm (monitor holder) may not fall down, even when the stopper pin might have sheared-off, we strongly recommend the retrofit of a mechanical locking device that excludes any possible risk even upon improper use.

Enclosed to this letter please find this necessary mechanical locking device as well as the respective installation guide therefore.

We ask for your respective feedback on the enclosed form sheet as soon as the retrofit of the OT lights (monitor holders) has been realized. This allows us to register the retrofit in our product records and to create a summary on all activities.

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Following our sales records you have received the products in question with the following documents:

Order Confirmation No.

8XXXXXXXX

Invoice No.

9XXXXXXXX

Please ensure that within your organization all users of the a.m. products as well as any other relevant persons involved are being informed about this urgent and important safety information.

Should you have delivered these OT light systems (monitor holders) to any third party, please make sure that this information including the enclosures is passed on to this party/user.

Please keep this information at least until this measure has been finalized.

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

## 2. Demontage Deckel von Z-Achse/Disassembly cover of Z-axis


Beide Befestigungsschrauben am Ende des Auslegers lösen (Abb. 2.2)  
Unscrew both fixing screws on the upper end of the boom. (Pic 2.2)

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Press both plastic covers apart and remove via (Pic 4.3) careful levering with one screwdriver (Pic 4.1/4.2) within both grooves.

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Through wiring carefully from the Pull spring arm hinge (Pic 5.1)

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Press plastic sleeve with a screwdriver flush (Pic 5.2) under the casting

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Pull out the plastic sleeve using the disassembly tools (Pic 5.3) from the spring arm.

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
Sleeve with side cutter.  
Unpick (Pic 5.5) and remove (Pic 5.6).

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
**Attention:** Install locking tab with radius down (Pic 6.2).

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Insert the Securing tab from bottom into the the spring arm, place it parallel to the upper surface, and hold it there (Pic 6.3/6.4).

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
Microencapsulated screw mount M 8,0x130,0 8.8 vz with lock washer and locking tab in the spring arm (Pic 6.5).

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Tighten cylinder bolt with a torque wrench 12,5 Nm (Pic 7.2).

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Mount panels and lid in reverse order, Make sure the correct position of the aperture plates in spring arm hinge. Then fit both plastic covers in sequence and press tight until they are heard to snap in place