

Urgent Field Safety Notice

MEDICAL DEVICE VOLUNTARY PRODUCT RECALL AND FIELD CORRECTION TO LASEREDGE® KNIVES, ACCORDING TO ATTACHED LIST OF LOT NUMBERS OF THE AFFECTED PRODUCTS

Thursday 23rd February 2017

Please forward this information to all users and staff who may use Laseredge® Knives -products manufactured by Angiotech (Surgical Specialties)


FSCA-identifier: CAC-2016-006 - LaserEdge® Knives / 15.02.2017

Type of Action: VOLUNTARY PRODUCT RECALL

Dear Valued Customer,

This is to inform you of a medical device voluntary product recall and field correction involving LaserEdge® knives, due to an increase in complaints of a dull knife edge.

See below an example of a LaserEdge® Knives product label for ease in identifying the product.



During the time period of January 2016 – January 2017, it has been determined through the complaint trending program that the LaserEdge® Surgical Knives may have demonstrated higher than normal complaints for dull knives. We are committed to ensuring that all of our products meet the highest standards of quality and take matters such as this very seriously, which is why we are taking this action.

Customer Services and Orders:



Visioncare
T 0845 602 2350
F 0845 602 2351

Surgical
T 020 8781 0000
F 020 8781 0001

Pharmaceutical
T 020 8781 2991
F 020 8781 0001

Aesthetics
T 0845 600 5212
F 0845 600 5215



If excessive force is required to push a dull knife through the cornea, this may result in:

- 1) Sub-optimal incision shape, such as short tunnel or lack of multi-plane beveling. The consequence may be incisions that are not watertight, requiring sutures, or inducement of corneal astigmatism.
- 2) Uncontrolled penetration through the cornea resulting in injury to anterior segment structures, such as iris, capsule, or lens.

It has come to our attention that some boxes of LaserEdge® Knives have not been as sharp as previous lots of this product. Please carefully review the notes outlined in this letter regarding your LaserEdge® Knives. This action represents a voluntary product recall and we have notified the appropriate authority of this voluntary recall.

According to our records, your facility may have a supply of LaserEdge® Knives according to attached list of Lot numbers of affected products.

Actions to be taken:

We ask that you please quarantine any unused boxes (full and partial) and take the following steps to return the product to Bausch + Lomb at our company's expense:

1. **Quarantine the product:** Please review your inventory and hold all unused (full and partial) boxes of LaserEdge® Knives according to attached a list of Lot numbers of affected products. Please note that only blades in the sealed blister can be returned.
2. **Return the product:** Please complete the enclosed **Recall and Field Correction Acknowledgement Form** and return it to the Customer Services Department at Bausch + Lomb, including all details for the pickup service by UPS (name of your contact person including email address, exact location in your facilities and estimated number of shipper boxes). Once we have received the completed **Recall and Field Correction Form**, we will be in contact to arrange collection of the identified products in your facilities.

This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

The decision to conduct this voluntary recall is part of our commitment to quality and customer satisfaction. We sincerely apologise for any inconvenience and assure you that we are working diligently to resolve this issue in a timely manner. Should you have any questions or concerns, please feel free to contact customer services on 020 8781 0000 or csd.backoffice@bausch.com.

Yours faithfully,

A handwritten signature in black ink, appearing to read 'Duncan Dow', written in a cursive style.

Duncan Dow

UK/I General Manager

Recall and Field Correction Acknowledgement Form

This is to acknowledge receipt of the above referenced recall and field correction notification dated **February 23rd, 2017**

Product Details:

LaserEdge® Knives (6/Box or individually)

A list of Lot numbers of affected products is attached to this letter

Please confirm inventory levels of the affected product at your facility with the 7-digit lot numbers:

Product	Lot #	# Received	# Used	# In Inventory	Responsible Person Initials

To obtain a Return Material Number (RMA) and arrange pick up of the identified products, please complete, sign and return this form to:

Fax: 020 8781 0001 / Email: csd.backoffice@bausch.com

Contact person for communication with UPS: _____

Contact person e-mail address: _____

Pickup location in your facilities: _____

Estimated shipper boxes to be returned: _____

- I hereby certify that I have quarantined the above listed product to prevent use and am awaiting pick up by a Valeant/Bausch + Lomb representative or agent.
- I hereby certify that I do not have inventory of the Lot numbers of affected products attached to this letter

Date

Name (Print)

Bausch + Lomb Account Number

Signature

Facility Name

Telephone Number

Confidentiality Agreement: The information contained in this message is privileged and confidential information intended for the use of the addressee listed above. If you are neither the intended recipient nor the employee or agent responsible for delivering this message to the intended recipient, you are hereby notified that any disclosure, copying, distribution or the taking of any action in reliance on the contents of the information is strictly prohibited. If you have received this in error, please immediately notify us by telephone to arrange for the return of the original document to us.