

January 12, 2017

To: Distributors

Subject: **URGENT MEDICAL DEVICE FIELD SAFETY NOTICE / INFORMATION**

Reference: **FA 2017-01 (ZFA 2016-164)**

Affected Product: **Revitan® Straight and Curved Distal Stems (List in attachment 2)**

Dear Sir or Madam,

This notice is to inform you that Zimmer GmbH is conducting a voluntary medical device Field Safety Notice for the Revitan Revision System. This Field Safety Notice is to provide information on the use of the Revitan Revision System and is jointly provided with the latest updated surgical techniques with references 06.01169.012X - Ed 2016.11 – rev. 03 for the Revitan Curved and 06.01109.012X - Ed 2016.11 – rev. 04 for the Revitan Straight¹.

Zimmer GmbH would like to emphasize that patient selection for modular hip revision systems, such as the Revitan system, is a key parameter. Modular stems have been used in revision surgery for many years and provide advantages for issues that are frequently encountered in a revision setting.² However, a modular system also comes with the risk for potential stem fractures.³ This is a rare complication of modular revision stems caused, most commonly, by large stresses to the modular junction when the proximal component lacks medial bone contact at the calcar.^{4,5} The risk of stem failure is not unique to any one modular system.^{6,7,8,9,10} Bone stock of adequate quality must be present and appraised at the time of surgery. For patients with severe proximal deficiency, a surgeon should consider surgical options to ensure proximal bone support (such as medial and/or lateral strut grafts) or switching to a monobloc revision stem, such as the Wagner SL Revision® Hip Stem. Please review the indications and contra indications in the associated surgical techniques and Instructions for Use of the Revitan Revision System.

¹ English version, other languages available.

² Murray, et al. Modularity in Revision Total Hips: The Use of the Distally-Fixed Stem. *Seminars in Arthroplasty*, 2008

³ Sekundiak, T., A Modular Hip System for Simplification of Revision Hip Arthroplasty. *American Journal of Orthopedics*, 2010; 39:7-12

⁴ Fink, B., Urbansky, K., & Schuster, P. (2014). Midterm results with the curved modular tapered, fluted titanium Revitan stem in revision hip replacement. *Bone Joint J*, 96-B(7), 889-895

⁵ Buttaro, M.A., Mayor, M.B., Van Citters, D. et al. Fatigue fracture of a proximally modular, distally tapered fluted implant with diaphyseal fixation. *Journal of Arthroplasty*, 22(2007), 780.

⁶ Kwong, L.M., Miller, A.J., Lubinus, P. A modular distal fixation option for proximal bone loss in revision total hip arthroplasty: a 2-to-6 year follow-up study. *Journal of Arthroplasty*, 18 (3 Suppl 1) (2003), S22.

⁷ Rodriguez, J., Fada, R. Murphy, S.B., et al. Two year to five year follow-up on femoral defects in femoral revision treated with the Link MP modular stem. *Journal of Arthroplasty*, 24 (2009), 94.

⁸ Slover, James, et al. "Fatigue failure of newer generation modular revision femoral stem following fracture healing: a case report." *Bulletin of the NYU Hospital for Joint Diseases* 73.1 (2015): 54.

⁹ Efe, T., Schmitt, J. Analyses of prosthesis stem failures in noncemented modular hip revision prostheses. *Journal of Arthroplasty* (2011); 665-667.

¹⁰ Skyttä, E.T., Eskeinen, A., Remes, V. Successful femoral reconstruction with a fluted and tapered modular distal fixation stem in revision total hip arthroplasty. *Scandinavian Journal of Surgery* (2012): 101; 222-226.

Potential Risks associated with the issue:

| Risks | | |
|---|---------------|---|
| Describe immediate health consequences (injuries or illness) that may result from use of or exposure to the product issue. | Most Probable | Worst Case |
| | Pain | Revision surgery due to breakage of the stem. Additional possible complications in order to achieve the removal of the broken prosthesis: Transfemoral approach necessary to remove fractured implant which leads to a significant prolongation of initial surgery time. An additional damage of tissue and bone quantity may be expected. |
| Describe long range health consequences (injuries or illness) that may result from use of or exposure to the product issue. | Most Probable | Worst Case |
| | Pain | Potential increased infection risk due to prolongation of surgery time and additionally a potential prolongation of healing time. The mobility phase of the patient will be decreased due to the fact that an immediate loading after a transfemoral approach is not applicable. An additional damage of tissue and bone quantity may be expected. |

Please note that the latest revision rates for the Revitan Revision System due to breakages, as monitored on regular basis by Zimmer GmbH are within acceptable rates.

Our records indicate you may have received one or more of the affected products.

Your Responsibilities

1. Review this notification and ensure affected team members are aware of the contents.
2. Inform your affected customers using the surgeon letter.
3. There are no specific patient monitoring instructions related to this Field Safety Notice that are recommended beyond your existing follow up schedule.
4. Complete Attachment 1 – Certificate of Acknowledgement.



- a. Return a digital copy to uk.fieldaction@zimmerbiomet.com or to your local Zimmer Biomet contact.
- b. Retain a copy of the Acknowledgement Form with your recall records in the event of a compliance audit of your documentation.

Other Information

This voluntary medical device Field Safety Notice was reported to all relevant Competent Authorities and the related Notified Body as required under the applicable regulations for Medical Devices.

Please keep Zimmer Biomet informed of any adverse events associated with this product or any other Zimmer Biomet product by emailing per.uk@zimmerbiomet.com or to your local Zimmer Biomet contact.

The undersigned confirms that this notice has been delivered to the appropriate Competent Authorities, as per MEDDEV 2.12-1 revision 8.

We would like to thank you for your co-operation in advance and regret any inconveniences caused by this Field Safety Notice.

Sincerely,

A handwritten signature in blue ink that reads 'Matthias Bürger'. The signature is written in a cursive style and is positioned above a horizontal line.

Matthias Bürger

Zimmer Biomet Vice President QARC EMEA



ZIMMER BIOMET

ATTACHMENT 1

Certificate of Acknowledgement

IMMEDIATE RESPONSE REQUIRED –TIME SENSITIVE ACTION NEEDED

Affected Product: _____ *[Enter affected product name from page 1 of the Field Safety Notice]*

Territory Number: _____ **Account Number:** _____

Account Name: _____

Account Address: _____

Certificate of Acknowledgement:

By signing below, I acknowledge that I have received, read, and understood the contents of this Field Safety Notice communication. All required activities are complete or are in the process of being completed.

Printed Name: _____ **Signature:** _____

Title: _____ **Tel:** () _____ - _____ x _____ **Date:** ___/___/___

Note: This form must be returned to Zimmer Biomet before this action can be considered closed for your account. It is important that you complete this form and email a copy to: uk.fieldaction@zimmerbiomet.com or to your local Zimmer Biomet contact: _____.



ATTACHMENT 2

ATTACHMENT 2- Affected Product List

| Description | Material |
|-------------------------------|--------------|
| REVITAN dist. straight 14/120 | 01.00405.014 |
| REVITAN dist. straight 14/140 | 01.00405.114 |
| REVITAN dist. straight 16/140 | 01.00405.116 |
| REVITAN dist. straight 18/140 | 01.00405.118 |
| REVITAN dist. straight 20/140 | 01.00405.120 |
| REVITAN dist. straight 22/140 | 01.00405.122 |
| REVITAN dist. straight 24/140 | 01.00405.124 |
| REVITAN dist. straight 14/200 | 01.00405.214 |
| REVITAN dist. straight 16/200 | 01.00405.216 |
| REVITAN dist. straight 18/200 | 01.00405.218 |
| REVITAN dist. straight 20/200 | 01.00405.220 |
| REVITAN dist. straight 22/200 | 01.00405.222 |
| REVITAN dist. straight 24/200 | 01.00405.224 |
| REVITAN dist. straight 26/200 | 01.00405.226 |
| REVITAN dist. straight 28/200 | 01.00405.228 |
| REVITAN dist. straight 16/260 | 01.00405.316 |
| REVITAN dist. straight 18/260 | 01.00405.318 |
| REVITAN dist. straight 20/260 | 01.00405.320 |
| REVITAN dist. straight 22/260 | 01.00405.322 |
| REVITAN dist. straight 24/260 | 01.00405.324 |
| REVITAN dist. straight 26/260 | 01.00405.326 |
| REVITAN dist. straight 28/260 | 01.00405.328 |
| REVITAN dist. curved 14/140 | 01.00406.114 |
| REVITAN dist. curved 16/140 | 01.00406.116 |
| REVITAN dist. curved 18/140 | 01.00406.118 |
| REVITAN dist. curved 20/140 | 01.00406.120 |
| REVITAN dist. curved 22/140 | 01.00406.122 |
| REVITAN dist. curved 24/140 | 01.00406.124 |
| REVITAN dist. curved 14/200 | 01.00406.214 |
| REVITAN dist. curved 16/200 | 01.00406.216 |
| REVITAN dist. curved 18/200 | 01.00406.218 |
| REVITAN dist. curved 20/200 | 01.00406.220 |
| REVITAN dist. curved 22/200 | 01.00406.222 |
| REVITAN dist. curved 24/200 | 01.00406.224 |
| REVITAN dist. curved 26/200 | 01.00406.226 |
| REVITAN dist. curved 28/200 | 01.00406.228 |
| REVITAN dist. curved 16/260 | 01.00406.316 |
| REVITAN dist. curved 18/260 | 01.00406.318 |
| REVITAN dist. curved 20/260 | 01.00406.320 |
| REVITAN dist. curved 22/260 | 01.00406.322 |
| REVITAN dist. curved 24/260 | 01.00406.324 |
| REVITAN dist. curved 26/260 | 01.00406.326 |
| REVITAN dist. curved 28/260 | 01.00406.328 |