DIRECTIONS FOR USE ANCHOR BOLT (LSB Series)

| | |
|---------|----|
| STERILE | EO |
| USA | |
| Rx ONLY | |

CE
2797

| | |
|----|-----|
| EC | REP |
|----|-----|




E.C. REP Ltd
5 Fitzwilliam Square East
Dublin 2, D02 R744, Ireland
Tel: +353 1 2 544 944

TABLE OF CONTENTS

| | |
|------------------|---------------|
| English | Pages 4 - 5 |
| Čeština | Pages 6 - 7 |
| Dansk | Pages 8 - 9 |
| Deutsch | Pages 10 - 11 |
| Español | Pages 12 - 13 |
| Ελληνικά | Pages 14 - 15 |
| Français | Pages 16 - 17 |
| Italiano | Pages 18 - 19 |
| Nederlands | Pages 20 - 21 |
| Turkish | Pages 22 - 23 |
| Portuguese | Pages 24 - 25 |

Please refer to Ad-Tech's Symbols Glossary @ www.adtechmedical.com

DIRECTIONS FOR USE ANCHOR BOLT (LSB Series)

To Insert:

- -
 -
 -
- -
 -
- -
 -
 -
 -
- -

To Remove:

-
-
-
-
-

DIRECTIONS FOR USE

ANCHOR BOLT

(LSB Series)

Purpose: The Ad-Tech Anchor Bolts are optional accessories for use with Depth Electrodes. The Anchor Bolts may be applied when it is desired to minimize concerns about potential cerebrospinal fluid (CSF) leakage and infection of the subdural space while stabilizing the electrode. Anchor Bolts are secured in the skull to provide an access point for Depth Electrodes.

Contraindications: Anchor Bolts should not be used on any patient whom the physician/surgeon considers at risk for infection or on whom the use cannot be performed safely.

The Anchor Bolt should not be used with patients that have softening of the skull or low skull bone density.

Use: The Anchor Bolt (LSB Series) is supplied STERILE. Sterility is guaranteed unless package is damaged or seal is broken.

This product should only be used by a physician/surgeon trained in the use of Anchor Bolts.

Do not apply excessive torque to the anchor bolt during placement. The anchor bolt should introduce easily and SHOULD NOT BE FORCED. Excess torque may result in a bent or broken anchor bolt.

Anchor Bolt is for Single Use Only. Do Not Re-Sterilize or Reuse the Anchor Bolt.

Not Intended for Implantation (21 CFR 860.3(d): ≥ 30 days). For Surgical Use Only. Do not use if packaging is damaged.

The LSB Anchor Bolt Placement and Removal Wrench is supplied sterile and may be re-used. Reference AT10054-B Cleaning and Sterilization Instructions Placement Tools and Accessories to support re-use.

The Anchor Bolts are compatible with the recommended Styles of Ad-Tech Depth Electrodes and recommended cranial drill bit diameter hole to support placement.

| Anchor Bolt Model (Gasket Color) | Recommended Depth Electrode Style | Recommended Cranial Drill Bit Diameter for Anchor Bolt Placement |
|-------------------------------------|--------------------------------------|---|
| LSBK1-AX-05 (Green) | SD | 2.4 mm |
| LSBK1-AX-06 (Clear) | BF, LD, MM | 2.4mm |
| LSBK1-BX-04 (Blue) | RD | 2.4 mm |
| LSBK1-BX-05 (Green) | SD | 2.4 mm |
| LSBK1-BX-06 (Clear) | BF, LD, MM | 2.4 mm |
| LSBK1-CX-04 (Blue) | RD | 2.4 mm |

DIRECTIONS FOR USE ANCHOR BOLT (LSB Series)

| Anchor Bolt Model (Gasket Color) | Recommended Depth Electrode Style | Recommended Cranial Drill Bit Diameter for Anchor Bolt Placement |
|-------------------------------------|--------------------------------------|---|
| LSBK1-CX-05 (Green) | SD | 2.4 mm |
| LSBK1-CX-06 (Clear) | BF, LD, MM | 2.4 mm |
| LSBK2-AX-04 (Blue) | RD | 2.4 mm |
| LSBK2-BX-04 (Blue) | RD | 2.4 mm |
| LSBK2-CX-04 (Blue) | RD | 2.4 mm |

In the event that an anchor bolt failure is experienced (i.e. breakage), proceed with the following steps:

1. Stop use immediately and monitor the patient.
2. Attempt to remove the retained portion of the bolt using the LSB Anchor Bolt Wrench.
3. If unable to use the wrench for removal, use other medical equipment available (e.g. clamp) to remove the retained portion.
4. If still unable to remove the retained portion of the bolt, additional surgical intervention may be required.

CAUTION: Do not use Ad-Tech Anchor Bolts if third party interface (e.g. frame) is damaged or flawed in any way.

CAUTION: Discontinue placement of anchor bolt if you experience any instrument (e.g. wrench, drill, etc.) bending during the placement procedure.

CAUTION: Discontinue use of anchor bolt if you experience bending or breakage of the anchor bolt during the placement procedure or during use.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

CAUTION: Reuse of this device is prohibited as it may malfunction and cause contamination and risk to the patient.

CAUTION: Disconnect from monitoring equipment during cardiac defibrillation.

CAUTION: Protect the Anchor Bolt from blunt impact that may occur if the patient has a seizure. Significant force may fracture the Anchor Bolt.

STORAGE: LSB Anchor Bolts need to be stored and used within typical hospital/office room ambient temperature and humidity conditions.

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400 WEST OAKVIEW PARKWAY, OAK CREEK, WI, 53154 U.S.A
Toll Free (U.S.A): 800.776.1555