

URGENT: FIELD SAFETY NOTICE
NanoClave Leaks
See Table 1 for Affected Product and Lot Numbers

02 July 2020

Dear Valued Customers:

Director of Risk Management
Director of Nursing
Director of Materials Management

ICU Medical, Inc. is issuing this Urgent Field Safety Notice to notify you of a potential lot-specific defect with the NanoClave connector. This Urgent Field Safety Notice letter details the issue and the required steps for you to perform.

Issue:

ICU Medical has identified the potential for a manufacturing defect on the internal surface of the NanoClave within specific lots of NanoClave sets, which may inhibit a proper seal with the NanoClave spike.

Potential Risk:

Inadequate sealing with the NanoClave spike may potentially cause fluid leak, blood loss, air ingress, or contamination. To date, ICU Medical has not received any reports of serious injury or death associated with this issue.

Affected Product:

Our records indicate that you may have received some of the affected products, which were distributed in Ireland in February 2020. The affected item and lot numbers are provided in Table 1.

Required Actions for Users:

- 1) Please discontinue the use and distribution of the affected product immediately. Check your inventory and quarantine all affected product at your facility.
- 2) Inform potential users of the product in your organization of this notification and complete the attached response form. Return the completed response form to the e-mail address on the form, even if you do not have the affected product.
- 3) Upon receipt of the completed response form and return of the affected product, ICU Medical will credit you for any product returned. You will only receive credit for product that you return. NOTE: Credits for product purchased through distributor will be credited by the distributor.
- 4) If you have distributed the product further, immediately notify your accounts that received the product identified in the Affected Product / Table 1 sections of this notification and ask them to complete a response form and return to the e-mail address on the form.

Follow up Actions by ICU Medical:

Please contact Customer Service using the information provided below for assistance reordering replacement product.

For further inquiries, please contact ICU Medical using the information provided below.

ICU Medical Contact	Contact Information	Areas of Support
Global Complaint Management	ProductComplaintsBucharest@icumed.com	To report adverse events or product complaints
ICU Customer Service	EMEA distributor-support@icumed.com	Additional information or assistance

The HPRA has been notified of this action.

ICU Medical is committed to patient safety and is focused on providing exceptional product reliability and the highest level of customer satisfaction. Thank you for your prompt support on this important matter. We appreciate your cooperation.

Sincerely,



Corine Broekhuizen
Director, Quality and Regulatory Affairs
ICU Medical BV

Enclosures:

- Affected Product and Lot Numbers
- Response Form

Table 1: Affected Product and Lot Numbers

Item No.	Description	Lot Numbers
011-A1000	NanoClave™ Connector	4476111

URGENT: FIELD SAFETY NOTICE RESPONSE FORM

NanoClave Leaks

See Table 1 for Affected Product and Lot Numbers

02 July 2020

Please check your inventory and complete the information below, even if you do not have the affected product. *Failure to complete all sections of this page may result in improper, delayed or denied credit.*

Please return the completed form to EMEA-Quality@icumed.com, EMEA distributor-support@icumed.com or your ICU Medical sales representative.

Hospital/Facility Name	
ICU Medical Customer # (if applicable)	
Address/City/ Postal code	
Contact Name/Title/Phone/E-mail Address	
Name and Title of Person Completing this Form	
Completed by: Printed Name/Signature/Date	

I have **NO** affected product (complete and return this form to the e-mail addresses above).

YES, I have affected product, I have followed the instruction provided to me and I am going to contact EMEA distributor-support@icumed.com to make arrangement to return the affected products.

If affected product is not being returned, please explain below:

- Have you distributed the product further to the retail level? YES___ NO___
 - If yes, have you notified your retail customers? YES___ NO___ (if no, explain below)

If you have distributed the product further, please provide the list of your retail customers, inclusive of customer name, address, city, state, zip code, telephone number and quantity of product distributed along with your completed response form to the contact information listed above so ICU Medical can verify effectiveness of the recall notification to the appropriate level.

Lot Number	Quantity to be returned	Wholesaler/Distributor Name If you purchased from Wholesalers/Distributors include name, address, city, state, zip, quantity from each, and invoice number. If you purchased directly from ICU Medical leave this section blank.	PO, debit memo or invoice
		1.	
		2.	

Adverse events and complaints associated with the use of these products should be reported and emailed to HPRA or to the ICU Medical at the contact information provided.