



## Urgent Field Safety Notice Device Correction

24<sup>th</sup> May 2019

Dear Sir/Madam

**Affected Product**

Product Code	Product Family	Serial Numbers
114870	PRISMAFLEX 7.XX ROW	All
955052	PRISMAFLEX 8.XX ROW	All
114489	PRISMAFLEX 6.10 ROW	All

**Problem Description**

Baxter Healthcare has received reports of Prismaflex 8.10 devices with inactive syringe pump during Continuous Renal Replacement Therapy (CRRT) treatment while using RCA (Regional Citrate Anticoagulation). It was determined that there is a potential for the calcium syringe pump to be inactive without the device alarming after completed change syringe procedure.

Baxter will be upgrading all Prismaflex devices with software version 8.10 to software version 8.20. The new software version will include an enhancement to ensure alarm generation when the issue occurs while performing therapy using Regional Citrate Anticoagulation (RCA)

**Hazard Involved**

An inactive pump may result in under delivery of calcium leading to hypocalcemia. Hypocalcemia may potentially result in serious adverse health consequences. There have been two reports of serious injury associated with this issue.

**Action to be taken by the user**

Baxter is kindly asking that you take the following actions:

1. Clinicians may continue to safely use the Prismaflex devices while utilizing additional caution to ensure that the syringe pump operates as intended after the change syringe procedure, until the software upgrade can be performed.
2. Clinicians should use appropriate syringes when using the Prismaflex device per the operator's manual, see section 15.6.2 Citrate – calcium method.

3. A local Baxter service representative will contact your facility to arrange for the software upgrade for all Prismaflex devices with current 8.10 software. Your facility will be receiving this software upgrade from Baxter at no charge. This will be performed during scheduled servicing.
4. **If you purchased this product directly from Baxter, complete the enclosed Baxter Customer Reply Form and return it to Baxter** by faxing it to 01 206 5577 or scanning and e-mailing it to qa\_dublin@baxter.com **even if you do not have any inventory.** Returning the customer reply form promptly will confirm your receipt of this notification and prevent you from receiving repeat notices.
5. If you purchased this product from a distributor, please note that the Baxter customer reply form is not applicable. If a reply form is provided by your distributor or wholesaler, please return it to the supplier per their instructions.
6. If you distribute this product to other facilities or departments within your institution, please forward a copy of this communication to them.
7. If you are a dealer, wholesaler, distributor/reseller, or original equipment manufacturer (OEM) that distributed any affected product to other facilities, please notify your customers of this Urgent Medical Device Correction in accordance with your customary procedures.

**Further  
information  
and support**

If you have additional questions, please contact your Baxter sales representative.

We apologize for any inconvenience this may cause you and your staff. Baxter's software version update will take additional measures to further ensure patient safety. Baxter is committed to ensuring our products and services consistently meet the highest standards of quality and safety for our patients and healthcare providers.

The Local MOH has been informed about this action.

Sincerely,



Sam Nickerson  
Business Unit Head, Acute Therapies, UK and Ireland  
Baxter Healthcare Ltd.



Enclosure: Baxter Customer Reply Form  
Attachment 1: Affected Products Table