Medtronic Ireland Limited

T Block 3090-3094
Lake Drive
Citywest Business Campus
Dublin
DN24 XN47
Ireland
Tel: 01 511 1400
Fax: 01 807 7220
www.medtronic.ie

URGENT FIELD SAFETY NOTICE

**Mahurkar™ Acute Triple Lumen Catheters and
Mahurkar™ Acute High Pressure Triple Lumen Catheters
Recall**

June 2023

Medtronic Reference: FA1333

EU MDR: EU Manufacturer Single Registration Number (SRN): US-MF-000028763

Dear Risk Manager/Healthcare Professional:

The purpose of this letter is to advise you that Medtronic is voluntarily initiating a recall for specific lots of **Mahurkar™ Acute Triple Lumen Catheters** and **Mahurkar™ Acute High Pressure Triple Lumen Catheters**. **Please note:** This recall does **not** include any Mahurkar **Elite** Catheters.

You are receiving this letter as Medtronic records indicate your facility may have at least one of the Mahurkar™ Acute Triple Lumen Catheters and/or at least one of the Mahurkar™ Acute High Pressure Triple Lumen Catheters outlined in Attachment B. Medtronic initiated this action to prevent the use of potentially affected product that may impact patients.

Issue Description:

During manufacturing related testing, the catheter center lumen was found to have an occlusion in the tip of the catheter. Investigation identified the source of the occlusion as excessive MDX, a silicone-based lubricant which coats the catheter tip. As of 10 June 2023, there have been zero (0) confirmed complaints. Additionally, there have been zero (0) reported adverse events and there have been zero (0) reported deaths.

Risk to Health:

An incorrect application of MDX to catheters may result in the hazardous situation whereby the catheter is occluded, partially or fully, and/or uncured or excessive MDX may dislodge from the catheter. An occurrence of the hazardous situation may lead to potential harms identified as full catheter obstruction resulting in delay to treatment and partial obstruction resulting in reduced flow or particulate dislodgement that may result in delay to treatment, hemolysis, embolism/embolus or thrombosis/thrombus.

Patient Recommendation:

Mahurkar™ Acute Triple Lumen Catheters and Mahurkar™ Acute High Pressure Triple Lumen Catheters are intended for short term use of up to 29 days. For patients with affected lot(s) currently in place, a replacement procedure is recommended. If a patient is found to have a catheter from an affected lot, the patient's medical team should assess the overall patient risk when considering the timing of a replacement. Clinicians should continue to follow current product Instructions For Use (IFU) along with facility specific policies and procedures for routine assessment of the hemodialysis access device for patency, function, and efficacy.

Required Actions:

1. Immediately quarantine and discontinue use of all unused Mahurkar™ Acute Triple Lumen Catheters and Mahurkar™ Acute High Pressure Triple Lumen Catheters referenced in Attachment B - List of affected Lot numbers (see Attachment A for guidance to identify impacted product).
2. To help you identify if you have affected product, please visit our website www.Medtronic.com/Mahurkar-Triple-Lumen-Catheter-Recall. Here you will find a tool to help you determine if the product you have is affected by this recall.

Please note: This recall does **not** include any Mahurkar **Elite** Catheters.

3. Return all unused affected product(s) to Medtronic. Your Medtronic Sales Representative can assist in returning any affected product.
4. Please complete the enclosed Customer Acknowledgement Form and email to rs.regulatoryuk-ire@medtronic.com
5. This notice should be passed on to all those who need to be aware within your organization or to any organization including but not limited to Nephrologists, Intensivists, physicians, renal nurses, critical care nurses, or other dialysis staff where the potentially affected devices have been transferred. Please maintain a copy of this notice in your records.

Additional Information:

Medtronic has notified the Competent Authority of your country of this action.

We regret any inconvenience this may cause. We are committed to patient safety and appreciate your prompt attention to this matter. If you have any questions regarding this communication, please contact our Representative directly or via Tel No: 01 511 1400.

Sincerely,

Bethany Moxon

Associate Regulatory Affairs Specialist - UK and Ireland

Enclosures:

Attachment A: Identifying Affected Product

Attachment B: List of Affected Lot Numbers

Customer Confirmation Form

Attachment A: Identifying Affected Product

The **MAHURKAR™* Triple Lumen Catheter** is a radiopaque, polyurethane tube with two clear silicone catheter extensions and three internal lumina distinguished by color-coded adapters on the extensions:

- Red adapter (arterial)
- Blue adapter (venous)
- Clear adapter (medial)

The 12 Fr/Ch catheter is available in various implantable lengths as shown on the box and device labels. The **MAHURKAR™* Triple Lumen Catheter** is intended for short term central venous access for hemodialysis, apheresis, and infusion.

The **MAHURKAR™* Acute High Pressure Triple Lumen Catheter** is a radiopaque, polyurethane tube with two clear silicone catheter extensions and one clear polyurethane Infusion lumen. The three internal lumina can be distinguished by the color-coded luer-lock adapters on the silicone rubber extensions:

- Red adapter (arterial)
- Blue adapter (venous)
- Clear adapter (medial)

The proximal lumen provides "arterial" outflow from the patient; the distal lumen provides "venous" return; the medial lumen is for infusion of fluid, blood products, medications, blood sampling, pressure injection of contrast media and central venous pressure monitoring.

Attachment A: Identifying Affected Product

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© 2011 Covidien. Made in Costa Rica.
MAHURKAR is a trademark of Saktharam D. Mahurkar, M.D., used under license.