

Smiths Medical International Ltd.  
St Crispin Way, Haslingden, Rossendale  
Lancashire, BB4 4PW, UK  
T: +44 (0)1706 233800  
F: +44 (0)1706 218834  
www.smiths-medical.com

Our Ref: CAR - 172

Date: 26<sup>th</sup> February 2008

## **IMPORTANT PRODUCT INFORMATION**

### **PRODUCT:**

Smiths Medical Logical, NovaTrans, and Transtar Invasive Pressure Monitoring Kits with Trigger Flush Device, Product Codes DPSXXXXXX / MX95XXXX / MX96XXXX / MXXXX / SXXXXXX

### **DETAILS OF AFFECTED DEVICES:**

All product codes above manufactured between October 2007 and February 2008

### **TYPE OF ACTION:**

Field Safety Corrective Action

### **ATTENTION:**

Risk/ Safety Managers, Distributors, Clinicians and other Customers

### **Dear Customer,**

Smiths Medical International has received reports that some pressure monitoring kits allow flow rates which exceed the 3ml/ hour nominal flow rate specification. In the affected devices, the flow rate can be between the nominal flow rate, of 3ml/ hour, and the flush rate. This increased flow rate allows for more fluid than intended to be delivered to the patient and the higher flow rate can cause a positive offset in the pressure measurement indicated on the patient monitor leading to an inaccurate pressure reading.

Smiths Medical International recommends that you take action(s) as detailed in the attached "Field Safety Corrective Action Notice".

If after reading the Field Safety Corrective Action Notice you still have questions or queries relating to this matter please do not hesitate to contact: -

**UK Customer Services** – Tel: 01923 241411

**International Customer Services** – Tel: +44 (0)1303 260551

**Quality Assurance Department**- +44 (0)1706 233831

Yours Sincerely

**Omar Khan**

**QA Manager**

**Smiths Medical International**

St. Crispin Way, Haslingden

Rosendale, Lancs. BB4 4PW

Tel: +44 (0)1706 233831

E-Mail: omar.khan@smiths-medical.com

Smiths Medical International Limited  
St. Crispin Way,  
Haslingden, Rossendale,  
Lancashire, BB4 4PW UK  
T: +44 (0)1706 233800  
F: +44 (0)1706 218834  
[www.smiths-medical.com](http://www.smiths-medical.com)

## Urgent Field Safety Notice

---

**Affected Devices:** Smiths Medical Logical, NovaTrans, and Transtar Invasive Pressure Monitoring Kits with Trigger Flush Device, Product Codes DPSxxxxxx / MX95xxxx / MX96xxxx / MXxxx / SXxxxxxx

**Type of Action:** Field Safety Corrective Action

**Manufacturers Reference:** CAR-172

**Date:** 21<sup>st</sup> February 2008

**Details on affected devices:** All product codes above manufactured between October 2007 and February 2008

**Complete and Fax promptly to:** Omar Khan, Quality Manager , Smiths Medical International Ltd.

Fax No.: +44 1706 211878

SECTION 1 -:

Distributor Name: \_\_\_\_\_  
Printed Name: \_\_\_\_\_  
Country: \_\_\_\_\_  
Tel No.: \_\_\_\_\_ Ext: \_\_\_\_\_ Fax No.: \_\_\_\_\_ E-mail: \_\_\_\_\_

Tick one of the  
Boxes:

We have not received any products covered by the Field Safety Notice

We have passed a copy of this notice on to all customers in receipt of items manufactured during the affected period

SECTION 2 - ANSWER ALL ITEMS BELOW AND SIGN.

Has Regulatory Authority been notified?  YES  NO  
If yes, please indicate name of regulatory authority:  
If no please indicate why notification hasn't occurred

Questions should be directed to Omar Khan at +44 1706 211878 (e-mail [omar.khan@smiths-medical.com](mailto:omar.khan@smiths-medical.com))

Submitted By ..... Date:...../...../.....  
(Name). Signature:

**For Smiths Medical Use Only**

Customer Account No.