



S  
A  
F  
E  
T  
Y

NOTICE

# Cardiac Science's Powerheart and CardioVive Automated External Defibrillators

**IMB Safety Notice: SN2012(11)**  
**Circulation Date: 03 August 2012**

**MANUFACTURER/SUPPLIER**

Cardiac Science

**TARGET GROUPS**

General Practitioners  
General Public  
Community First Responder schemes  
Risk Managers  
Hospitals  
Nursing Homes  
Schools  
Sports clubs

**Please bring this safety notice to the attention of all who need to be aware of it.**

**ISSUE**

The Cardiac Science Powerheart and CardioVive AEDs affected by two separate issues may not deliver therapy.

**BACKGROUND**

The manufacturer Cardiac Science is currently conducting two Field Safety Corrective Actions (FSCAs) on the Irish market relating to their Powerheart and CardioVive automated external defibrillators (AEDs). The issues were communicated separately in Field Safety Notices (FSN) distributed in November 2009 and January 2012 (see attached).

[FSCA November 2009 - Cardiac Science Ref MDD09.070](#)

The first issue relates to Powerheart 9300A, 9300C, 9300D, 9300E, 9300P, 9390A, 9390E, CardioVive 92531, 92532 and 92533 automated external defibrillators (AEDs) manufactured between August 2003 and August 2009. These AEDs have resistors which may fail and the failure will not be detected by the AED's periodic self-test. If the issue were to occur the AED may not be able to deliver therapy during a rescue attempt.

# SAFETY NOTICE

The manufacturer distributed a follow up FSN in June 2010 to notify users that a mandatory software update to address this issue was available for all affected models. To date, despite communications distributed by Cardiac Science and their representatives, not all customers have carried out the software update on the affected devices in Ireland.

## FSCA January 2012 – Cardiac Science Ref MDD12.005

The second issue relates to Powerheart 9300A, 9300E, 9300P, 9390A, and 9390E AEDs manufactured between July 1, 2011 and December 30, 2011. In addition, certain CardioVive 92532 and 92533 model AEDs serviced during this time are also affected. The affected AEDs contain a circuit board manufactured with a component that may fail unexpectedly due to a supplier manufacturing defect. If the component were to fail during a rescue attempt, the AED may not deliver therapy.

Cardiac Science informed the IMB that they distributed the communications to all of their direct sales customers in Ireland. Please see accompanying FSNs for additional information. As some users in Ireland may have been supplied with affected devices indirectly, they may not be aware of these two issues. Cardiac Sciences advises that such customers should contact them directly if they identify that they have such a device.

### **ACTION OR RECOMMENDATIONS**

1. Ensure that the relevant personnel in your organisation are made aware of the two issues.
2. Determine if you have purchased a device that is affected by either of the above issues by reviewing the attached FSNs and lists of affected AED model numbers and serial numbers.
3. If you have an affected device, please follow the instructions provided by Cardiac Science in the relevant FSN.

### **ENQUIRIES**

Enquiries to the manufacturer should be addressed to:

Bob Odell  
Cardiac Science  
3303 Monte Villa Parkway Bothell  
WA 98021-8969  
Telephone: +1 425 402 2482  
Email: [aed210@cardiacscience.com](mailto:aed210@cardiacscience.com)

Or to Cardiac Science's European Authorised Representative:

MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany  
Tel.: +49 511 6262 8630  
Fax: +49 511 6262 8633  
Email: [vigilance@mdss.com](mailto:vigilance@mdss.com)

# S A F E T Y

## NOTICE

All adverse incidents relating to a medical device should be reported to the:

Irish Medicines Board  
Kevin O'Malley House  
Earlsfort Centre  
Earlsfort Terrace  
Dublin 2

Telephone: +353-1-6764971  
Fax: +353-1-6344033  
E-mail: [vigilance@imb.ie](mailto:vigilance@imb.ie)  
Website: [www.imb.ie](http://www.imb.ie)