

# Safety Notice

## Medical Devices

---

### Life-Point Pro AED / bPLUS / mPLUS defibrillator devices and Adult Electrode Pads & Paediatric Electrode Pads

#### Priority 1 – Immediate Action

HPRA Safety Notice: SN2016(46)

Issue Date: 21<sup>st</sup> December 2016

MANUFACTURER / SUPPLIER	HPRA CASE REFERENCE
METsis Medikal / LifePoint	V29146 / V30286

#### ISSUE

METsis Medikal has become aware of an issue whereby the green status indicator on the Life-Point Pro AED / bPLUS / mPLUS may not assure that the device is ready to use.

A defibrillator in this condition may be inoperable and unable to deliver a shock. This may result in therapy not being delivered to a patient.

The HPRA has also become aware that some METsis Medikal electrode pads produced by a specific supplier (Medeks Medikal) in 2015 may degrade prematurely, before their expiry date. Electrode Pads manufactured from the 1<sup>st</sup> of January 2015 to the 1<sup>st</sup> of June 2016 are being recalled as a result of this issue.

#### ACTION OR RECOMMENDATIONS

The HPRA advises the following;

All users/owners of Life-Point Pro AED / bPLUS / mPLUS AEDs should:

- 1) Ensure that you read and follow the instructions provided in the manufacturer's field safety notice (FSN) (see attached).
- (2) Ensure that you download the software upgrade and updated user manual for these devices from the manufacturer's website or contact your distributor to arrange the software upgrade.
- (3) Acknowledge receipt of the FSN if you have not already done so.
- (4) Report any concerns regarding these devices to the manufacturer and the HPRA.

All users/owners of lifePOINT Adult Electrode Pads / Paediatric Electrode Pads should:

- (5) Identify affected products by referring to the attached FSN.
- (6) Cease use and/or distribution of stock from the affected product batches, quarantine and return to manufacturer immediately.
- (7) Acknowledge receipt of the FSN if you have not already done so.
- (8) Forward this HPRA safety notice and attached FSNs to all those who need to be made aware within your organisation or to any organisation/person where these devices have been transferred.
- (9) Report any concerns regarding these devices to the manufacturer and the HPRA.

**The HPRA advises users to exercise extreme caution when using Life-Point Pro AED / bPLUS / mPLUS defibrillator devices. The HPRA is continuing to investigate both of these issues with the manufacturer.**

TARGET GROUPS	
Community First Responder Schemes County Councils Clinics Educations Centres Emergency First Responders Emergency Medical Technicians Fire Services Hospitals HSE Ambulance Services Medical Directors	Nursing homes Paramedics / Advanced Paramedics Private Ambulance Services Private Medical Practitioners Risk Managers Schools Sports clubs Supplies Officers / Managers Town Councils Voluntary / Auxiliary Ambulance Services

## BACKGROUND

METsis Medikal has issued a FSN advising customers that the green status indicator on the Life-Point Pro AED / bPLUS / mPLUS defibrillator devices may not signify that the device is ready to use. The status indicator of the device may be fixed "GREEN" as a result of mechanical and electrical failures.

METsis Medikal has advised in the FSN that the status indicator on the device is an advisory indicator and that the device should be checked manually by the user to determine if it is operable / functional by pushing the On/Off button as frequently as possible.

METsis Medikal, to address this issue, has revised the user manual and also developed a software upgrade which can be downloaded from the manufacturer's website free of charge.

The HPRA has also been advised that certain electrode pads produced by a specific supplier (Medeks Medikal) are being recalled due to the deterioration of certain components, which may lead to a serious operational problem. Please see the manufacturer's FSN (attached) for advice on how to identify impacted products.

#### MANUFACTURER / DISTRIBUTOR CONTACT INFORMATION

Enquiries to in relation to **FSN 1** should be addressed to:

METsis Medikal	Telephone:	+90 (312) 3862190
2B-1 ODTU	E-mail:	<a href="mailto:info@metsismedikal.com">info@metsismedikal.com</a>
Teknokent Binasi	Website:	<a href="http://www.metsismedikal.com">www.metsismedikal.com</a>
Ostim		
Yenimahalle / Ankara		
Turkey		

Enquiries to in relation to **FSN 2** should be addressed to:

LifePoint	Telephone:	+34 619 121 295
	E-mail:	<a href="mailto:lifepoint@lifepoint.es">lifepoint@lifepoint.es</a>
	Website:	<a href="http://www.lifepoint.es">www.lifepoint.es</a>

#### HPRA CONTACT INFORMATION

All **adverse incidents** relating to a medical device should be reported to:

Health Products Regulatory Authority	Telephone:	+353-1-6764971
Kevin O'Malley House	Fax:	+353-1-6344033
Earlsfort Centre	E-mail:	<a href="mailto:devicesafety@hpra.ie">devicesafety@hpra.ie</a>
Earlsfort Terrace	Website:	<a href="http://www.hpra.ie">www.hpra.ie</a>
Dublin 2		