

Safety Notice

Medical Devices

ARCHITECT c4000 Mixer ARCHITECT c8000 Mixer ARCHITECT c16000 Mixer

Priority 2 – Warning

HPRA Safety Notice: SN2019(14)

Issue Date: 17 May 2019

MANUFACTURER / SUPPLIER	HPRA CASE REFERENCE
Abbott Laboratories	V39936

ISSUE

Abbott has issued a field safety notice (FSN) to raise awareness of the potential for the mixer blade to separate from the mixer of ARCHITECT c4000, c8000 and c16000 analysers. A missing mixer blade will result in inadequate mixing of the reaction mixtures which could potentially lead to incorrect results (failed assay calibration, Quality Control (QC) out of range, and/or incorrect patient results).

Abbot advises users of these analysers to inspect the mixer blade as part of routine maintenance. Users who identify a mixer with a detached blade should consult their medical director regarding previously generated results.

Please refer to the accompanying FSN for further details.

ACTION OR RECOMMENDATIONS

The HPRA advises that healthcare professionals:

1. Refer to the accompanying FSN and follow the instructions provided.

2. Acknowledge receipt of the FSN if you have not already done so.
3. Forward a copy of this safety notice and the FSN to all relevant personnel within your organisation and to any other organisations/persons to which/whom these devices have been transferred.
4. Report any concerns regarding these devices and this issue to the manufacturer and the HPRA.

TARGET GROUPS

Biochemistry Departments	Laboratory Technicians
Clinical Chemistry Departments	Medical Scientists
Clinical Directors	Purchasing Managers
Hospital Managers / CEOs	Risk Managers
Laboratory Managers	Supplies Managers
Laboratory Staff	

BACKGROUND

Abbott has issued an FSN to provide further information to enable users to inspect the mixer blade and identify error codes which signal the need for an immediate inspection of the mixer. A procedure for inspecting the mixer is included alongside further precautions for the weekly maintenance of the mixer

Abbott has indicated that the investigation is ongoing and that updated information on these procedures will be provided once the investigation is completed.

Please refer to the accompanying FSN for further details.

The HPRA is issuing this safety notice to raise awareness of this issue.

MANUFACTURER CONTACT INFORMATION

Enquiries to the **manufacturer** should be addressed to:

Abbott	Abbott Helpline
Max-Planck-Ring 2	Telephone: +353-1-4691568
Wiesbaden 65205	E-mail: AbbottCareIrl@Abbott.com
Germany	

HPRA CONTACT INFORMATION

All **adverse incidents** relating to a medical device should be reported to:

Health Products Regulatory Authority	Telephone: +353-1-6764971
Kevin O'Malley House	Fax: +353-1-6344033
Earlsfort Centre	E-mail: devicesafety@hpra.ie
Earlsfort Terrace	Website: www.hpra.ie
Dublin 2	