

July 1, 2015

Important Information about Animas[®] Vibe[®] Insulin Pump

Dear Valued Health Care Professional:

Animas takes our commitment to notify patients of issues regarding their pumps seriously. We have received complaints concerning faded, dim or discolored displays on Animas[®] Vibe[®] Insulin Pumps manufactured prior to February 2013. The display may fade or become dim over time as a result of exposure to extremes of temperature and humidity.

We need to inform you that we are notifying patients who are still using an Animas[®] Vibe[®] Insulin Pump produced before February 2013. If you have any questions on the notification, please do speak to your Animas Representative or reach out to us on the contact details at the top of this letter.

What is the issue?

The display on the Animas[®] Vibe[®] Insulin Pump may fade over time, making it difficult to read the information on the screen.

What are the health and safety risks that may result?

If the display is faded dim or discolored, the patient may not be able to read the screen, which could prevent the patient from using the pump correctly.

The Owner's Booklet includes instructions for patients to contact Animas Technical Assistance if they suspect the pump has been damaged.

We appreciate your prompt attention to this issue.

Sincerely,

Animas Corporation