



## Urgent Field Safety Notice

### Medtronic MiniMed Infusion Sets - Recall of specific Lot Numbers Potential over-delivery of insulin

September 2017

Medtronic reference: FA784

Dear Distributor,

The purpose of this letter is to notify you that we are voluntarily recalling specific lots of infusion sets used with Medtronic insulin pumps.

Medtronic has become aware of recent reports of potential over-delivery of insulin shortly after an infusion set change that have required medical intervention. Our investigation has shown this can be caused by fluid blocking the infusion set membrane during the priming/fill-tubing process. A membrane blocked by fluid most likely occurs if insulin, alcohol, or water is spilled on the top of the insulin reservoir which then could prevent the infusion set from working properly. Infusion sets currently being shipped by Medtronic contain an enhanced membrane material that significantly reduces this risk.

We request you to inform your patients of this recall, using the attached letter, to ask them to do the following:

- A. Patients are instructed to go to [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look) and determine if they have recalled infusion sets. The website will prompt them to enter the REF and LOT numbers for all infusion sets in their possession. The website will then tell them which infusion sets are part of this recall and which are not.

The REF and LOT numbers are listed on the labels as shown in the below examples:



- B. Medtronic recommends that the recalled infusion sets are not used by patients.
  - a. If they have new and enhanced infusion sets that are not part of this recall, they should use only those sets starting with their next set change. As a reminder, we have enclosed Key Steps regarding the priming/fill-tubing process.
  - b. If they only have recalled infusion sets right now, it is very important that they carefully follow the Key Steps.
- C. Patients are instructed to discard their recalled infusion. Medtronic will replace the recalled infusion sets free of charge.


# Medtronic

In order to check MiniMed Infusion Sets in your inventory, please visit [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look) to check if you have any affected lot inventory. If you do, please cut out the LOT number labels of the affected boxes and discard recalled infusion sets and contact your Medtronic Representative who will discuss replacement stock and collate your LOT numbers labels.

The Competent Authority of your country has been notified of this action.

Medtronic considers patient safety and customer satisfaction our primary priorities. We appreciate your time and attention in reading this important notification.

Sincerely,



Samantha Baxter  
Regulatory Affairs Manager UK and Ireland

**Enclosures:**

1. Pump user letter
2. HCP letter V2
3. Key Steps: Infusion Set Priming / Fill-tubing Process

Note: Medtronic Ireland will also contact affected patients that we hold contact details for directly.



## Urgent Field Safety Notice

### Medtronic MiniMed Infusion Sets - Recall of specific Lot Numbers Potential over-delivery of insulin

September 2017

Medtronic reference: FA784

Dear Physician, Healthcare Professional,

The purpose of this letter is to notify you that we are voluntarily recalling specific lots of infusion sets used with Medtronic insulin pumps.

Medtronic has become aware of recent reports of potential over-delivery of insulin shortly after an infusion set change that have required medical intervention. Our investigation has shown this can be caused by fluid blocking the infusion set membrane during the priming/fill-tubing process. A membrane blocked by fluid most likely occurs if insulin, alcohol, or water is spilled on the top of the insulin reservoir which then could prevent the infusion set from working properly. Infusion sets currently being shipped by Medtronic contain an enhanced membrane material that significantly reduces this risk.

Because Medtronic does not have your patients' records on file\*, we ask you to inform your patients of this recall, using the attached letter, to ask them to do the following:

- A. Patients are instructed to go to [www.mmc.medtronic-diabetes.com/link](http://www.mmc.medtronic-diabetes.com/link) and determine if they have recalled infusion sets. The website will prompt them to enter the REF and LOT numbers for all infusion sets in their possession. The website will then tell them which infusion sets are part of this recall and which are not.

The REF and LOT numbers are listed on the labels as shown in the below examples:



- B. Medtronic recommends that the recalled infusion sets are not used by patients.
  - a. If they have new and enhanced infusion sets that are not part of this recall, they should use only those sets starting with their next set change. As a reminder, we have enclosed Key Steps regarding the priming/fill-tubing process.
  - b. If they only have recalled infusion sets right now, it is very important that they carefully follow the Key Steps.
- C. Patients are instructed to discard their recalled infusion. Medtronic will replace the recalled infusion sets free of charge.

# Medtronic

In order to check MiniMed Infusion Sets in your hospital inventory, please visit [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look) to check if you have any affected lot inventory, if you do please discard recalled infusion sets and contact your Medtronic Representative who will order you replacement stock free of charge.

The Competent Authority of your country has been notified of this action.

Medtronic considers patient safety and customer satisfaction our primary priorities. We appreciate your time and attention in reading this important notification.

Sincerely,



Samantha Baxter  
Regulatory Affairs Manager UK and Ireland

**\*Note: Medtronic Ireland will also contact affected patients that we hold contact details for directly.**

**Enclosures:**

1. Pump User Letter
2. Key Steps: Infusion Set Priming / Fill-tubing Process



## Medtronic MiniMed Infusion Sets Potential over-delivery of insulin

September 2017

Medtronic reference: FA784

Dear pump user,

Because the safety of our customers is our top priority, we are voluntarily recalling specific lots of infusion sets used with Medtronic insulin pumps.

### Explanation of the Issue

Medtronic has become aware of recent reports of potential over-delivery of insulin shortly after an infusion set change. Over-delivery of insulin can cause hypoglycemia and in extreme cases, death. Medtronic has received reports of hypoglycemia requiring medical intervention potentially related to this issue.

Our investigation has shown this can be caused by fluid blocking the infusion set membrane during the priming/fill-tubing process. A membrane blocked by fluid most likely occurs if insulin, alcohol, or water is spilled on top of the insulin reservoir which then could prevent the infusion set from working properly. Infusion sets currently being shipped by Medtronic contain a new and enhanced membrane material that significantly reduces this risk.

### Actions Required by You

- A. Go to [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look) determine if you have recalled infusion sets. You will be prompted to enter the Ref and Lot number for all infusion set boxes in your possession. The website will then tell you which infusion sets are part of this recall and which are not.

Your REF and LOT numbers are listed on the labels as shown in the below examples:



- B. Medtronic recommends you not use recalled infusion sets.
  - If you have new and enhanced infusion sets that are not part of this recall, use only those new and improved sets starting with your next set change.

# Medtronic

- If you only have recalled infusion sets, **it is very important to carefully follow the instructions for use** regarding the priming/fill-tubing process. You will find key steps enclosed.


C. Discard recalled infusion sets when you have new and improved infusion sets and follow instructions on the website [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look). Medtronic will replace the recalled infusion sets free of charge.

## **What if I have more questions?**

Follow the process on the website at [www.mmc.medtronic-diabetes.com/look](http://www.mmc.medtronic-diabetes.com/look). If you have additional questions call Medtronic at 01511 1444.

Medtronic considers patient safety and customer satisfaction our top priorities. We appreciate your time and attention in reading this important notification.

Sincerely,



Samantha Baxter  
Regulatory Affairs Manager UK and Ireland

## **Enclosure:**

Key Steps: Infusion Set Priming / Fill-tubing Process


## KEY STEPS

### Infusion Set Priming / Fill-tubing Process

To avoid potential over-delivery of insulin shortly after an infusion set change


It is very important to follow the steps below to prevent fluid from getting on the infusion set membrane during the priming/fill-tubing process:

- a) After filling the reservoir, make sure the vial of insulin is held upright when removing the reservoir from the blue transfer guard. This prevents insulin from accidentally getting on the top of the reservoir, which could be transferred onto the infusion set membrane.



*Hold insulin vial upright when removing reservoir.*

- b) If any liquid (such as insulin, isopropyl alcohol, or water) gets on the top of the reservoir or inside the tubing connector, it can block the infusion set membrane. This may lead to increased pressure inside the pump's reservoir chamber during the priming/fill-tubing process. This may potentially lead to over-delivery of insulin shortly after the infusion set change. If liquid is visible on top of the reservoir or inside the tubing connector, start over with a new reservoir and infusion set.



*Make sure these are dry when connecting.*

**If you notice anything unusual after the priming/fill-tubing process, such as insulin continuing to drip or squirt from the infusion set cannula, do not insert. Start over with a new reservoir and infusion set.**



**Best Practices for Changing Your Infusion Set Include:**

- Do not change your infusion set before going to sleep so that you are able to monitor your glucose levels.
- As an extra precaution, check your blood glucose at 1 hour after your infusion set change in addition to your routine monitoring.

Following these best practices will allow you to identify potential hypoglycemia and hyperglycemia so you can take necessary action.



### FA784 Affected Lots Distributed in Ireland

Product Number	Lot Serial Number	Expiry Date
MMT-368	5018235	15/01/18
MMT-368	5046554	15/11/18
MMT-368	5135966	15/01/21
MMT-368	5119768	15/08/20
MMT-368	5131873	15/11/20
MMT-368	5156825	15/07/21
MMT-368	5161553	01/09/21
MMT-399	5154358	01/09/19
MMT-368	5004175	15/08/17
MMT-368	5009099	15/10/17
MMT-368	5011831	15/10/17
MMT-368	5012330	15/10/17
MMT-368	5002988	15/10/17
MMT-368	5019152	15/01/18
MMT-368	5022898	15/02/18
MMT-368	5028303	15/03/18
MMT-368	5030861	15/04/18
MMT-368	5031734	15/04/18
MMT-368	5029838	15/04/18
MMT-368	5037923	15/07/18
MMT-368	5051351	15/10/18
MMT-368	5056651	15/12/18
MMT-368	5059299	15/01/19
MMT-368	5064182	15/03/19
MMT-368	5064849	15/03/19
MMT-368	5064848	15/03/19
MMT-368	5066434	15/06/19
MMT-368	5078887	15/08/19
MMT-368	5075369	15/08/19
MMT-368	5084142	15/09/19
MMT-368	5083714	15/01/20
MMT-368	5100278	15/02/20

MMT-368	5103796	15/03/20
MMT-368	5098611	15/04/20
MMT-368	5102050	15/05/20
MMT-368	5108594	15/05/20
MMT-368	5113196	15/06/20
MMT-368	5108027	15/07/20
MMT-368	5129181	15/10/20
MMT-368	5129182	15/11/20
MMT-368	5144273	15/03/21
MMT-368	5147467	15/04/21
MMT-368	5138510	15/04/21
MMT-368	5152671	15/07/21
MMT-368	5160444	01/08/21
MMT-368	5169287	01/10/21
MMT-368	5168199	01/10/21
MMT-368	5167003	01/10/21
MMT-368	5173454	01/12/21
MMT-368	5175616	01/01/22
MMT-368	5067603	15/05/19
MMT-377	5011832	15/10/17
MMT-377	5022226	15/01/18
MMT-377	5026475	15/03/18
MMT-377	5026444	15/03/18
MMT-377	5030967	15/04/18
MMT-377	5034236	15/06/18
MMT-377	5042576	15/08/18
MMT-377	5043393	15/08/18
MMT-377	5051354	15/11/18
MMT-377	5047491	15/01/19
MMT-377	5066395	15/03/19
MMT-377	5070833	15/05/19
MMT-377	5084755	15/09/19
MMT-377	5083719	15/01/20
MMT-377	5101576	15/02/20
MMT-377	5098612	15/03/20
MMT-377	5103798	15/04/20
MMT-377	5106197	15/04/20
MMT-377	5132536	15/11/20
MMT-377	5134841	15/12/20
MMT-377	5137380	15/02/21
MMT-377	5140283	15/02/21
MMT-377	5147170	15/04/21

MMT-377	5156826	15/07/21
MMT-377	5160446	01/08/21
MMT-377	5168201	01/10/21
MMT-377	5175025	01/01/22
MMT-378	5017473	15/01/18
MMT-378	5004197	15/09/17
MMT-378	5004711	15/09/17
MMT-378	5011141	15/10/17
MMT-378	5013129	15/11/17
MMT-378	5017087	15/12/17
MMT-378	5022901	15/01/18
MMT-378	5017472	15/01/18
MMT-378	5017486	15/03/18
MMT-378	5025151	15/03/18
MMT-378	5032760	15/05/18
MMT-378	5032670	15/06/18
MMT-378	5041847	15/08/18
MMT-378	5041856	15/08/18
MMT-378	5043399	15/09/18
MMT-378	5048409	15/09/18
MMT-378	5059302	15/01/19
MMT-378	5059301	15/01/19
MMT-378	5063531	15/02/19
MMT-378	5063530	15/02/19
MMT-378	5063532	15/03/19
MMT-378	5065534	15/03/19
MMT-378	5066633	15/04/19
MMT-378	5072352	15/06/19
MMT-378	5084756	15/09/19
MMT-378	5089245	15/11/19
MMT-378	5088660	15/11/19
MMT-378	5088659	15/11/19
MMT-378	5095742	15/01/20
MMT-378	5095744	15/01/20
MMT-378	5099024	15/02/20
MMT-378	5097152	15/03/20
MMT-378	5106198	15/04/20
MMT-378	5103801	15/04/20
MMT-378	5102250	15/05/20
MMT-378	5109101	15/05/20
MMT-378	5108599	15/05/20
MMT-378	5113965	15/06/20

MMT-378	5113966	15/07/20
MMT-378	5117518	15/07/20
MMT-378	5119910	15/09/20
MMT-378	5119911	15/09/20
MMT-378	5122378	15/09/20
MMT-378	5129185	15/10/20
MMT-378	5117805	15/10/20
MMT-378	5131876	15/11/20
MMT-378	5134844	15/12/20
MMT-378	5140427	15/02/21
MMT-378	5134443	15/03/21
MMT-378	5147172	15/04/21
MMT-378	5148148	15/04/21
MMT-378	5144320	15/04/21
MMT-378	5149823	15/04/21
MMT-378	5156556	15/06/21
MMT-378	5152673	15/06/21
MMT-378	5148939	15/07/21
MMT-378	5165249	01/10/21
MMT-378	5168518	01/11/21
MMT-378	5168200	01/11/21
MMT-378	5174169	01/12/21
MMT-378	5175023	01/01/22
MMT-378	5177468	01/01/22
MMT-378	5007690	15/09/17
MMT-381	5035709	15/07/18
MMT-381	5038849	15/07/18
MMT-381	5042562	15/09/18
MMT-381	5066397	15/04/19
MMT-381	5150925	15/07/21
MMT-381	5119246	15/08/20
MMT-381	5032657	15/05/18
MMT-381	5099025	15/02/20
MMT-381	5155453	15/06/21
MMT-381	5109405	15/06/20
MMT-381	5022228	15/01/18
MMT-381	5064857	15/03/19
MMT-381	5018243	15/01/18
MMT-381	5027641	15/04/18
MMT-381	5065536	15/03/19
MMT-381	5170949	01/12/21
MMT-381	5037925	15/06/18

MMT-381	5055906	15/02/19
MMT-381	5034237	15/06/18
MMT-381	5012918	15/11/17
MMT-381	5015509	15/12/17
MMT-381	5033425	15/05/18
MMT-381	5059940	15/01/19
MMT-381	5067596	15/04/19
MMT-381	5067488	15/06/19
MMT-381	5148813	15/04/21
MMT-381	5115661	15/07/20
MMT-1001	0212069586	15/07/19
MMT-381	5117102	15/07/20
MMT-381	5134847	15/01/21
MMT-381	5161551	01/09/21
MMT-381	5009103	15/09/17
MMT-381	5003718	15/09/17
MMT-381	5003719	15/09/17
MMT-381	5010306	15/09/17
MMT-381	5007685	15/09/17
MMT-381	5012300	15/10/17
MMT-381	5011171	15/10/17
MMT-381	5012917	15/11/17
MMT-381	5016921	15/12/17
MMT-381	5023618	15/02/18
MMT-381	5022908	15/02/18
MMT-381	5022906	15/02/18
MMT-381	5023616	15/02/18
MMT-381	5026447	15/03/18
MMT-381	5028262	15/03/18
MMT-381	5028263	15/04/18
MMT-381	5037927	15/06/18
MMT-381	5037924	15/06/18
MMT-381	5051363	15/11/18
MMT-381	5053728	15/12/18
MMT-381	5064856	15/03/19
MMT-381	5066396	15/04/19
MMT-381	5066645	15/05/19
MMT-381	5070169	15/05/19
MMT-381	5079923	15/08/19
MMT-381	5083739	15/10/19
MMT-381	5088662	15/10/19
MMT-381	5088666	15/10/19

MMT-381	5089250	15/12/19
MMT-381	5095745	15/01/20
MMT-381	5095746	15/01/20
MMT-381	5179978	01/02/20
MMT-381	5179858	01/02/20
MMT-381	5093624	15/02/20
MMT-381	5098613	15/03/20
MMT-381	5103363	15/03/20
MMT-381	5100604	15/04/20
MMT-381	5105374	15/04/20
MMT-381	5109107	15/05/20
MMT-381	5109103	15/05/20
MMT-381	5109106	15/05/20
MMT-381	5109105	15/05/20
MMT-381	5113201	15/06/20
MMT-381	5113968	15/07/20
MMT-381	5119642	15/09/20
MMT-381	5119643	15/09/20
MMT-381	5122169	15/09/20
MMT-381	5117798	15/09/20
MMT-381	5119651	15/10/20
MMT-381	5129192	15/11/20
MMT-381	5132941	15/12/20
MMT-381	5136402	15/01/21
MMT-381	5136400	15/01/21
MMT-381	5130633	15/01/21
MMT-381	5137381	15/02/21
MMT-381	5140285	15/02/21
MMT-381	5132811	15/02/21
MMT-381	5144322	15/03/21
MMT-381	5140433	15/03/21
MMT-381	5144279	15/03/21
MMT-381	5139060	15/03/21
MMT-381	5141806	15/04/21
MMT-381	5147470	15/04/21
MMT-381	5149826	15/04/21
MMT-381	5148814	15/04/21
MMT-381	5148950	15/06/21
MMT-381	5148951	15/06/21
MMT-381	5158224	15/08/21
MMT-381	5163053	01/09/21
MMT-381	5161569	01/09/21

MMT-381	5169288	01/10/21
MMT-381	5164403	01/10/21
MMT-381	5171144	01/11/21
MMT-381	5170676	01/11/21
MMT-381	5170950	01/12/21
MMT-381	5177843	01/02/22
MMT-381	5034320	15/06/18
MMT-382	5004185	15/08/17
MMT-382	5022229	15/01/18
MMT-382	5022909	15/02/18
MMT-382	5030870	15/04/18
MMT-382	5030978	15/04/18
MMT-382	5034322	15/05/18
MMT-382	5052917	15/11/18
MMT-382	5063534	15/03/19
MMT-382	5063533	15/03/19
MMT-382	5063629	15/04/19
MMT-382	5067610	15/05/19
MMT-382	5073130	15/10/19
MMT-382	5088669	15/11/19
MMT-382	5093900	15/12/19
MMT-382	5095747	15/01/20
MMT-382	5093607	15/03/20
MMT-382	5102054	15/05/20
MMT-382	5110549	15/06/20
MMT-382	5118536	15/08/20
MMT-382	5119912	15/09/20
MMT-382	5129193	15/10/20
MMT-382	5131894	15/11/20
MMT-382	5149830	15/04/21
MMT-382	5150523	15/05/21
MMT-382	5152670	15/07/21
MMT-382	5160447	01/08/21
MMT-382	5165251	01/10/21
MMT-382	5173457	01/12/21
MMT-383	5017516	15/01/18
MMT-383	5028261	15/03/18
MMT-383	5030979	15/04/18
MMT-383	5028970	15/05/18
MMT-383	5037929	15/07/18
MMT-383	5048111	15/09/18
MMT-383	5050061	15/10/18

MMT-383	5046570	15/11/18
MMT-383	5056656	15/12/18
MMT-383	5059306	15/01/19
MMT-383	5057790	15/01/19
MMT-383	5059943	15/02/19
MMT-383	5066652	15/04/19
MMT-383	5067385	15/06/19
MMT-383	5079924	15/08/19
MMT-383	5084149	15/09/19
MMT-383	5103366	15/03/20
MMT-383	5105376	15/04/20
MMT-383	5115662	15/07/20
MMT-383	5119290	15/08/20
MMT-383	5129194	15/10/20
MMT-383	5132539	15/11/20
MMT-383	5135980	15/01/21
MMT-383	5132392	15/02/21
MMT-383	5157554	15/07/21
MMT-383	5152669	15/07/21
MMT-383	5168202	01/10/21
MMT-383	5033532	15/07/18
MMT-384	5010137	15/11/17
MMT-384	5015518	15/12/17
MMT-384	5038025	15/08/18
MMT-384	5043411	15/08/18
MMT-386	5085615	15/01/18
MMT-386	5108135	15/08/18
MMT-386	5130913	15/12/18
MMT-386	5139081	15/04/19
MMT-386	5141579	15/05/19
MMT-386	5164682	01/10/19
MMT-386	5164685	01/11/19
MMT-386	5164686	01/11/19
MMT-386	5173833	01/12/19
MMT-386	5179666	01/02/20
MMT-387	5109316	15/09/18
MMT-387	5131828	15/12/18
MMT-387	5144452	15/06/19
MMT-387	5164595	01/01/20
MMT-387	5131823	15/11/18
MMT-387	5144451	15/05/19
MMT-387	5144453	15/06/19



MMT-387	5153091	15/07/19
MMT-387	5153095	01/09/19
MMT-387	5164592	01/11/19
MMT-387	5164593	01/12/19
MMT-387	5176046	01/01/20
MMT-387	5176048	01/02/20
MMT-394	5122571	15/10/18
MMT-394	5079863	15/08/17
MMT-394	5079865	15/08/17
MMT-394	5079868	15/10/17
MMT-394	5079871	15/12/17
MMT-394	5094272	15/01/18
MMT-394	5094273	15/02/18
MMT-394	5094274	15/03/18
MMT-394	5094276	15/04/18
MMT-394	5094278	15/04/18
MMT-394	5109328	15/05/18
MMT-394	5109330	15/06/18
MMT-394	5109334	15/07/18
MMT-394	5109336	15/08/18
MMT-394	5122185	15/09/18
MMT-394	5122186	15/09/18
MMT-394	5122557	15/09/18
MMT-394	5122569	15/10/18
MMT-394	5122572	15/11/18
MMT-394	5132385	15/12/18
MMT-394	5139112	15/02/19
MMT-394	5139115	15/03/19
MMT-394	5144055	15/03/19
MMT-394	5144500	15/05/19
MMT-394	5144502	15/06/19
MMT-394	5154158	15/07/19
MMT-394	5144505	15/07/19
MMT-394	5154160	01/08/19
MMT-394	5154161	01/08/19
MMT-394	5154165	01/09/19
MMT-394	5164693	01/10/19
MMT-394	5164697	01/11/19
MMT-394	5164694	01/11/19
MMT-394	5164700	01/12/19
MMT-394	5164699	01/12/19
MMT-394	5174475	01/01/20

MMT-394	5174479	01/02/20
MMT-394	5109332	15/06/18
MMT-396	5079874	15/09/17
MMT-396	5086244	15/12/17
MMT-396	5093716	15/03/18
MMT-396	5105644	15/06/18
MMT-396	5122190	15/09/18
MMT-396	5144487	15/04/19
MMT-396	5144541	01/09/19
MMT-396	5139137	15/02/19
MMT-396	5078044	15/08/17
MMT-396	5079873	15/09/17
MMT-396	5086240	15/09/17
MMT-396	5093713	15/01/18
MMT-396	5093715	15/02/18
MMT-396	5093718	15/03/18
MMT-396	5093717	15/03/18
MMT-396	5105638	15/04/18
MMT-396	5105637	15/04/18
MMT-396	5105640	15/05/18
MMT-396	5105641	15/05/18
MMT-396	5105642	15/05/18
MMT-396	5108161	15/08/18
MMT-396	5108162	15/08/18
MMT-396	5122578	15/10/18
MMT-396	5137416	15/01/19
MMT-396	5139139	15/03/19
MMT-396	5139138	15/03/19
MMT-396	5144489	15/04/19
MMT-396	5144490	15/04/19
MMT-396	5144491	15/04/19
MMT-396	5144488	15/04/19
MMT-396	5144492	15/04/19
MMT-396	5144495	15/07/19
MMT-396	5144496	15/08/19
MMT-396	5144539	01/09/19
MMT-396	5144544	01/10/19
MMT-396	5144542	01/10/19
MMT-396	5171604	01/11/19
MMT-396	5144545	01/11/19
MMT-396	5171605	01/11/19
MMT-396	5171607	01/12/19

MMT-396	5171606	01/12/19
MMT-396	5173458	01/12/19
MMT-396	5171608	01/12/19
MMT-396	5173459	01/01/20
MMT-397	5093680	15/01/18
MMT-397	5101211	15/04/18
MMT-397	5122054	15/10/18
MMT-397	5100666	15/04/18
MMT-397	5108255	15/06/18
MMT-397	5105580	15/05/18
MMT-397	5131595	15/12/18
MMT-397	5143480	15/05/19
MMT-397	5122042	15/10/18
MMT-397	5131582	15/12/18
MMT-397	5143484	15/05/19
MMT-397	5093679	15/01/18
MMT-397	5108254	15/06/18
MMT-397	5165416	01/10/19
MMT-397	5097169	15/03/18
MMT-397	5105581	15/05/18
MMT-397	5105569	15/05/18
MMT-397	5108267	15/06/18
MMT-397	5108259	15/06/18
MMT-397	5114567	15/08/18
MMT-397	5114570	15/08/18
MMT-397	5139929	15/03/19
MMT-397	5143469	15/04/19
MMT-397	5143498	15/06/19
MMT-397	5154333	01/08/19
MMT-923	5161461	01/08/19
MMT-397	5176002	01/02/20
MMT-397	5078273	15/08/17
MMT-397	5078283	15/08/17
MMT-397	5078284	15/08/17
MMT-397	5088111	15/10/17
MMT-397	5088115	15/10/17
MMT-397	5078299	15/10/17
MMT-397	5088657	15/12/17
MMT-397	5093708	15/02/18
MMT-397	5097160	15/02/18
MMT-397	5100659	15/03/18
MMT-397	5097171	15/03/18

MMT-397	5100661	15/03/18
MMT-397	5100660	15/03/18
MMT-397	5101209	15/04/18
MMT-397	5105567	15/05/18
MMT-397	5108258	15/06/18
MMT-397	5108244	15/06/18
MMT-397	5108245	15/06/18
MMT-397	5108281	15/08/18
MMT-397	5114568	15/08/18
MMT-397	5108279	15/08/18
MMT-397	5120649	15/09/18
MMT-397	5120653	15/09/18
MMT-397	5122035	15/09/18
MMT-397	5120651	15/09/18
MMT-397	5122180	15/09/18
MMT-397	5122055	15/10/18
MMT-397	5122056	15/10/18
MMT-397	5131591	15/12/18
MMT-397	5137492	15/01/19
MMT-397	5138221	15/02/19
MMT-397	5139842	15/02/19
MMT-397	5139841	15/02/19
MMT-397	5143448	15/03/19
MMT-397	5143459	15/03/19
MMT-397	5139932	15/03/19
MMT-397	5143453	15/03/19
MMT-397	5143449	15/03/19
MMT-397	5143457	15/03/19
MMT-397	5143462	15/04/19
MMT-397	5143482	15/05/19
MMT-397	5143491	15/05/19
MMT-397	5143486	15/05/19
MMT-397	5143494	15/05/19
MMT-397	5143481	15/05/19
MMT-397	5143502	15/06/19
MMT-397	5143500	15/06/19
MMT-397	5151295	15/07/19
MMT-397	5152595	15/07/19
MMT-397	5152597	15/07/19
MMT-397	5154332	01/08/19
MMT-397	5154329	01/08/19
MMT-397	5154335	01/08/19

MMT-397	5162856	01/09/19
MMT-397	5162860	01/09/19
MMT-397	5165415	01/09/19
MMT-397	5154340	01/09/19
MMT-397	5165413	01/09/19
MMT-397	5165429	01/10/19
MMT-397	5165417	01/10/19
MMT-397	5165442	01/11/19
MMT-397	5165439	01/11/19
MMT-397	5165437	01/11/19
MMT-397	5173510	01/12/19
MMT-397	5173509	01/12/19
MMT-397	5171649	01/12/19
MMT-397	5171655	01/12/19
MMT-397	5173523	01/01/20
MMT-397	5173518	01/01/20
MMT-397	5175996	01/01/20
MMT-397	5176011	01/02/20
MMT-397	5176003	01/02/20
MMT-398	5123562	15/10/18
MMT-398	5131253	15/03/19
MMT-398	5076583	15/08/17
MMT-398	5079887	15/09/17
MMT-398	5079885	15/09/17
MMT-398	5079892	15/11/17
MMT-398	5094989	15/01/18
MMT-398	5103356	15/05/18
MMT-398	5108320	15/05/18
MMT-398	5108325	15/06/18
MMT-398	5108327	15/07/18
MMT-398	5123561	15/10/18
MMT-398	5109518	15/10/18
MMT-398	5108336	15/10/18
MMT-398	5123565	15/12/18
MMT-398	5131247	15/02/19
MMT-398	5131250	15/02/19
MMT-398	5144206	15/03/19
MMT-398	5131252	15/03/19
MMT-398	5144208	15/04/19
MMT-398	5144211	15/05/19
MMT-398	5144214	15/06/19
MMT-398	5146686	15/06/19

MMT-398	5146690	01/08/19
MMT-398	5164604	01/10/19
MMT-398	5164612	01/12/19
MMT-398	5164615	01/01/20
MMT-398	5164618	01/02/20
MMT-399	5114613	15/09/18
MMT-399	5139145	15/01/19
MMT-399	5139944	15/03/19
MMT-399	5164664	01/11/19
MMT-399	5131675	15/01/19
MMT-399	5146714	15/06/19
MMT-399	5123591	15/11/18
MMT-399	5144226	15/04/19
MMT-399	5079032	15/08/17
MMT-399	5101289	15/04/18
MMT-399	5154356	01/09/19
MMT-399	5087095	15/10/17
MMT-399	5146713	15/05/19
MMT-399	5154360	01/09/19
MMT-399	5170863	01/12/19
MMT-399	5164645	01/09/19
MMT-399	5101291	15/04/18
MMT-399	5123600	15/11/18
MMT-399	5131660	15/01/19
MMT-399	5131662	15/01/19
MMT-399	5139854	15/03/19
MMT-399	5144216	15/03/19
MMT-399	5144215	15/03/19
MMT-399	5144218	15/04/19
MMT-399	5144223	15/04/19
MMT-399	5146707	15/05/19
MMT-399	5146709	15/05/19
MMT-399	5146723	15/06/19
MMT-399	5154350	01/08/19
MMT-399	5154354	01/09/19
MMT-1001	0212104481	15/07/19
MMT-862	5159805	01/08/21
MMT-399	5170872	01/12/19
MMT-1001	0211097217	15/01/19
MMT-1001	0212104646	15/07/19
MMT-399	5146716	15/06/19
MMT-399	5079045	15/09/17

MMT-399	5087099	15/10/17
MMT-399	5079051	15/10/17
MMT-399	5079049	15/10/17
MMT-399	5079050	15/10/17
MMT-399	5088367	15/11/17
MMT-399	5088366	15/11/17
MMT-399	5088368	15/12/17
MMT-399	5088370	15/12/17
MMT-399	5088373	15/01/18
MMT-399	5088371	15/01/18
MMT-399	5098581	15/03/18
MMT-399	5098579	15/03/18
MMT-399	5098580	15/03/18
MMT-399	5105505	15/04/18
MMT-399	5101292	15/04/18
MMT-399	5105509	15/05/18
MMT-399	5105518	15/06/18
MMT-399	5105517	15/06/18
MMT-399	5105521	15/06/18
MMT-399	5108378	15/06/18
MMT-399	5108387	15/07/18
MMT-399	5108394	15/08/18
MMT-399	5114614	15/09/18
MMT-399	5114622	15/09/18
MMT-399	5114611	15/09/18
MMT-399	5123575	15/10/18
MMT-399	5123568	15/10/18
MMT-399	5123599	15/11/18
MMT-399	5123597	15/11/18
MMT-399	5123576	15/11/18
MMT-399	5131654	15/12/18
MMT-399	5131673	15/01/19
MMT-399	5131668	15/01/19
MMT-399	5139851	15/02/19
MMT-399	5139853	15/03/19
MMT-399	5139938	15/03/19
MMT-399	5139945	15/03/19
MMT-399	5139941	15/03/19
MMT-399	5144224	15/04/19
MMT-399	5144220	15/04/19
MMT-399	5146697	15/04/19
MMT-399	5146706	15/05/19

MMT-399	5146717	15/06/19
MMT-399	5146718	15/06/19
MMT-399	5146727	15/06/19
MMT-399	5146732	15/07/19
MMT-399	5146735	15/07/19
MMT-399	5154348	01/08/19
MMT-399	5154349	01/08/19
MMT-399	5154357	01/09/19
MMT-399	5154352	01/09/19
MMT-399	5164654	01/10/19
MMT-399	5164648	01/10/19
MMT-399	5164656	01/10/19
MMT-399	5164657	01/10/19
MMT-399	5164658	01/10/19
MMT-399	5164653	01/10/19
MMT-399	5164649	01/10/19
MMT-399	5164650	01/10/19
MMT-399	5164660	01/10/19
MMT-399	5164666	01/11/19
MMT-399	5164665	01/11/19
MMT-399	5164663	01/11/19
MMT-399	5164668	01/11/19
MMT-399	5170861	01/12/19
MMT-399	5170860	01/12/19
MMT-399	5170859	01/12/19
MMT-399	5170969	01/01/20
MMT-399	5170881	01/01/20
MMT-399	5170970	01/01/20
MMT-399	5170977	01/02/20
MMT-399	5105520	15/06/18
MMT-399	5131667	15/01/19
MMT-1001	0211117021	15/01/19
MMT-399	5131658	15/01/19
MMT-399	5093759	15/02/18
MMT-399	5101290	15/04/18
MMT-399	5123592	15/11/18
MMT-399	5131665	15/01/19
MMT-399	5146738	15/07/19
MMT-399	5114618	15/09/18
MMT-862	5038783	15/08/18
MMT-862	5118537	15/08/20
MMT-862	5043415	15/09/18



MMT-862	5072333	15/07/19
MMT-862	5127927	15/10/20
MMT-862	5056680	15/12/18
MMT-862	5047168	15/12/18
MMT-862	5119920	15/08/20
MMT-862	5046638	15/09/18
MMT-862	5119226	15/08/20
MMT-862	5063694	15/02/19
MMT-862	5087581	15/10/19
MMT-862	5117946	15/07/20
MMT-862	5149834	15/05/21
MMT-862	5089775	15/11/19
MMT-862	5029016	15/04/18
MMT-862	5030895	15/04/18
MMT-862	5040157	15/07/18
MMT-862	5043414	15/08/18
MMT-862	5046639	15/10/18
MMT-862	5052903	15/11/18
MMT-862	5057953	15/01/19
MMT-862	5068516	15/05/19
MMT-862	5073682	15/06/19
MMT-862	5085471	15/09/19
MMT-862	5086259	15/10/19
MMT-862	5088971	15/11/19
MMT-862	5095083	15/12/19
MMT-862	5103784	15/03/20
MMT-862	5101637	15/03/20
MMT-862	5087768	15/04/20
MMT-862	5103785	15/04/20
MMT-862	5105356	15/05/20
MMT-862	5107545	15/05/20
MMT-862	5109337	15/07/20
MMT-862	5123218	15/09/20
MMT-862	5129156	15/11/20
MMT-862	5131848	15/11/20
MMT-862	5134865	15/12/20
MMT-862	5138176	15/01/21
MMT-862	5137408	15/01/21
MMT-862	5140266	15/02/21
MMT-862	5134866	15/03/21
MMT-862	5143881	15/03/21
MMT-862	5143880	15/03/21

MMT-862	5147532	15/04/21
MMT-862	5150583	15/05/21
MMT-862	5156260	15/06/21
MMT-862	5152957	15/06/21
MMT-862	5150565	15/07/21
MMT-862	5161076	01/08/21
MMT-862	5165191	01/09/21
MMT-862	5161874	01/09/21
MMT-862	5165190	01/09/21
MMT-862	5169314	01/11/21
MMT-862	5165395	01/12/21
MMT-862	5175703	01/12/21
MMT-862	5176806	01/01/22
MMT-862	5178110	01/01/22
MMT-862	5180074	01/02/22
MMT-862	5140377	15/02/21
MMT-862	5033394	15/06/18
MMT-864	5118411	15/07/20
MMT-864	5065089	15/03/19
MMT-864	5047171	15/11/18
MMT-864	5118185	15/08/20
MMT-864	5108397	15/05/20
MMT-864	5149836	15/04/21
MMT-864	5034327	15/06/18
MMT-864	5106839	15/05/20
MMT-864	5118295	15/08/20
MMT-864	5139031	15/02/21
MMT-864	5153932	15/06/21
MMT-864	5004572	15/11/17
MMT-864	5007429	15/11/17
MMT-864	5016922	15/12/17
MMT-864	5020537	15/01/18
MMT-864	5022955	15/02/18
MMT-864	5023463	15/02/18
MMT-864	5021556	15/02/18
MMT-864	5023465	15/03/18
MMT-864	5026739	15/03/18
MMT-864	5030896	15/04/18
MMT-864	5028965	15/04/18
MMT-864	5030899	15/05/18
MMT-864	5035649	15/07/18
MMT-864	5036356	15/07/18

MMT-864	5035690	15/08/18
MMT-864	5040466	15/08/18
MMT-864	5046395	15/09/18
MMT-864	5047582	15/09/18
MMT-864	5050035	15/10/18
MMT-864	5043106	15/11/18
MMT-864	5054558	15/12/18
MMT-864	5060948	15/02/19
MMT-864	5057231	15/03/19
MMT-864	5064142	15/03/19
MMT-864	5067589	15/05/19
MMT-864	5070664	15/05/19
MMT-864	5072716	15/06/19
MMT-864	5079683	15/07/19
MMT-864	5077462	15/07/19
MMT-864	5077107	15/09/19
MMT-864	5086971	15/10/19
MMT-864	5085438	15/10/19
MMT-864	5089778	15/11/19
MMT-864	5093164	15/12/19
MMT-864	5083273	15/12/19
MMT-864	5096974	15/01/20
MMT-864	5097021	15/01/20
MMT-864	5104623	15/04/20
MMT-864	5109089	15/05/20
MMT-864	5114171	15/06/20
MMT-864	5113369	15/06/20
MMT-864	5117947	15/07/20
MMT-864	5119915	15/08/20
MMT-864	5118296	15/08/20
MMT-864	5123292	15/09/20
MMT-864	5123295	15/10/20
MMT-864	5129148	15/11/20
MMT-864	5131850	15/11/20
MMT-864	5134869	15/12/20
MMT-864	5134555	15/12/20
MMT-864	5137410	15/01/21
MMT-864	5140268	15/02/21
MMT-864	5143886	15/03/21
MMT-864	5143885	15/03/21
MMT-864	5147533	15/04/21
MMT-864	5147559	15/04/21

MMT-864	5148883	15/04/21
MMT-864	5151824	15/05/21
MMT-864	5150585	15/05/21
MMT-864	5152314	15/06/21
MMT-864	5155456	15/06/21
MMT-864	5158839	01/08/21
MMT-864	5159784	01/08/21
MMT-864	5162136	01/08/21
MMT-864	5166005	01/09/21
MMT-864	5167800	01/10/21
MMT-864	5165398	01/10/21
MMT-864	5167794	01/10/21
MMT-864	5167799	01/10/21
MMT-864	5172573	01/11/21
MMT-864	5172577	01/12/21
MMT-864	5177297	01/01/22
MMT-864	5177299	01/01/22
MMT-864	5180091	01/02/22
MMT-864	5180066	01/02/22
MMT-866	5167125	01/10/21
MMT-866	5003932	15/08/17
MMT-866	5017733	15/12/17
MMT-866	5017653	15/02/18
MMT-866	5026883	15/03/18
MMT-866	5043420	15/08/18
MMT-866	5046397	15/09/18
MMT-866	5047583	15/10/18
MMT-866	5050037	15/10/18
MMT-866	5047175	15/12/18
MMT-866	5059551	15/01/19
MMT-866	5061497	15/02/19
MMT-866	5066515	15/04/19
MMT-866	5059523	15/05/19
MMT-866	5067592	15/05/19
MMT-866	5085473	15/09/19
MMT-866	5085432	15/09/19
MMT-866	5086972	15/10/19
MMT-866	5088975	15/10/19
MMT-866	5090814	15/01/20
MMT-866	5096975	15/01/20
MMT-866	5104628	15/04/20
MMT-866	5109094	15/05/20

MMT-866	5113371	15/06/20
MMT-866	5109527	15/06/20
MMT-866	5115466	15/07/20
MMT-866	5114173	15/07/20
MMT-866	5117180	15/07/20
MMT-866	5118412	15/07/20
MMT-866	5118646	15/08/20
MMT-866	5118648	15/09/20
MMT-866	5129144	15/10/20
MMT-866	5131851	15/11/20
MMT-866	5133978	15/12/20
MMT-866	5134873	15/12/20
MMT-866	5138179	15/01/21
MMT-866	5137411	15/01/21
MMT-866	5140722	15/02/21
MMT-866	5139033	15/02/21
MMT-866	5143889	15/03/21
MMT-866	5143888	15/03/21
MMT-866	5148885	15/04/21
MMT-866	5150410	15/05/21
MMT-866	5150588	15/05/21
MMT-866	5150411	15/05/21
MMT-866	5152991	15/06/21
MMT-866	5155458	15/06/21
MMT-866	5158847	01/07/21
MMT-866	5158846	01/07/21
MMT-866	5152331	15/07/21
MMT-866	5162246	01/08/21
MMT-866	5159825	01/08/21
MMT-866	5158842	01/09/21
MMT-866	5161866	01/10/21
MMT-866	5175697	01/12/21
MMT-866	5180088	01/02/22
MMT-866	5181484	01/02/22
MMT-866	5184073	01/03/22
MMT-866	5184072	01/03/22
MMT-866	5079689	15/08/19
MMT-866	5020541	15/01/18
MMT-866	5019182	15/01/18
MMT-874	5040471	15/08/18
MMT-874	5002817	15/08/17
MMT-874	5002816	15/08/17

MMT-874	5009115	15/09/17
MMT-874	5017734	15/12/17
MMT-874	5020544	15/02/18
MMT-874	5040470	15/07/18
MMT-874	5055465	15/04/19
MMT-874	5072719	15/06/19
MMT-874	5077464	15/07/19
MMT-874	5081723	15/10/19
MMT-874	5089781	15/11/19
MMT-874	5093165	15/12/19
MMT-874	5098841	15/02/20
MMT-874	5105360	15/04/20
MMT-874	5105361	15/05/20
MMT-874	5100585	15/05/20
MMT-874	5113372	15/06/20
MMT-874	5114175	15/07/20
MMT-874	5120630	15/09/20
MMT-874	5128418	15/12/20
MMT-874	5135125	15/12/20
MMT-874	5140270	15/02/21
MMT-874	5146999	15/03/21
MMT-874	5143892	15/03/21
MMT-874	5151827	15/05/21
MMT-874	5156821	15/06/21
MMT-874	5152960	15/06/21
MMT-874	5150896	15/07/21
MMT-874	5162247	01/08/21
MMT-874	5161079	01/08/21
MMT-874	5166031	01/09/21
MMT-874	5165196	01/09/21
MMT-874	5167797	01/10/21
MMT-874	5174592	01/12/21
MMT-874	5180069	01/01/22
MMT-874	5177301	01/01/22
MMT-876	5003938	15/08/17
MMT-876	5047179	15/01/19
MMT-876	5059553	15/01/19
MMT-876	5093166	15/12/19
MMT-876	5095828	15/01/20
MMT-876	5107546	15/05/20
MMT-876	5109532	15/06/20
MMT-876	5109531	15/06/20

MMT-876	5108435	15/07/20
MMT-876	5117948	15/07/20
MMT-876	5119918	15/08/20
MMT-876	5123223	15/10/20
MMT-876	5138181	15/01/21
MMT-876	5143971	15/03/21
MMT-876	5148146	15/04/21
MMT-876	5143972	15/04/21
MMT-876	5149838	15/05/21
MMT-876	5151828	15/05/21
MMT-876	5153935	15/06/21
MMT-876	5161080	01/08/21
MMT-876	5163929	01/09/21
MMT-876	5167128	01/10/21
MMT-876	5175702	01/12/21
MMT-876	5165968	01/12/21
MMT-876	5176807	01/12/21
MMT-876	5178117	01/01/22
MMT-876	5178488	01/01/22
MMT-876	5186039	01/03/22
MMT-884	5039586	15/01/19
MMT-884	5085435	15/09/19
MMT-886	5105647	15/06/20
MMT-886	5068557	15/05/19
MMT-886	5154652	15/06/21
MMT-886	5156823	15/07/21
MMT-886	5150905	01/08/21
MMT-886	5163930	01/09/21
MMT-886	5180070	01/01/22
MMT-905	5119529	15/08/18
MMT-905	5119525	15/08/18
MMT-905	5114235	15/08/18
MMT-905	5114246	15/08/18
MMT-905	5114264	15/10/18
MMT-905	5148891	15/04/19
MMT-905	5156265	15/07/19
MMT-905	5119527	15/08/18
MMT-921	5085187	15/09/17
MMT-921	5097241	15/01/18
MMT-921	5103409	15/05/18
MMT-921	5150591	15/05/19
MMT-921	5084460	15/10/17

MMT-921	5091305	15/11/17
MMT-921	5095069	15/01/18
MMT-921	5119233	15/08/18
MMT-921	5131854	15/11/18
MMT-1001	0212124157	15/07/19
MMT-921	5079697	15/08/17
MMT-921	5082268	15/08/17
MMT-921	5079698	15/08/17
MMT-921	5082622	15/09/17
MMT-921	5082773	15/09/17
MMT-921	5086974	15/10/17
MMT-921	5084464	15/12/17
MMT-921	5101297	15/04/18
MMT-921	5104632	15/04/18
MMT-921	5104247	15/05/18
MMT-921	5109789	15/06/18
MMT-921	5117935	15/07/18
MMT-921	5115456	15/07/18
MMT-921	5109702	15/08/18
MMT-921	5109707	15/09/18
MMT-921	5122398	15/09/18
MMT-921	5123224	15/10/18
MMT-921	5123228	15/10/18
MMT-921	5127932	15/10/18
MMT-921	5131255	15/11/18
MMT-921	5135127	15/12/18
MMT-921	5133982	15/12/18
MMT-921	5135128	15/01/19
MMT-921	5135176	15/01/19
MMT-921	5140384	15/02/19
MMT-921	5132018	15/03/19
MMT-921	5143983	15/03/19
MMT-921	5143984	15/04/19
MMT-921	5147485	15/04/19
MMT-921	5141588	15/05/19
MMT-921	5151829	15/06/19
MMT-921	5154653	15/06/19
MMT-921	5156088	15/07/19
MMT-921	5150551	15/07/19
MMT-921	5158191	01/08/19
MMT-921	5163949	01/09/19
MMT-921	5162250	01/09/19



MMT-921	5162886	01/10/19
MMT-921	5167139	01/10/19
MMT-921	5152526	01/10/19
MMT-921	5167818	01/10/19
MMT-921	5167819	01/11/19
MMT-921	5169320	01/11/19
MMT-921	5169321	01/11/19
MMT-921	5174510	01/12/19
MMT-921	5176581	01/01/20
MMT-921	5180040	01/02/20
MMT-921	5180041	01/02/20
MMT-921	5178033	01/02/20
MMT-921	5101624	15/03/18
MMT-1001	0212287272	15/07/19
MMT-923	5122400	15/09/18
MMT-923	5096419	15/01/18
MMT-923	5122399	15/09/18
MMT-923	5178180	01/01/20
MMT-923	5115454	15/07/18
MMT-923	5085476	15/09/17
MMT-923	5087587	15/10/17
MMT-923	5120633	15/08/18
MMT-923	5084462	15/09/17
MMT-923	5128448	15/10/18
MMT-923	5139009	15/01/19
MMT-923	5151832	15/06/19
MMT-923	5175707	01/12/19
MMT-923	5131856	15/11/18
MMT-923	5165202	01/10/19
MMT-923	5088982	15/10/17
MMT-923	5109709	15/06/18
MMT-923	5134433	15/01/19
MMT-923	5158194	01/08/19
MMT-923	5162225	01/08/19
MMT-923	5174512	01/12/19
MMT-923	5167802	01/10/19
MMT-923	5163247	01/09/19
MMT-923	5119934	15/08/18
MMT-923	5084126	15/09/17
MMT-923	5095196	15/01/18
MMT-923	5109792	15/06/18
MMT-923	5109711	15/09/18

MMT-923	5107548	15/05/18
MMT-923	5151830	15/06/19
MMT-923	5169323	01/11/19
MMT-923	5140386	15/02/19
MMT-923	5154655	15/06/19
MMT-923	5147487	15/04/19
MMT-1001	0212523379	01/08/19
MMT-923	5078337	15/08/17
MMT-923	5084127	15/09/17
MMT-923	5088984	15/10/17
MMT-923	5091306	15/11/17
MMT-923	5091307	15/11/17
MMT-923	5088612	15/12/17
MMT-923	5095071	15/01/18
MMT-923	5095072	15/01/18
MMT-923	5103315	15/03/18
MMT-923	5103313	15/03/18
MMT-923	5104633	15/04/18
MMT-923	5104634	15/04/18
MMT-923	5105363	15/04/18
MMT-923	5100590	15/04/18
MMT-923	5109796	15/06/18
MMT-923	5105437	15/06/18
MMT-923	5117939	15/07/18
MMT-923	5115441	15/07/18
MMT-923	5117936	15/08/18
MMT-923	5118541	15/08/18
MMT-923	5117938	15/08/18
MMT-923	5127871	15/09/18
MMT-923	5123227	15/09/18
MMT-923	5119663	15/10/18
MMT-923	5128411	15/10/18
MMT-923	5130065	15/11/18
MMT-923	5128409	15/11/18
MMT-923	5128413	15/11/18
MMT-923	5130067	15/11/18
MMT-923	5131857	15/11/18
MMT-923	5130066	15/11/18
MMT-923	5134549	15/12/18
MMT-923	5137570	15/01/19
MMT-923	5140387	15/02/19
MMT-923	5143988	15/03/19

MMT-923	5147486	15/04/19
MMT-923	5141512	15/04/19
MMT-923	5149496	15/05/19
MMT-923	5149494	15/05/19
MMT-923	5150471	15/05/19
MMT-923	5151831	15/06/19
MMT-923	5153938	15/06/19
MMT-923	5152963	15/06/19
MMT-923	5154654	15/06/19
MMT-923	5157423	15/07/19
MMT-923	5157372	15/07/19
MMT-923	5156447	15/07/19
MMT-923	5161815	01/08/19
MMT-923	5162224	01/08/19
MMT-923	5162251	01/09/19
MMT-923	5163246	01/09/19
MMT-923	5162252	01/09/19
MMT-923	5163932	01/09/19
MMT-923	5163258	01/09/19
MMT-923	5164574	01/09/19
MMT-923	5167806	01/11/19
MMT-923	5172578	01/11/19
MMT-923	5174511	01/12/19
MMT-923	5175708	01/12/19
MMT-923	5176471	01/01/20
MMT-923	5177288	01/01/20
MMT-923	5176470	01/01/20
MMT-923	5180046	01/02/20
MMT-923	5103889	15/04/18
MMT-925	5127872	15/09/18
MMT-925	5138172	15/01/19
MMT-925	5143413	15/03/19
MMT-925	5156832	15/07/19
MMT-925	5167807	01/11/19
MMT-925	5081819	15/08/17
MMT-925	5084130	15/09/17
MMT-925	5078343	15/10/17
MMT-925	5088986	15/10/17
MMT-925	5094245	15/12/17
MMT-925	5084466	15/01/18
MMT-925	5115457	15/07/18
MMT-925	5122403	15/09/18

MMT-925	5130069	15/11/18
MMT-925	5130068	15/11/18
MMT-925	5131258	15/11/18
MMT-925	5129134	15/12/18
MMT-925	5136422	15/01/19
MMT-925	5135134	15/02/19
MMT-925	5143990	15/04/19
MMT-925	5147488	15/04/19
MMT-925	5141470	15/04/19
MMT-925	5151302	15/05/19
MMT-925	5153939	15/06/19
MMT-925	5157425	15/07/19
MMT-925	5154270	01/08/19
MMT-925	5158196	01/08/19
MMT-925	5150558	01/10/19
MMT-925	5169324	01/11/19
MMT-925	5174516	01/12/19
MMT-925	5165957	01/12/19
MMT-925	5180047	01/02/20
MMT-941	5151303	15/05/19
MMT-941	5109799	15/05/18
MMT-941	5135136	15/01/19
MMT-941	5140396	15/02/19
MMT-941	5148154	15/04/19
MMT-941	5148888	15/04/19
MMT-941	5157426	15/07/19
MMT-941	5139011	15/02/19
MMT-941	5129131	15/10/18
MMT-941	5151833	15/06/19
MMT-941	5151834	15/06/19
MMT-941	5170779	01/11/19
MMT-941	5077112	15/08/17
MMT-941	5081818	15/08/17
MMT-941	5079704	15/08/17
MMT-941	5082776	15/08/17
MMT-941	5086975	15/10/17
MMT-941	5090810	15/12/17
MMT-941	5091308	15/12/17
MMT-941	5105444	15/04/18
MMT-941	5117218	15/07/18
MMT-941	5105446	15/07/18
MMT-941	5117940	15/08/18

MMT-941	5120634	15/08/18
MMT-941	5118543	15/08/18
MMT-941	5122582	15/09/18
MMT-941	5123241	15/10/18
MMT-941	5117246	15/10/18
MMT-941	5131259	15/11/18
MMT-941	5129132	15/11/18
MMT-941	5133985	15/12/18
MMT-941	5135135	15/12/18
MMT-941	5137370	15/01/19
MMT-941	5140272	15/02/19
MMT-941	5143414	15/03/19
MMT-941	5137876	15/03/19
MMT-941	5146899	15/04/19
MMT-941	5141601	15/04/19
MMT-941	5150474	15/05/19
MMT-941	5153940	15/06/19
MMT-941	5152965	15/06/19
MMT-941	5154656	15/06/19
MMT-941	5156833	15/07/19
MMT-941	5151276	15/07/19
MMT-941	5157370	15/07/19
MMT-941	5158198	01/08/19
MMT-941	5161814	01/08/19
MMT-941	5162253	01/09/19
MMT-941	5163951	01/09/19
MMT-941	5163950	01/09/19
MMT-941	5167808	01/10/19
MMT-941	5170780	01/11/19
MMT-941	5174517	01/12/19
MMT-941	5175709	01/12/19
MMT-941	5174518	01/12/19
MMT-941	5176577	01/01/20
MMT-941	5178186	01/01/20
MMT-941	5176568	01/01/20
MMT-941	5180048	01/02/20
MMT-941	5101627	15/02/18
MMT-943	5167134	01/10/19
MMT-943	5177292	01/01/20
MMT-943	5178026	01/01/20
MMT-943	5084471	15/10/17
MMT-943	5140401	15/02/19

MMT-943	5151305	15/05/19
MMT-943	5131863	15/11/18
MMT-943	5163934	01/09/19
MMT-943	5117254	15/08/18
MMT-943	5137371	15/01/19
MMT-943	5152556	01/08/19
MMT-943	5147490	15/04/19
MMT-943	5120635	15/08/18
MMT-943	5162205	01/09/19
MMT-943	5170782	01/11/19
MMT-943	5169337	01/11/19
MMT-943	5096425	15/01/18
MMT-943	5119936	15/08/18
MMT-943	5132966	15/12/18
MMT-943	5137373	15/01/19
MMT-943	5105366	15/04/18
MMT-943	5082272	15/08/17
MMT-943	5084469	15/09/17
MMT-943	5085190	15/09/17
MMT-943	5082615	15/09/17
MMT-943	5085191	15/09/17
MMT-943	5084128	15/09/17
MMT-943	5086976	15/09/17
MMT-943	5082442	15/09/17
MMT-943	5093159	15/12/17
MMT-943	5091310	15/12/17
MMT-943	5088605	15/12/17
MMT-943	5095076	15/12/17
MMT-943	5095078	15/01/18
MMT-943	5095075	15/01/18
MMT-943	5105177	15/04/18
MMT-943	5106189	15/04/18
MMT-943	5105161	15/05/18
MMT-943	5109803	15/05/18
MMT-943	5109720	15/06/18
MMT-943	5109729	15/07/18
MMT-943	5118545	15/08/18
MMT-943	5117941	15/08/18
MMT-943	5118544	15/08/18
MMT-943	5122409	15/09/18
MMT-943	5123236	15/09/18
MMT-943	5122408	15/09/18

MMT-943	5118959	15/10/18
MMT-943	5130074	15/11/18
MMT-943	5131262	15/11/18
MMT-943	5130072	15/11/18
MMT-943	5130073	15/11/18
MMT-943	5130075	15/11/18
MMT-943	5129479	15/11/18
MMT-943	5132085	15/12/18
MMT-943	5132081	15/12/18
MMT-943	5135140	15/01/19
MMT-943	5140400	15/02/19
MMT-943	5137889	15/02/19
MMT-943	5140275	15/02/19
MMT-943	5138482	15/02/19
MMT-943	5138483	15/02/19
MMT-943	5143415	15/03/19
MMT-943	5147491	15/04/19
MMT-943	5146895	15/04/19
MMT-943	5148157	15/04/19
MMT-943	5146894	15/04/19
MMT-943	5150476	15/05/19
MMT-943	5151304	15/05/19
MMT-943	5150478	15/05/19
MMT-943	5150477	15/05/19
MMT-943	5150479	15/05/19
MMT-943	5149499	15/05/19
MMT-943	5152966	15/06/19
MMT-943	5155466	15/06/19
MMT-943	5155465	15/06/19
MMT-943	5154668	15/06/19
MMT-943	5151836	15/06/19
MMT-943	5156449	15/07/19
MMT-943	5156836	15/07/19
MMT-943	5156835	15/07/19
MMT-943	5157428	15/07/19
MMT-943	5156834	15/07/19
MMT-943	5157430	15/07/19
MMT-943	5157429	15/07/19
MMT-943	5161465	01/08/19
MMT-943	5161070	01/08/19
MMT-943	5162236	01/08/19
MMT-943	5163250	01/09/19

MMT-943	5163935	01/09/19
MMT-943	5163249	01/09/19
MMT-943	5162833	01/09/19
MMT-943	5162249	01/09/19
MMT-943	5167132	01/10/19
MMT-943	5165211	01/10/19
MMT-943	5165210	01/10/19
MMT-943	5167135	01/10/19
MMT-943	5167820	01/10/19
MMT-943	5165215	01/10/19
MMT-943	5170783	01/11/19
MMT-943	5169338	01/11/19
MMT-943	5174523	01/12/19
MMT-943	5176578	01/01/20
MMT-943	5180052	01/02/20
MMT-943	5180053	01/02/20
MMT-943	5132087	15/12/18
MMT-943	5137372	15/01/19
MMT-1001	0212524900	01/08/19
MMT-945	5139013	15/01/19
MMT-945	5154659	15/06/19
MMT-945	5175705	01/12/19
MMT-945	5174524	01/12/19
MMT-945	5085192	15/09/17
MMT-945	5085477	15/09/17
MMT-945	5086710	15/12/17
MMT-945	5091312	15/12/17
MMT-945	5097971	15/02/18
MMT-945	5098833	15/02/18
MMT-945	5103412	15/04/18
MMT-945	5109809	15/05/18
MMT-945	5115448	15/06/18
MMT-945	5109810	15/06/18
MMT-945	5117220	15/07/18
MMT-945	5117943	15/07/18
MMT-945	5120637	15/08/18
MMT-945	5122410	15/09/18
MMT-945	5122583	15/09/18
MMT-945	5129126	15/10/18
MMT-945	5129125	15/10/18
MMT-945	5131867	15/11/18
MMT-945	5129127	15/11/18



MMT-945	5132968	15/12/18
MMT-945	5134551	15/12/18
MMT-945	5136424	15/01/19
MMT-945	5147493	15/04/19
MMT-945	5150481	15/05/19
MMT-945	5140360	15/05/19
MMT-945	5157365	15/07/19
MMT-945	5156451	15/07/19
MMT-945	5157433	15/07/19
MMT-945	5162254	01/08/19
MMT-945	5150908	01/08/19
MMT-945	5164578	01/09/19
MMT-945	5163260	01/09/19
MMT-945	5151297	01/10/19
MMT-945	5170785	01/11/19
MMT-945	5166417	01/01/20
MMT-945	5180056	01/02/20
MMT-965	5117046	15/08/18
MMT-965	5132969	15/11/18
MMT-965	5138175	15/01/19
MMT-965	5143416	15/03/19
MMT-965	5156452	15/07/19
MMT-965	5140403	15/03/19
MMT-965	5148159	15/04/19
MMT-965	5082273	15/08/17
MMT-965	5084131	15/09/17
MMT-965	5095079	15/01/18
MMT-965	5101632	15/02/18
MMT-965	5103887	15/03/18
MMT-965	5109812	15/06/18
MMT-965	5105460	15/06/18
MMT-965	5109740	15/07/18
MMT-965	5109742	15/09/18
MMT-965	5123232	15/10/18
MMT-965	5123231	15/10/18
MMT-965	5123238	15/11/18
MMT-965	5135142	15/12/18
MMT-965	5136428	15/01/19
MMT-965	5139014	15/01/19
MMT-965	5144000	15/03/19
MMT-965	5148890	15/04/19
MMT-965	5153943	15/06/19

MMT-965	5152968	15/06/19
MMT-965	5155467	15/07/19
MMT-965	5157367	15/07/19
MMT-965	5157434	15/07/19
MMT-965	5161463	01/08/19
MMT-965	5152564	01/09/19
MMT-965	5165217	01/10/19
MMT-965	5175720	01/12/19
MMT-965	5176570	01/01/20
MMT-965	5178197	01/02/20
MMT-965	5095189	15/03/18
MMT-975	5130079	15/11/18
MMT-975	5148162	15/04/19
MMT-975	5172584	01/12/19
MMT-975	5085478	15/09/17
MMT-975	5085480	15/09/17
MMT-975	5091225	15/11/17
MMT-975	5086708	15/12/17
MMT-975	5093161	15/12/17
MMT-975	5095081	15/12/17
MMT-975	5101241	15/03/18
MMT-975	5103321	15/03/18
MMT-975	5103416	15/04/18
MMT-975	5105145	15/04/18
MMT-975	5106836	15/05/18
MMT-975	5115460	15/07/18
MMT-975	5115450	15/07/18
MMT-975	5109745	15/07/18
MMT-975	5118549	15/08/18
MMT-975	5109758	15/08/18
MMT-975	5119237	15/08/18
MMT-975	5119238	15/08/18
MMT-975	5122177	15/09/18
MMT-975	5118932	15/09/18
MMT-975	5128918	15/10/18
MMT-975	5134426	15/01/19
MMT-975	5139015	15/01/19
MMT-975	5135144	15/01/19
MMT-975	5137568	15/01/19
MMT-975	5140276	15/02/19
MMT-975	5140409	15/02/19
MMT-975	5144005	15/03/19

MMT-975	5147495	15/04/19
MMT-975	5148160	15/04/19
MMT-975	5143667	15/04/19
MMT-975	5154661	15/06/19
MMT-975	5150439	15/07/19
MMT-975	5152338	15/07/19
MMT-975	5152354	01/08/19
MMT-975	5158205	15/08/19
MMT-975	5162239	01/09/19
MMT-975	5163938	01/09/19
MMT-975	5165226	01/10/19
MMT-975	5169332	01/11/19
MMT-975	5176574	01/01/20
MMT-975	5178191	01/01/20
MMT-1000	0209918204	15/03/18
MMT-1000	0210805376	15/10/18
MMT-1000	0209796050	15/02/18
MMT-1000	0209800546	15/02/18
MMT-1000	0209806385	15/02/18
MMT-1000	0209825072	15/02/18
MMT-1000	0209830199	15/02/18
MMT-1000	0209831816	15/02/18
MMT-1000	0209836347	15/02/18
MMT-1000	0209842963	15/02/18
MMT-1000	0209903829	15/03/18
MMT-1000	0209907372	15/03/18
MMT-1000	0209907374	15/03/18
MMT-1000	0209907377	15/03/18
MMT-1000	0209936029	15/03/18
MMT-1000	0209936055	15/03/18
MMT-1000	0209969900	15/03/18
MMT-1000	0209972008	15/03/18
MMT-1000	0210031884	15/04/18
MMT-1000	0210049681	15/04/18
MMT-1000	0210068544	15/04/18
MMT-1000	0210077688	15/04/18
MMT-1000	0210109050	15/03/18
MMT-1001	0211723685	15/04/19
MMT-1001	0211729178	15/04/19
MMT-1001	0211907525	15/05/19
MMT-1001	0211912287	15/06/19
MMT-1001	0213085718	01/10/19

MMT-1001	0210051873	15/08/17
MMT-1001	0210058751	15/08/17
MMT-1001	0210101132	15/08/17
MMT-1001	0210101133	15/08/17
MMT-1001	0210101134	15/08/17
MMT-1001	0210128899	15/09/17
MMT-1001	0210128940	15/08/17
MMT-1001	0210147247	15/06/18
MMT-1001	0210172656	15/05/18
MMT-1001	0210245918	15/06/18
MMT-1001	0210246107	15/07/18
MMT-1001	0210246110	15/07/18
MMT-1001	0210250628	15/07/18
MMT-1001	0210290356	15/07/18
MMT-1001	0210290357	15/05/18
MMT-1001	0210291302	15/05/18
MMT-1001	0210294912	15/06/18
MMT-1001	0210328771	15/07/18
MMT-1001	0210364483	15/07/18
MMT-1001	0210373047	15/08/18
MMT-1001	0210445107	15/08/18
MMT-1001	0210451698	15/08/18
MMT-1001	0210452643	15/08/18
MMT-1001	0210460991	15/09/18
MMT-1001	0210468571	15/09/18
MMT-1001	0210474221	15/09/18
MMT-1001	0210482613	15/09/18
MMT-1001	0210488073	15/09/18
MMT-1001	0210563848	15/09/18
MMT-1001	0210587012	15/10/18
MMT-1001	0210587016	15/10/18
MMT-1001	0210594151	15/10/18
MMT-1001	0210597823	15/10/18
MMT-1001	0210601639	15/10/18
MMT-1001	0210607422	15/10/18
MMT-1001	0210672557	15/10/18
MMT-1001	0210672558	15/11/18
MMT-1001	0210672559	15/11/18
MMT-1001	0210675761	15/10/18
MMT-1001	0210681111	15/11/18
MMT-1001	0210738311	15/11/18
MMT-1001	0210741232	15/11/18

MMT-1001	0210763716	15/11/18
MMT-1001	0210852240	15/11/18
MMT-1001	0210852242	15/11/18
MMT-1001	0210853928	15/12/18
MMT-1001	0210854260	15/12/18
MMT-1001	0210854276	15/11/18
MMT-1001	0210854277	15/12/18
MMT-1001	0210889185	15/12/18
MMT-1001	0210890079	15/12/18
MMT-1001	0210960831	15/12/18
MMT-1001	0210981468	15/12/18
MMT-1001	0210995950	15/12/18
MMT-1001	0211013745	15/12/18
MMT-1001	0211039389	15/12/18
MMT-1001	0211039395	15/11/18
MMT-1001	0211058319	15/01/19
MMT-1001	0211072724	15/01/19
MMT-1001	0211073799	15/01/19
MMT-1001	0211089296	15/01/19
MMT-1001	0211117175	15/01/19
MMT-1001	0211128858	15/01/19
MMT-1001	0211128859	15/01/19
MMT-1001	0211154587	15/01/19
MMT-1001	0211157890	15/01/19
MMT-1001	0211212998	15/02/19
MMT-1001	0211213321	15/01/19
MMT-1001	0211285629	15/01/19
MMT-1001	0211285630	15/01/19
MMT-1001	0211289429	15/09/18
MMT-1001	0211293490	15/02/19
MMT-1001	0211321650	15/02/19
MMT-1001	0211321658	15/02/19
MMT-1001	0211332173	15/02/19
MMT-1001	0211336547	15/11/18
MMT-1001	0211338650	15/03/19
MMT-1001	0211338658	15/02/19
MMT-1001	0211355090	15/03/19
MMT-1001	0211370645	15/03/19
MMT-1001	0211376472	15/03/19
MMT-1001	0211405241	15/03/19
MMT-1001	0211460088	15/03/19
MMT-1001	0211460092	15/03/19

MMT-1001	0211473612	15/01/19
MMT-1001	0211528599	15/04/19
MMT-1001	0211536613	15/04/19
MMT-1001	0211536662	15/04/19
MMT-1001	0211565606	15/04/19
MMT-1001	0211565669	15/04/19
MMT-1001	0211667295	15/04/19
MMT-1001	0211718145	15/04/19
MMT-1001	0211718149	15/04/19
MMT-1001	0211718152	15/04/19
MMT-1001	0211768946	15/04/19
MMT-1001	0211769499	15/05/19
MMT-1001	0211782026	15/05/19
MMT-1001	0211801301	15/05/19
MMT-1001	0211809093	15/05/19
MMT-1001	0211824668	15/05/19
MMT-1001	0211907524	15/04/19
MMT-1001	0211912283	15/04/19
MMT-1001	0211912285	15/05/19
MMT-1001	0211912286	15/04/19
MMT-1001	0211912288	15/04/19
MMT-1001	0211913166	15/06/19
MMT-1001	0211913174	15/06/19
MMT-1001	0211927795	15/06/19
MMT-1001	0211935660	15/06/19
MMT-1001	0211977004	15/06/19
MMT-1001	0212038299	15/06/19
MMT-1001	0212066366	15/06/19
MMT-1001	0212524901	01/08/19
MMT-397	5165428	01/10/19
MMT-399	5164667	01/11/19
MMT-1001	0212532534	01/08/19
MMT-1001	0212109137	15/06/19
MMT-1001	0212110435	15/06/19
MMT-1001	0212110436	15/06/19
MMT-1001	0212112860	15/06/19
MMT-1001	0212124152	15/05/19
MMT-1001	0212124700	15/03/19
MMT-1001	0212124703	15/04/19
MMT-1001	0212256884	15/04/19
MMT-399	5164670	01/12/19
MMT-1001	0212537121	01/08/19

MMT-1001	0212296902	15/04/19
MMT-1001	0212538820	01/08/19
MMT-1001	0212834770	01/10/19
MMT-1001	0212856144	01/10/19
MMT-1001	0212532533	01/08/19
MMT-1001	0212970837	01/10/19
MMT-1001	0212561644	01/08/19
MMT-1001	0212628818	01/09/19
MMT-1001	0212717612	01/10/19
MMT-1001	0212829305	01/10/19
MMT-399	5170882	01/01/20
MMT-1001	0212981916	01/10/19
MMT-1001	0213055130	01/10/19
MMT-1001	0212965943	01/09/19
MMT-399	5170968	01/01/20
MMT-1001	0213064208	01/10/19
MMT-1001	0213058331	01/10/19
MMT-1001	0213064339	01/10/19
MMT-1001	0213208843	01/10/19
MMT-1001	0213395713	01/10/19
MMT-1001	0213094816	01/10/19
MMT-1001	0213098295	01/10/19
MMT-1001	0213103684	01/10/19
MMT-1001	0213426643	01/10/19
MMT-864	5180448	01/03/22
MMT-1001	0213524667	01/10/19
MMT-1001	0213380679	01/12/19
MMT-1001	0213524668	01/10/19
MMT-1001	0213524669	01/10/19
MMT-1001	0213524748	01/10/19
MMT-1001	0213584272	01/10/19
MMT-1001	0213584275	01/10/19
MMT-1001	0213721081	01/10/19
MMT-1001	0213721082	01/10/19
MMT-1001	0213352081	01/12/19
MMT-1001	0213361272	01/12/19
MMT-1001	0213379768	01/12/19
MMT-1001	0211213320	15/01/19
MMT-1001	0211285628	15/01/19