

«Hospital_Name»
«Users_Name» - «Department»
«Customer_Address»
«Zip_Code» «City» - «Country_name»

Reference: 91033431-FA

28 August 2015

Urgent Field Safety Notice Important Medical Device Information

WATCHMAN[®] Left Atrial Appendage Closure Device with Delivery System and WATCHMAN[®] Access System Access Sheath with Dilator

Dear «Users_Name»,

Boston Scientific (BSC) is initiating a voluntary Field Safety Notice for the WATCHMAN technology, which includes the WATCHMAN Left Atrial Appendage (LAA) Closure Device with Delivery System and the WATCHMAN Access System.

BSC has determined cross-threading of the hemostasis valve may occur if the valve is tightened with the dilator in place, potentially preventing subsequent sealing of the valve when desired. This Field Safety Notice reinforces existing Directions for Use (DFU) and provides further guidance regarding the correct use of the hemostasis valve in order to avoid cross-threading and to securely seal the valve, minimizing the potential for undesirable blood leakage.

No product is being recalled and you are not required to return product to Boston Scientific.

Note that there is no impact to previously implanted devices.

Action Required

It is very important that you read this entire Field Safety Notice and ensure that all users of the WATCHMAN LAA Closure Device with Delivery System and WATCHMAN Access System are aware of this Field Safety Notice (please refer to the affected list of products attached). **You must also complete the enclosed Acknowledgment Form and return it to Boston Scientific, indicating that you have received, read, and understood the important information contained in this Field Safety Notice.**

Updated Directions for Use

The following sections of the WATCHMAN LAA Closure Device with Delivery System and the WATCHMAN Access System DFUs clarify correct use of the hemostasis valve and interaction between the access sheath and dilator. The additions to the current DFU are highlighted in **blue**:

Applicable to WATCHMAN LAA Closure Device with Delivery System and the WATCHMAN Access System DFUs:

3. Prepare WATCHMAN Access System.

Note: Inspect sterile package and WATCHMAN Access System prior to use. If sterile barrier, labeling, packaging, or device have been compromised in any way, DO NOT USE.

- A. Remove access sheath and dilator from package under sterile conditions.
- B. Inspect prior to use to ensure no damage.
- C. Flush access sheath and dilator with sterile saline prior to use.
- D. Insert dilator into hemostasis valve of access sheath until the two snap together.

Note: Do not tighten the hemostasis valve while the dilator is inserted in the WATCHMAN Access System. The dilator by itself will occlude the lumen of the WATCHMAN Access System creating hemostasis. Tightening the valve onto the dilator may damage the valve threads, which can lead to subsequent difficulty in closing the valve and an incomplete seal, once the dilator is removed.

4. Advance WATCHMAN Access System over guidewire into left atrium (LA). As access sheath nears center of LA, unsnap the access sheath from the dilator, hold dilator and advance access sheath into initial position in LA or ostium of Left Upper Pulmonary Vein (LUPV).

PRECAUTION: Use caution when introducing the WATCHMAN Access System to prevent damage to cardiac structures.

5. Remove dilator and guidewire, leaving access sheath in LA or LUPV. Allow back bleed to minimize potential for introducing air before tightening the hemostasis valve. Flush the access sheath with saline.

If continued back bleed is observed from the valve after the dilator is removed despite attempting to close it, loosen the valve cap (counter-clockwise rotation) until the cap spins freely. Then re-attempt closure of the valve while exerting gentle forward pressure on the valve cap during closure (clockwise rotation) to ensure proper engagement of the valve thread. While these steps are being undertaken, manual occlusion of the valve opening using a gloved finger is recommended to minimize blood loss.

Note: These steps may be repeated if necessary. However, if this does not mitigate the blood leak, the user should remove and replace the WATCHMAN Access Sheath before proceeding with the procedure.

(Continue to follow existing DFU instructions)

Applicable to the WATCHMAN LAA Closure Device with Delivery System DFU only:

8. Loosen hemostasis valve of the WATCHMAN Access Sheath allowing back bleed before inserting the WATCHMAN Delivery System. *Note: Hemostasis valve should spin freely (fully open).*

Note: Tightening the valve onto the WATCHMAN Delivery System may damage the valve threads, which can lead to subsequent difficulty in closing the valve and an incomplete seal, once the WATCHMAN Delivery System is removed.

(Continue to follow existing DFU instructions)

Our records indicate that your facility received one or more of the affected products subject to this Field Safety Notice. **The table below provides a complete list of all affected products**, including Product Description, Material Number (UPN), Lot/Batch numbers and Expiration date. **Please note that only the material listed in the table below is affected.**
No other Boston Scientific product is involved by this Field Safety Notice.

Product Description	Material/UPN	Lot/Batch Number	Expiration Date
WATCHMAN® Left Atrial Appendage Closure Device with Delivery System	M635WS21060	See attached list	
	M635WS24060		
	M635WS27060		
	M635WS30060		
	M635WS33060		
WATCHMAN® Access System	M635TS10060		
	M635TS20060		
	M635TS40060		

INSTRUCTIONS:

- Please read carefully the Field Safety Notice and immediately post this information in a visible location near the product to ensure this information is easily accessible to all users of the device.**
- Please complete the attached Acknowledgement Form even if you do not have any affected product.**
- When completed, please return the Acknowledgement Form to your local Boston Scientific office to the attention of «Customer_Service_Fax_Number» on or before **14 September 2015.****
- Please pass on this notice to any health professional of your organization that need to be aware and to any organization where the potentially affected devices have been transferred (if appropriate). Please provide Boston Scientific with details of any affected devices that have been transferred to other organizations (if appropriate).

Although Boston Scientific is not physically recalling any product, your Competent Authority is being notified of this Field Safety Notice.

We regret any inconvenience that this action may cause and we appreciate your understanding as we take action to ensure patient safety and customer satisfaction.

If you have any questions or would like assistance with this Field Safety Notice, please contact your local Sales Representative.

Yours sincerely,



Marie Pierre Barlangua
 Quality Department
 Boston Scientific International S.A.

Attachments: - Acknowledgement Form
 - List of affected product



Please complete the form & send it to Your Local Office:
«Customer_Service_Fax_Number»

«Sold_to» - «Hospital_Name» - «City» - «Country_name»

Acknowledgement Form – Important Medical Device Information

**WATCHMAN® Left Atrial Appendage Closure Device with Delivery System
and
WATCHMAN® Access System Access Sheath with Dilator**

91033431-FA

I acknowledge receipt of the Boston Scientific Field Safety Notice

dated 28 August 2015 for the

**WATCHMAN® Left Atrial Appendage Closure Device
with Delivery System**

and

WATCHMAN® Access System Access Sheath with Dilator

and took action as required in the “Instructions” of the letter.

NAME* _____ **Title** _____

Telephone _____ **Department** _____

Customer' SIGNATURE* _____ **DATE*** _____

* Required field

DD/MM/YYYY

91033431-FA WATCHMAN® Left Atrial Appendage Closure Device with Delivery System and WATCHMAN®
 Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

Material	Product Description	Batch	Catalog #	Expiration Date
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15711285	TS-1006	28-Nov-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15722505	TS-1006	3-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15716329	TS-1006	4-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15733228	TS-1006	6-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15742846	TS-1006	11-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15746783	TS-1006	12-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15750090	TS-1006	13-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15760116	TS-1006	18-Dec-2015
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15785032	TS-1006	8-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15802490	TS-1006	15-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15808197	TS-1006	17-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15821967	TS-1006	23-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15825660	TS-1006	24-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15833585	TS-1006	28-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15839839	TS-1006	30-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524103	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524104	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524105	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524106	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524107	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524108	TS-1006	31-Jan-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15870735	TS-1006	11-Feb-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15906798	TS-1006	26-Feb-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524109	TS-1006	29-Feb-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524404	TS-1006	29-Feb-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15942740	TS-1006	11-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15948999	TS-1006	13-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15952034	TS-1006	14-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15957983	TS-1006	18-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15961293	TS-1006	19-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524113	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524115	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524118	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524400	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17524402	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17539282	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17539283	TS-1006	31-Mar-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	15992933	TS-1006	2-Apr-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16028407	TS-1006	16-Apr-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16035580	TS-1006	18-Apr-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17538716	TS-1006	30-Apr-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17538719	TS-1006	30-Apr-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16135267	TS-1006	30-May-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16145323	TS-1006	4-Jun-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16164714	TS-1006	12-Jun-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16168810	TS-1006	13-Jun-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16172942	TS-1006	16-Jun-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16179406	TS-1006	18-Jun-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16242184	TS-1006	16-Jul-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16256223	TS-1006	22-Jul-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16277475	TS-1006	30-Jul-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16409463	TS-1006	24-Sep-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16422803	TS-1006	30-Sep-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16423442	TS-1006	30-Sep-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16576447	TS-1006	1-Dec-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16627421	TS-1006	26-Dec-2016
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16691145	TS-1006	26-Jan-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16702549	TS-1006	29-Jan-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16706056	TS-1006	30-Jan-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17538717	TS-1006	31-Jan-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17539281	TS-1006	31-Jan-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16719766	TS-1006	5-Feb-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16743849	TS-1006	13-Feb-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16764290	TS-1006	23-Feb-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16773649	TS-1006	25-Feb-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16773650	TS-1006	25-Feb-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16794371	TS-1006	5-Mar-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16794372	TS-1006	5-Mar-2017

91033431-FA WATCHMAN® Left Atrial Appendage Closure Device with Delivery System and WATCHMAN®
 Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16840109	TS-1006	23-Mar-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16842573	TS-1006	24-Mar-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16852558	TS-1006	27-Mar-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16865853	TS-1006	1-Apr-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16875735	TS-1006	6-Apr-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16930186	TS-1006	28-Apr-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16983795	TS-1006	19-May-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	16987369	TS-1006	20-May-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17014384	TS-1006	1-Jun-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17170094	TS-1006	3-Aug-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17246303	TS-1006	1-Sep-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17265821	TS-1006	8-Sep-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	17313481	TS-1006	24-Sep-2017
M635TS10060	WATCHMAN Access Single Curve-14F (SA)	18158641	TS-1006	31-Jul-2018

Material	Product Description	Batch	Catalog #	Expiration Date
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00543403	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00543404	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00543409	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00543410	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00543411	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544801	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544802	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544804	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544805	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544806	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544901	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544902	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544903	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544904	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00544905	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546002	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546004	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546005	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546006	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546007	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546008	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546011	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00546012	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547702	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547704	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547705	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547706	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547710	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547711	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547712	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547713	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00547715	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548001	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548003	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548004	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548005	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548006	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548007	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548008	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548009	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548010	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548011	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548012	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548014	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	F00548015	WS-2406	31-Aug-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15533157	WS-2406	14-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15533159	WS-2406	14-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15534567	WS-2406	14-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15560537	WS-2406	25-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15560538	WS-2406	25-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15563850	WS-2406	26-Sep-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15584346	WS-2406	4-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15601788	WS-2406	11-Oct-2015

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15613894	WS-2406	17-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15624071	WS-2406	22-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15627321	WS-2406	23-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15630804	WS-2406	24-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15639756	WS-2406	29-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15645827	WS-2406	31-Oct-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15650908	WS-2406	1-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15656641	WS-2406	5-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15660122	WS-2406	6-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15664112	WS-2406	7-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15670042	WS-2406	8-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15675296	WS-2406	12-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15678800	WS-2406	13-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15694583	WS-2406	20-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15712768	WS-2406	28-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17510339	WS-2406	30-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17513899	WS-2406	30-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17514330	WS-2406	30-Nov-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15725135	WS-2406	4-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15732982	WS-2406	6-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15752943	WS-2406	14-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15757059	WS-2406	17-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15761457	WS-2406	18-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15764438	WS-2406	19-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15767310	WS-2406	21-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15768509	WS-2406	26-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17510421	WS-2406	31-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17514339	WS-2406	31-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17514544	WS-2406	31-Dec-2015
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15785030	WS-2406	8-Jan-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15802897	WS-2406	15-Jan-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15828835	WS-2406	25-Jan-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15837276	WS-2406	29-Jan-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15839943	WS-2406	30-Jan-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15873851	WS-2406	12-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15880319	WS-2406	14-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15892274	WS-2406	20-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15904005	WS-2406	25-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15905274	WS-2406	25-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15907468	WS-2406	26-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17514540	WS-2406	29-Feb-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15917731	WS-2406	1-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15923968	WS-2406	4-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15927422	WS-2406	5-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15942658	WS-2406	11-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15945694	WS-2406	12-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15952033	WS-2406	14-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15972750	WS-2406	24-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15980371	WS-2406	26-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15982572	WS-2406	27-Mar-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15986988	WS-2406	1-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	15996021	WS-2406	3-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16008893	WS-2406	8-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16009246	WS-2406	8-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16012590	WS-2406	9-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16025225	WS-2406	15-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16028569	WS-2406	16-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16031642	WS-2406	17-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16049990	WS-2406	24-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16050155	WS-2406	24-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16053803	WS-2406	25-Apr-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16065499	WS-2406	1-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16075310	WS-2406	6-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16095358	WS-2406	14-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16110279	WS-2406	20-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16113543	WS-2406	21-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16128522	WS-2406	28-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16135270	WS-2406	30-May-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16138985	WS-2406	2-Jun-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16141220	WS-2406	3-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16155838	WS-2406	9-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16161519	WS-2406	11-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16179407	WS-2406	18-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16187284	WS-2406	20-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16196931	WS-2406	25-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16201410	WS-2406	26-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16204071	WS-2406	27-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17516997	WS-2406	30-Jun-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16210174	WS-2406	1-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16225593	WS-2406	9-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16231399	WS-2406	11-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16239135	WS-2406	15-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16245795	WS-2406	17-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16249321	WS-2406	18-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16260101	WS-2406	23-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16264696	WS-2406	24-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16274329	WS-2406	29-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17516991	WS-2406	31-Jul-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16282440	WS-2406	1-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16284578	WS-2406	4-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16294524	WS-2406	7-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16304635	WS-2406	12-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16316948	WS-2406	18-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16326971	WS-2406	21-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16338000	WS-2406	26-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16342273	WS-2406	27-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17517172	WS-2406	31-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17517421	WS-2406	31-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17517701	WS-2406	31-Aug-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16405297	WS-2406	23-Sep-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16408162	WS-2406	24-Sep-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16436680	WS-2406	6-Oct-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16443230	WS-2406	8-Oct-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16447292	WS-2406	9-Oct-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16457877	WS-2406	14-Oct-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16511475	WS-2406	4-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16524446	WS-2406	10-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16529379	WS-2406	11-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16543353	WS-2406	17-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16551893	WS-2406	19-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16551899	WS-2406	19-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16559456	WS-2406	21-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16566248	WS-2406	25-Nov-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16619148	WS-2406	17-Dec-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16621808	WS-2406	18-Dec-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16623958	WS-2406	19-Dec-2016
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16635565	WS-2406	5-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16638767	WS-2406	6-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16649394	WS-2406	9-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16654988	WS-2406	12-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16675496	WS-2406	20-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16676236	WS-2406	20-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16684502	WS-2406	22-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16691148	WS-2406	26-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17588071	WS-2406	26-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16701580	WS-2406	29-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17504979	WS-2406	31-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17505620	WS-2406	31-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17505627	WS-2406	31-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17522737	WS-2406	31-Jan-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16716474	WS-2406	4-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16723383	WS-2406	6-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16735187	WS-2406	11-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16740111	WS-2406	12-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16752662	WS-2406	18-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16764291	WS-2406	23-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16775462	WS-2406	26-Feb-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16790634	WS-2406	4-Mar-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16797991	WS-2406	6-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16799659	WS-2406	6-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16803524	WS-2406	9-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16805225	WS-2406	9-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16812946	WS-2406	12-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16823549	WS-2406	17-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16842574	WS-2406	24-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16857215	WS-2406	30-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16858004	WS-2406	30-Mar-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16897924	WS-2406	14-Apr-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16900811	WS-2406	15-Apr-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16902464	WS-2406	15-Apr-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16943540	WS-2406	4-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16953117	WS-2406	7-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16967573	WS-2406	13-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16987245	WS-2406	20-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	16999910	WS-2406	26-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17003018	WS-2406	27-May-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17015191	WS-2406	1-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17020296	WS-2406	3-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17035415	WS-2406	9-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17056856	WS-2406	17-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17075341	WS-2406	24-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17091393	WS-2406	30-Jun-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17101485	WS-2406	6-Jul-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17139718	WS-2406	20-Jul-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17142697	WS-2406	21-Jul-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17164903	WS-2406	30-Jul-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17167874	WS-2406	31-Jul-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17175419	WS-2406	4-Aug-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17199427	WS-2406	13-Aug-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17232655	WS-2406	26-Aug-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17246388	WS-2406	1-Sep-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17268931	WS-2406	9-Sep-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17308289	WS-2406	23-Sep-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17323937	WS-2406	29-Sep-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17347993	WS-2406	7-Oct-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17360281	WS-2406	31-Oct-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17400535	WS-2406	31-Oct-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17404703	WS-2406	31-Oct-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17426248	WS-2406	30-Nov-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17445090	WS-2406	30-Nov-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17462355	WS-2406	30-Nov-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17513733	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17547678	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17549055	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17562117	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17563262	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17567816	WS-2406	31-Dec-2017
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17741851	WS-2406	14-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17585954	WS-2406	31-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17598806	WS-2406	31-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17607108	WS-2406	31-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17628229	WS-2406	31-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17646696	WS-2406	31-Jan-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17662280	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17673314	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17684461	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17687772	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17695516	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17734422	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17734431	WS-2406	28-Feb-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17757446	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17759451	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17780348	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17801349	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17815706	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17847860	WS-2406	31-Mar-2018
M635WS24060	FG WATCHMAN LAA Closure Device,24mm,OUS	17852577	WS-2406	30-Apr-2018

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00546912	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547201	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547203	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547204	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547205	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547206	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547209	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547210	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547301	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547302	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547303	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547304	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547305	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547306	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547307	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547308	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547309	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547310	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547401	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547404	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00547405	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549501	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549503	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549504	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549505	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549506	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549507	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549508	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	F00549509	WS-3306	31-Aug-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15533431	WS-3306	14-Sep-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15555660	WS-3306	24-Sep-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15556765	WS-3306	24-Sep-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15558791	WS-3306	25-Sep-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15576412	WS-3306	2-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15579548	WS-3306	3-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15586576	WS-3306	5-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15608022	WS-3306	15-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15610675	WS-3306	16-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15637569	WS-3306	26-Oct-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15650909	WS-3306	1-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15678801	WS-3306	13-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15684755	WS-3306	15-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15697854	WS-3306	21-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15711401	WS-3306	28-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17510337	WS-3306	30-Nov-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15732981	WS-3306	6-Dec-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15739114	WS-3306	10-Dec-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15770115	WS-3306	28-Dec-2015
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15772889	WS-3306	2-Jan-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15775655	WS-3306	3-Jan-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17513887	WS-3306	29-Feb-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	15917732	WS-3306	1-Mar-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16022223	WS-3306	14-Apr-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16045735	WS-3306	23-Apr-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16072787	WS-3306	5-May-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16081851	WS-3306	8-May-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16102079	WS-3306	16-May-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16102965	WS-3306	16-May-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16119791	WS-3306	23-May-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16141221	WS-3306	3-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16148063	WS-3306	5-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16165647	WS-3306	12-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16179408	WS-3306	18-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16187286	WS-3306	20-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16204068	WS-3306	27-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17517168	WS-3306	30-Jun-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16208090	WS-3306	1-Jul-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16210176	WS-3306	1-Jul-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16218850	WS-3306	7-Jul-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16242465	WS-3306	16-Jul-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16249320	WS-3306	18-Jul-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16271045	WS-3306	28-Jul-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16282651	WS-3306	1-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16291980	WS-3306	6-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16307822	WS-3306	13-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16319990	WS-3306	19-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16334236	WS-3306	25-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16342278	WS-3306	27-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17512672	WS-3306	31-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17517703	WS-3306	31-Aug-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16421082	WS-3306	29-Sep-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16524448	WS-3306	10-Nov-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16605572	WS-3306	11-Dec-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16626167	WS-3306	25-Dec-2016
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16638768	WS-3306	6-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16638769	WS-3306	6-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16642108	WS-3306	7-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17588074	WS-3306	7-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16655802	WS-3306	12-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16655806	WS-3306	12-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16676239	WS-3306	20-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16678436	WS-3306	21-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16684504	WS-3306	22-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16691076	WS-3306	26-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16706055	WS-3306	30-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17505940	WS-3306	31-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17517705	WS-3306	31-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17523741	WS-3306	31-Jan-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16709409	WS-3306	2-Feb-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16719786	WS-3306	5-Feb-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16743851	WS-3306	13-Feb-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16755490	WS-3306	19-Feb-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16787948	WS-3306	3-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16797993	WS-3306	6-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16799942	WS-3306	6-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16811005	WS-3306	11-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16816971	WS-3306	13-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16835335	WS-3306	20-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16845892	WS-3306	25-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16861290	WS-3306	31-Mar-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16886782	WS-3306	9-Apr-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16898028	WS-3306	14-Apr-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16904755	WS-3306	16-Apr-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16909513	WS-3306	20-Apr-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16909514	WS-3306	20-Apr-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16949829	WS-3306	6-May-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16983797	WS-3306	19-May-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16983799	WS-3306	19-May-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	16994017	WS-3306	22-May-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17027981	WS-3306	5-Jun-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17060695	WS-3306	18-Jun-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17060696	WS-3306	18-Jun-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17078717	WS-3306	25-Jun-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17091552	WS-3306	30-Jun-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17107949	WS-3306	8-Jul-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17152260	WS-3306	24-Jul-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17176021	WS-3306	4-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17185583	WS-3306	7-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214494	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214495	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214496	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214497	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214498	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17214499	WS-3306	19-Aug-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17272647	WS-3306	10-Sep-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17276356	WS-3306	11-Sep-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17282618	WS-3306	14-Sep-2017
M635WS33060	FG WATCHMAN LAA Closure Device, 33mm,OUS	17286603	WS-3306	15-Sep-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548604	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548605	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548606	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548607	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548608	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548609	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548610	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548701	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548702	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548704	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548705	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548706	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548708	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548709	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00548710	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549601	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549602	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549603	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549605	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549606	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549607	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549608	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	F00549609	WS-3006	31-Aug-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15533414	WS-3006	14-Sep-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15534507	WS-3006	14-Sep-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15554778	WS-3006	24-Sep-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15559112	WS-3006	25-Sep-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15567409	WS-3006	27-Sep-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15576410	WS-3006	2-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15513004	WS-3006	3-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15583465	WS-3006	4-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15586522	WS-3006	5-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15596777	WS-3006	10-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15597797	WS-3006	10-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15597799	WS-3006	10-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15601789	WS-3006	11-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15634207	WS-3006	25-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15636698	WS-3006	26-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15641980	WS-3006	30-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17510336	WS-3006	31-Oct-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15639840	WS-3006	15-Nov-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15691176	WS-3006	19-Nov-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15704330	WS-3006	26-Nov-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15709031	WS-3006	27-Nov-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17510338	WS-3006	30-Nov-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15726789	WS-3006	4-Dec-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15742845	WS-3006	11-Dec-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15768979	WS-3006	26-Dec-2015
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15775654	WS-3006	3-Jan-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15780644	WS-3006	7-Jan-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15785551	WS-3006	8-Jan-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15866546	WS-3006	8-Feb-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15899484	WS-3006	22-Feb-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15904700	WS-3006	25-Feb-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15921283	WS-3006	3-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15938317	WS-3006	10-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15945695	WS-3006	12-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15961294	WS-3006	19-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15976299	WS-3006	25-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15980377	WS-3006	26-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15980378	WS-3006	26-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	15982757	WS-3006	27-Mar-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16005082	WS-3006	7-Apr-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16018189	WS-3006	11-Apr-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16035581	WS-3006	18-Apr-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16043830	WS-3006	22-Apr-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16049994	WS-3006	24-Apr-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16066321	WS-3006	1-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16069158	WS-3006	2-May-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16092736	WS-3006	13-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16099085	WS-3006	15-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16119789	WS-3006	23-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16128521	WS-3006	28-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16131917	WS-3006	29-May-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16144975	WS-3006	4-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16150819	WS-3006	6-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16172941	WS-3006	16-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17587903	WS-3006	16-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16176188	WS-3006	17-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16194170	WS-3006	24-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16200327	WS-3006	26-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17489863	WS-3006	30-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17489864	WS-3006	30-Jun-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16213051	WS-3006	2-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16218855	WS-3006	7-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16235346	WS-3006	14-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16253200	WS-3006	21-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16256396	WS-3006	22-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16260102	WS-3006	23-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16274330	WS-3006	29-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16277970	WS-3006	30-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17523742	WS-3006	31-Jul-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16282535	WS-3006	1-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16301056	WS-3006	11-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16304511	WS-3006	12-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16311248	WS-3006	14-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16317993	WS-3006	18-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16319993	WS-3006	19-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16334234	WS-3006	25-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16342277	WS-3006	27-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17512673	WS-3006	31-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17517702	WS-3006	31-Aug-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16418477	WS-3006	29-Sep-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16552008	WS-3006	19-Nov-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16605571	WS-3006	11-Dec-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16611631	WS-3006	15-Dec-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16621807	WS-3006	18-Dec-2016
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16629591	WS-3006	1-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16645865	WS-3006	8-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16649395	WS-3006	9-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16676238	WS-3006	20-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16698284	WS-3006	28-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17485956	WS-3006	31-Jan-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16713085	WS-3006	3-Feb-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16727779	WS-3006	9-Feb-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16740110	WS-3006	12-Feb-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16773087	WS-3006	25-Feb-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16797992	WS-3006	6-Mar-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16799941	WS-3006	6-Mar-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16816970	WS-3006	13-Mar-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16833973	WS-3006	20-Mar-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16835333	WS-3006	20-Mar-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16863192	WS-3006	1-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16879810	WS-3006	7-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16883020	WS-3006	8-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16911586	WS-3006	21-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16915196	WS-3006	22-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16930188	WS-3006	28-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16931943	WS-3006	28-Apr-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16964404	WS-3006	12-May-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	16979303	WS-3006	18-May-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17013672	WS-3006	1-Jun-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17015194	WS-3006	1-Jun-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17049932	WS-3006	15-Jun-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17064564	WS-3006	19-Jun-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17068256	WS-3006	22-Jun-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17103307	WS-3006	6-Jul-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17125140	WS-3006	14-Jul-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17138031	WS-3006	20-Jul-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17152067	WS-3006	24-Jul-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17155705	WS-3006	27-Jul-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17182781	WS-3006	6-Aug-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17235441	WS-3006	26-Aug-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17240941	WS-3006	28-Aug-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17249478	WS-3006	2-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17274686	WS-3006	10-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17274689	WS-3006	10-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17284571	WS-3006	15-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17315901	WS-3006	25-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17327649	WS-3006	30-Sep-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17418948	WS-3006	30-Nov-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17435610	WS-3006	30-Nov-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17449151	WS-3006	30-Nov-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17474986	WS-3006	30-Nov-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17505149	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17505151	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17512095	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17516438	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17555483	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17563263	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17563434	WS-3006	31-Dec-2017
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17580530	WS-3006	31-Jan-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17598807	WS-3006	31-Jan-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17607463	WS-3006	31-Jan-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17628583	WS-3006	31-Jan-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17643288	WS-3006	31-Jan-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17766401	WS-3006	31-Mar-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17785911	WS-3006	31-Mar-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17829286	WS-3006	31-Mar-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17838466	WS-3006	31-Mar-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17858126	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17866074	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17881676	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17904999	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17914603	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17935423	WS-3006	30-Apr-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17969722	WS-3006	31-May-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17989461	WS-3006	31-May-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	17994394	WS-3006	31-May-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18005577	WS-3006	31-May-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18059018	WS-3006	30-Jun-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18075560	WS-3006	30-Jun-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18164716	WS-3006	31-Jul-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18187095	WS-3006	31-Jul-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18224198	WS-3006	31-Jul-2018
M635WS30060	FG WATCHMAN LAA Closure Device, 30mm,OUS	18268761	WS-3006	31-Aug-2018

Material	Product Description	Batch	Catalog #	Expiration Date
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544001	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544102	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544103	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544104	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544105	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544109	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544110	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544111	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544602	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544603	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544604	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544605	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544609	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544610	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00544611	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00545001	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00545006	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00545007	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00545009	WS-2706	31-Aug-2015

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548303	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548304	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548305	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548306	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548309	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548310	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548311	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548312	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548313	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	F00548314	WS-2706	31-Aug-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15534565	WS-2706	14-Sep-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15534566	WS-2706	14-Sep-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15560536	WS-2706	25-Sep-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15567410	WS-2706	27-Sep-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15573288	WS-2706	1-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15601787	WS-2706	11-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15610674	WS-2706	16-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15613896	WS-2706	17-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15617630	WS-2706	18-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15617819	WS-2706	18-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15641459	WS-2706	30-Oct-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15624074	WS-2706	1-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15649001	WS-2706	1-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15656642	WS-2706	5-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15659939	WS-2706	6-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15663165	WS-2706	7-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15670041	WS-2706	8-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15675294	WS-2706	12-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15691175	WS-2706	19-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15704329	WS-2706	26-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17506121	WS-2706	30-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17506122	WS-2706	30-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17510420	WS-2706	30-Nov-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15722027	WS-2706	3-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15726314	WS-2706	4-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15726790	WS-2706	4-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15735937	WS-2706	7-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15750089	WS-2706	13-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15757080	WS-2706	17-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15758621	WS-2706	17-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17593205	WS-2706	17-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15768508	WS-2706	26-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17489860	WS-2706	31-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17489861	WS-2706	31-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17489862	WS-2706	31-Dec-2015
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15780922	WS-2706	7-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15788129	WS-2706	9-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15795135	WS-2706	11-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15802898	WS-2706	15-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15825661	WS-2706	24-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15833586	WS-2706	28-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15843113	WS-2706	31-Jan-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15854664	WS-2706	5-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15862600	WS-2706	7-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15876340	WS-2706	13-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15907469	WS-2706	26-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15911364	WS-2706	27-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15913992	WS-2706	28-Feb-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15931040	WS-2706	6-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15945696	WS-2706	12-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15955252	WS-2706	17-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15961295	WS-2706	19-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15968028	WS-2706	21-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15976620	WS-2706	25-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15980375	WS-2706	26-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17510423	WS-2706	31-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17513885	WS-2706	31-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517870	WS-2706	31-Mar-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517873	WS-2706	31-Mar-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	15992404	WS-2706	2-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16000476	WS-2706	4-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16008894	WS-2706	8-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16014993	WS-2706	10-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16025226	WS-2706	15-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16031644	WS-2706	17-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16039471	WS-2706	21-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16043831	WS-2706	22-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16049993	WS-2706	24-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16050156	WS-2706	24-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16053802	WS-2706	25-Apr-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16072487	WS-2706	5-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16078703	WS-2706	7-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16099084	WS-2706	15-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16102961	WS-2706	16-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16116555	WS-2706	22-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16125462	WS-2706	27-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16131916	WS-2706	29-May-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16138324	WS-2706	2-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16144976	WS-2706	4-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16148062	WS-2706	5-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16155839	WS-2706	9-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16158712	WS-2706	10-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16165648	WS-2706	12-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16168783	WS-2706	13-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16176189	WS-2706	17-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16190832	WS-2706	23-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16194169	WS-2706	24-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16196932	WS-2706	25-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17505633	WS-2706	30-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517874	WS-2706	30-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17522736	WS-2706	30-Jun-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16210175	WS-2706	1-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16213053	WS-2706	2-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16218853	WS-2706	7-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16239132	WS-2706	15-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16245794	WS-2706	17-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16253201	WS-2706	21-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16256395	WS-2706	22-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16271046	WS-2706	28-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16277162	WS-2706	30-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17514334	WS-2706	31-Jul-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16289073	WS-2706	5-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16291559	WS-2706	6-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16311107	WS-2706	14-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16324449	WS-2706	20-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16330056	WS-2706	22-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16341152	WS-2706	27-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16342274	WS-2706	27-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17516999	WS-2706	31-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17522729	WS-2706	31-Aug-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16414731	WS-2706	26-Sep-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16450765	WS-2706	10-Oct-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16511476	WS-2706	4-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16524447	WS-2706	10-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16525613	WS-2706	10-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16529939	WS-2706	11-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16531474	WS-2706	12-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16539324	WS-2706	14-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16552003	WS-2706	19-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16559519	WS-2706	21-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16567395	WS-2706	25-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517704	WS-2706	30-Nov-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16576443	WS-2706	1-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16583155	WS-2706	3-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16591524	WS-2706	5-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16591525	WS-2706	5-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16599033	WS-2706	9-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16599034	WS-2706	9-Dec-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16599035	WS-2706	9-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16604972	WS-2706	11-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16604974	WS-2706	11-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16621809	WS-2706	18-Dec-2016
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16642110	WS-2706	7-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16642112	WS-2706	7-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16649393	WS-2706	9-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16676237	WS-2706	20-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16684503	WS-2706	22-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16698601	WS-2706	28-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16701581	WS-2706	29-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517706	WS-2706	31-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17517707	WS-2706	31-Jan-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16709412	WS-2706	2-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16716475	WS-2706	4-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16731674	WS-2706	10-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16738583	WS-2706	12-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16748486	WS-2706	17-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16759988	WS-2706	20-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16768235	WS-2706	24-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16778857	WS-2706	27-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17505628	WS-2706	28-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17523740	WS-2706	28-Feb-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16794222	WS-2706	5-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16799940	WS-2706	6-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16803521	WS-2706	9-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16805224	WS-2706	9-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16807374	WS-2706	10-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16820798	WS-2706	16-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16830259	WS-2706	19-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16838592	WS-2706	23-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16845891	WS-2706	25-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16858005	WS-2706	30-Mar-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16867900	WS-2706	2-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16870754	WS-2706	3-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16890153	WS-2706	10-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16918018	WS-2706	23-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16922573	WS-2706	24-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16930291	WS-2706	28-Apr-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16957199	WS-2706	8-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16970994	WS-2706	14-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16974814	WS-2706	15-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	16990338	WS-2706	21-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17006538	WS-2706	28-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17010685	WS-2706	29-May-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17015193	WS-2706	1-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17023537	WS-2706	4-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17058567	WS-2706	17-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17060697	WS-2706	18-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17072362	WS-2706	23-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17078413	WS-2706	25-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17083017	WS-2706	26-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17086914	WS-2706	29-Jun-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17111598	WS-2706	9-Jul-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17124358	WS-2706	14-Jul-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17148055	WS-2706	23-Jul-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17170093	WS-2706	3-Aug-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17179249	WS-2706	5-Aug-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17246304	WS-2706	1-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17261865	WS-2706	7-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17280718	WS-2706	14-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17288434	WS-2706	16-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17303753	WS-2706	22-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17320842	WS-2706	28-Sep-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17335800	WS-2706	2-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17345077	WS-2706	6-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17355465	WS-2706	9-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17361618	WS-2706	31-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17367565	WS-2706	31-Oct-2017

91033431-FA WATCHMAN® Left Atrial Appendage Closure Device with Delivery System and WATCHMAN®
 Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17385002	WS-2706	31-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17400534	WS-2706	31-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17404702	WS-2706	31-Oct-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17431206	WS-2706	30-Nov-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17440837	WS-2706	30-Nov-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17456606	WS-2706	30-Nov-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17505566	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17513728	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17549051	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17549052	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17557930	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17563433	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17563837	WS-2706	31-Dec-2017
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17742337	WS-2706	15-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17572195	WS-2706	31-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17594047	WS-2706	31-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17611275	WS-2706	31-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17623145	WS-2706	31-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17638865	WS-2706	31-Jan-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17684460	WS-2706	28-Feb-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17687771	WS-2706	28-Feb-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17700105	WS-2706	28-Feb-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17743276	WS-2706	28-Feb-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17750413	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17776215	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17789857	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17795361	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17810105	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17820605	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17829287	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17845078	WS-2706	31-Mar-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17852460	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17861963	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17866075	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17894268	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17910306	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17935421	WS-2706	30-Apr-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17951926	WS-2706	31-May-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	17960378	WS-2706	31-May-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18079952	WS-2706	30-Jun-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18139617	WS-2706	30-Jun-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18142003	WS-2706	30-Jun-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18196863	WS-2706	31-Jul-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18212341	WS-2706	31-Jul-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18220251	WS-2706	31-Jul-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18224620	WS-2706	31-Jul-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18255944	WS-2706	31-Aug-2018
M635WS27060	FG WATCHMAN LAA Closure Device, 27mm,OUS	18284349	WS-2706	31-Aug-2018

Material	Product Description	Batch	Catalog #	Expiration Date
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541702	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541703	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541707	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541708	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541709	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541903	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541904	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541905	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541906	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541907	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541908	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541909	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541913	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541914	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00541916	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544503	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544504	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544505	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544508	WS-2106	31-Aug-2015

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544509	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544511	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544703	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544704	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544706	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544707	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544709	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00544710	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00545201	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547501	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547502	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547505	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547506	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547507	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547508	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547509	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547510	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547511	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547512	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00547514	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548202	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548204	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548205	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548206	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548207	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548210	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548211	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548213	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	F00548215	WS-2106	31-Aug-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15513009	WS-2106	6-Sep-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15556842	WS-2106	24-Sep-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15559607	WS-2106	25-Sep-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15559608	WS-2106	25-Sep-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15573287	WS-2106	1-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15576411	WS-2106	2-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15579547	WS-2106	3-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15583022	WS-2106	4-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15586521	WS-2106	5-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15597796	WS-2106	10-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15608021	WS-2106	15-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15634208	WS-2106	25-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15639752	WS-2106	29-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15645828	WS-2106	31-Oct-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15648179	WS-2106	1-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15650910	WS-2106	1-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15664113	WS-2106	7-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15688025	WS-2106	16-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15694584	WS-2106	20-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15711219	WS-2106	28-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15719185	WS-2106	30-Nov-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15722025	WS-2106	3-Dec-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15746165	WS-2106	12-Dec-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15760115	WS-2106	18-Dec-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15767309	WS-2106	21-Dec-2015
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15791346	WS-2106	10-Jan-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15798829	WS-2106	14-Jan-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15802896	WS-2106	15-Jan-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15802971	WS-2106	15-Jan-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17514543	WS-2106	31-Jan-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15846552	WS-2106	1-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15851185	WS-2106	4-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15859123	WS-2106	6-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15870174	WS-2106	11-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15883167	WS-2106	15-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15889808	WS-2106	19-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15906797	WS-2106	26-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15913993	WS-2106	28-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15918189	WS-2106	29-Feb-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17505941	WS-2106	29-Feb-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15921282	WS-2106	3-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15927421	WS-2106	5-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15931041	WS-2106	6-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15934213	WS-2106	7-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15935418	WS-2106	7-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15942659	WS-2106	11-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15955253	WS-2106	17-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15957982	WS-2106	18-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15965734	WS-2106	20-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15972756	WS-2106	24-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15980369	WS-2106	26-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17514542	WS-2106	31-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17516996	WS-2106	31-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17517162	WS-2106	31-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17517166	WS-2106	31-Mar-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15989303	WS-2106	1-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15990230	WS-2106	1-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15992405	WS-2106	2-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	15996022	WS-2106	3-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16018190	WS-2106	11-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16022222	WS-2106	14-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16039472	WS-2106	21-Apr-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16063043	WS-2106	1-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16069159	WS-2106	2-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16081852	WS-2106	8-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16095359	WS-2106	14-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16102740	WS-2106	16-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16125461	WS-2106	27-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16135272	WS-2106	30-May-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16138323	WS-2106	2-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16151320	WS-2106	6-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16168784	WS-2106	13-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16172940	WS-2106	16-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16180118	WS-2106	18-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16190833	WS-2106	23-Jun-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16208087	WS-2106	1-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16228219	WS-2106	10-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16239134	WS-2106	15-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16242466	WS-2106	16-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16245796	WS-2106	17-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16271044	WS-2106	28-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16274328	WS-2106	29-Jul-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16284579	WS-2106	4-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16297407	WS-2106	8-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16326970	WS-2106	21-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16330055	WS-2106	22-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16337479	WS-2106	26-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16342271	WS-2106	27-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17505943	WS-2106	31-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17522735	WS-2106	31-Aug-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16410910	WS-2106	25-Sep-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16457878	WS-2106	14-Oct-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16459272	WS-2106	14-Oct-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16511473	WS-2106	4-Nov-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16528429	WS-2106	11-Nov-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16535544	WS-2106	13-Nov-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16547561	WS-2106	18-Nov-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16566477	WS-2106	25-Nov-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16586724	WS-2106	4-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16601179	WS-2106	10-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16605567	WS-2106	11-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16611630	WS-2106	15-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16615630	WS-2106	16-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16626168	WS-2106	25-Dec-2016
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16669688	WS-2106	16-Jan-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16676235	WS-2106	20-Jan-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16684501	WS-2106	22-Jan-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16691147	WS-2106	26-Jan-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16706057	WS-2106	30-Jan-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17505944	WS-2106	31-Jan-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16709410	WS-2106	2-Feb-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16714136	WS-2106	3-Feb-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16719785	WS-2106	5-Feb-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16743850	WS-2106	13-Feb-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16772082	WS-2106	25-Feb-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16783761	WS-2106	2-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16794223	WS-2106	5-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16799658	WS-2106	6-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16807375	WS-2106	10-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16811004	WS-2106	11-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16812948	WS-2106	12-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16820799	WS-2106	16-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16826393	WS-2106	18-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16835331	WS-2106	20-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16858003	WS-2106	30-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16860706	WS-2106	31-Mar-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16875092	WS-2106	6-Apr-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16919733	WS-2106	23-Apr-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16927173	WS-2106	27-Apr-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16931941	WS-2106	28-Apr-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16946166	WS-2106	5-May-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16957819	WS-2106	8-May-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	16984103	WS-2106	19-May-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17015190	WS-2106	1-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17016764	WS-2106	2-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17031819	WS-2106	8-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17053901	WS-2106	16-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17088542	WS-2106	29-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17588070	WS-2106	29-Jun-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17158555	WS-2106	28-Jul-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17161831	WS-2106	29-Jul-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17162162	WS-2106	29-Jul-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17162163	WS-2106	29-Jul-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17162164	WS-2106	29-Jul-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17174212	WS-2106	4-Aug-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17190029	WS-2106	10-Aug-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17193307	WS-2106	11-Aug-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17256409	WS-2106	4-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17265820	WS-2106	8-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17271026	WS-2106	9-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17274694	WS-2106	10-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17281586	WS-2106	14-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17282611	WS-2106	14-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17308887	WS-2106	23-Sep-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17331530	WS-2106	1-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17337336	WS-2106	2-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17368985	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17380187	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17384737	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17384738	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17384739	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17385000	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17389788	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17400531	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17404701	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17408444	WS-2106	31-Oct-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17421244	WS-2106	30-Nov-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17424065	WS-2106	30-Nov-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17428086	WS-2106	30-Nov-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17442611	WS-2106	30-Nov-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17462354	WS-2106	30-Nov-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17544822	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17553573	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17563836	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17566509	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17566874	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17566875	WS-2106	31-Dec-2017
M635WS21060	FG WATCHMAN LAA Closure Device, 21mm,OUS	17566876	WS-2106	31-Dec-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553707	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553708	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553709	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553710	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553711	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553712	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553713	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553714	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553715	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553716	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	F00553717	TS-2006	31-Aug-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15523924	TS-2006	11-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15534295	TS-2006	14-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15547790	TS-2006	20-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15547791	TS-2006	20-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15547792	TS-2006	20-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15549844	TS-2006	20-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15557261	TS-2006	24-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15557262	TS-2006	24-Sep-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15569373	TS-2006	1-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15573927	TS-2006	1-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15558971	TS-2006	2-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15577456	TS-2006	2-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15577564	TS-2006	2-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15577565	TS-2006	2-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15582831	TS-2006	4-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15591639	TS-2006	8-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15596436	TS-2006	10-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15597413	TS-2006	10-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15604737	TS-2006	12-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15607853	TS-2006	15-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15610775	TS-2006	16-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15614197	TS-2006	17-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15614198	TS-2006	17-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15630136	TS-2006	24-Oct-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15624240	TS-2006	1-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15652977	TS-2006	2-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15640384	TS-2006	5-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15657830	TS-2006	5-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15660313	TS-2006	6-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15640385	TS-2006	7-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15647162	TS-2006	9-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15630139	TS-2006	12-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15675108	TS-2006	12-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15675110	TS-2006	12-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15676735	TS-2006	12-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15647163	TS-2006	13-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15679745	TS-2006	13-Nov-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15679746	TS-2006	3-Dec-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15691631	TS-2006	5-Dec-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15684861	TS-2006	6-Dec-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15691634	TS-2006	11-Dec-2015
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15775097	TS-2006	3-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15778019	TS-2006	4-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15780647	TS-2006	7-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15788130	TS-2006	9-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15791345	TS-2006	10-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15795851	TS-2006	11-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15799915	TS-2006	14-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15802493	TS-2006	15-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15818893	TS-2006	22-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15828834	TS-2006	25-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15837151	TS-2006	29-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524411	TS-2006	31-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524416	TS-2006	31-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524417	TS-2006	31-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524418	TS-2006	31-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524720	TS-2006	31-Jan-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524721	TS-2006	31-Jan-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15846553	TS-2006	1-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15851186	TS-2006	4-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15854603	TS-2006	5-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15858077	TS-2006	6-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15862218	TS-2006	7-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15866281	TS-2006	8-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15870175	TS-2006	11-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15876700	TS-2006	13-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15883165	TS-2006	15-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15889807	TS-2006	19-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15892273	TS-2006	20-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15895940	TS-2006	21-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15899488	TS-2006	22-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15904004	TS-2006	25-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15910338	TS-2006	27-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15913991	TS-2006	28-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524722	TS-2006	29-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524723	TS-2006	29-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524724	TS-2006	29-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524725	TS-2006	29-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524726	TS-2006	29-Feb-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15917331	TS-2006	1-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15921281	TS-2006	3-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15923967	TS-2006	4-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15931042	TS-2006	6-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15934214	TS-2006	7-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15938318	TS-2006	10-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15945693	TS-2006	12-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15955251	TS-2006	17-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15965197	TS-2006	20-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15968027	TS-2006	21-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15972749	TS-2006	24-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15976621	TS-2006	25-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15980066	TS-2006	26-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15982571	TS-2006	27-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524727	TS-2006	31-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524729	TS-2006	31-Mar-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15986987	TS-2006	1-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15989302	TS-2006	1-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15992406	TS-2006	2-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15992934	TS-2006	2-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	15997550	TS-2006	3-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16005081	TS-2006	7-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16008695	TS-2006	8-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16015248	TS-2006	10-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16018005	TS-2006	11-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16022224	TS-2006	14-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16025224	TS-2006	15-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16031641	TS-2006	17-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16039474	TS-2006	21-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16042993	TS-2006	22-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16045736	TS-2006	23-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16049761	TS-2006	24-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16052878	TS-2006	25-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524731	TS-2006	30-Apr-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16056581	TS-2006	1-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16059512	TS-2006	2-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16063042	TS-2006	3-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16072486	TS-2006	5-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16075309	TS-2006	6-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16075311	TS-2006	6-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16078702	TS-2006	7-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16085500	TS-2006	9-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16089346	TS-2006	12-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16092734	TS-2006	13-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16095740	TS-2006	14-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16099082	TS-2006	15-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16102078	TS-2006	16-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16102676	TS-2006	16-May-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16107064	TS-2006	19-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16110480	TS-2006	20-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16113541	TS-2006	21-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16116554	TS-2006	22-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16119788	TS-2006	23-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16125460	TS-2006	27-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16128359	TS-2006	28-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16131915	TS-2006	29-May-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16138326	TS-2006	2-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16141222	TS-2006	3-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16144572	TS-2006	4-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16145326	TS-2006	4-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16150818	TS-2006	6-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16155310	TS-2006	9-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16158711	TS-2006	10-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16161518	TS-2006	11-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16168782	TS-2006	13-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16176187	TS-2006	17-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16183060	TS-2006	19-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16186299	TS-2006	20-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16193942	TS-2006	24-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16196930	TS-2006	25-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16200326	TS-2006	26-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16204070	TS-2006	27-Jun-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16207559	TS-2006	1-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16210178	TS-2006	1-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16213052	TS-2006	2-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16218852	TS-2006	7-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16222332	TS-2006	8-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16225384	TS-2006	9-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16228241	TS-2006	10-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16231880	TS-2006	11-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16235347	TS-2006	14-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16239221	TS-2006	15-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16245797	TS-2006	17-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16249323	TS-2006	18-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16253202	TS-2006	21-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16256863	TS-2006	22-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16262989	TS-2006	24-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16266635	TS-2006	25-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16270213	TS-2006	28-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16273633	TS-2006	29-Jul-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16282152	TS-2006	1-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16284577	TS-2006	4-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16287159	TS-2006	5-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16289074	TS-2006	5-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16294255	TS-2006	7-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16297274	TS-2006	8-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16300507	TS-2006	11-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16311011	TS-2006	14-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16311108	TS-2006	14-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16316943	TS-2006	18-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16319992	TS-2006	19-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16323506	TS-2006	20-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16326972	TS-2006	21-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16330057	TS-2006	22-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16334237	TS-2006	25-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16338001	TS-2006	26-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16341153	TS-2006	27-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524732	TS-2006	31-Aug-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16344102	TS-2006	3-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16347281	TS-2006	3-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16410911	TS-2006	25-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16418478	TS-2006	29-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16422805	TS-2006	30-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16423443	TS-2006	30-Sep-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16455414	TS-2006	13-Oct-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16483226	TS-2006	23-Oct-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16519540	TS-2006	6-Nov-2016

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16519541	TS-2006	6-Nov-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16547565	TS-2006	18-Nov-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16548609	TS-2006	18-Nov-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524733	TS-2006	30-Nov-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16579911	TS-2006	2-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16583156	TS-2006	3-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16590621	TS-2006	5-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16598142	TS-2006	9-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16601174	TS-2006	10-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16602910	TS-2006	10-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16608075	TS-2006	12-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16611632	TS-2006	15-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16615627	TS-2006	16-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16619147	TS-2006	17-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16621806	TS-2006	18-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16623957	TS-2006	19-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16626169	TS-2006	25-Dec-2016
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16629592	TS-2006	1-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16638766	TS-2006	6-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16642109	TS-2006	7-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16642111	TS-2006	7-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16649392	TS-2006	9-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16654987	TS-2006	12-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16655276	TS-2006	12-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16655277	TS-2006	12-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16655278	TS-2006	12-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16678438	TS-2006	21-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16679009	TS-2006	21-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16682438	TS-2006	22-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16695334	TS-2006	27-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16698285	TS-2006	28-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16701585	TS-2006	29-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16702550	TS-2006	29-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524734	TS-2006	31-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524735	TS-2006	31-Jan-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16709414	TS-2006	2-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16713086	TS-2006	3-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16716476	TS-2006	4-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16719767	TS-2006	5-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16719776	TS-2006	5-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16727778	TS-2006	9-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16735188	TS-2006	11-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16738581	TS-2006	12-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16748485	TS-2006	17-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16752079	TS-2006	18-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16768234	TS-2006	24-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16772083	TS-2006	25-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16775461	TS-2006	26-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17524736	TS-2006	28-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17538718	TS-2006	28-Feb-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16787947	TS-2006	3-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16790633	TS-2006	4-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16799406	TS-2006	6-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16803525	TS-2006	9-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16805226	TS-2006	9-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16807376	TS-2006	10-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16812947	TS-2006	12-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16816972	TS-2006	13-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16820960	TS-2006	16-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16823548	TS-2006	17-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16826394	TS-2006	18-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16830360	TS-2006	19-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16833974	TS-2006	20-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16849203	TS-2006	26-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16850145	TS-2006	27-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16852559	TS-2006	27-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16857214	TS-2006	30-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16860704	TS-2006	31-Mar-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16863193	TS-2006	1-Apr-2017

Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16865854	TS-2006	1-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16867902	TS-2006	2-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16870755	TS-2006	3-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16875093	TS-2006	6-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16875734	TS-2006	6-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16900812	TS-2006	15-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16911587	TS-2006	21-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16915195	TS-2006	22-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16918019	TS-2006	23-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16922572	TS-2006	24-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16927174	TS-2006	27-Apr-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16946167	TS-2006	5-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16949828	TS-2006	6-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16953116	TS-2006	7-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16957400	TS-2006	8-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16957818	TS-2006	8-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16964405	TS-2006	12-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16967574	TS-2006	13-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16970995	TS-2006	14-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16974815	TS-2006	15-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16979304	TS-2006	18-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16984104	TS-2006	19-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16990814	TS-2006	21-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16994018	TS-2006	22-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	16999911	TS-2006	26-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17003660	TS-2006	27-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17006537	TS-2006	28-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17010728	TS-2006	29-May-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17013673	TS-2006	1-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17020297	TS-2006	3-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17023538	TS-2006	4-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17027982	TS-2006	5-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17032220	TS-2006	8-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17035414	TS-2006	9-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17046073	TS-2006	12-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17053900	TS-2006	16-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17056855	TS-2006	17-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17058566	TS-2006	17-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17064562	TS-2006	19-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17068255	TS-2006	22-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17071798	TS-2006	23-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17075340	TS-2006	24-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17078716	TS-2006	25-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17083015	TS-2006	26-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17086912	TS-2006	29-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17088543	TS-2006	29-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17091394	TS-2006	30-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17092084	TS-2006	30-Jun-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17101486	TS-2006	6-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17103308	TS-2006	6-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17107951	TS-2006	8-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17111599	TS-2006	9-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17115175	TS-2006	10-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17118418	TS-2006	13-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17142698	TS-2006	21-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17148056	TS-2006	23-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17152068	TS-2006	24-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17155707	TS-2006	27-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17158556	TS-2006	28-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17161832	TS-2006	29-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17164904	TS-2006	30-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17167875	TS-2006	31-Jul-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17174211	TS-2006	4-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17175418	TS-2006	4-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17179248	TS-2006	5-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17182780	TS-2006	6-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17185582	TS-2006	7-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17190030	TS-2006	10-Aug-2017
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	17193308	TS-2006	11-Aug-2017

91033431-FA WATCHMAN® Left Atrial Appendage Closure Device with Delivery System and WATCHMAN®
 Access System Access Sheath with Dilator
List of affected UPN and Batch numbers

M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036461	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036462	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036463	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036467	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036468	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036469	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036470	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036471	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18036472	TS-2006	31-May-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18075936	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18075937	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18075938	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18110489	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18115034	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18115035	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18120793	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18154221	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18154222	TS-2006	30-Jun-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18158642	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18158643	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18158644	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18158645	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18158646	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18204079	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18204400	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18204401	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233257	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233259	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233260	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233261	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233262	TS-2006	31-Jul-2018
M635TS20060	FG WATCHMAN Access Sys Double Curve,OUS	18233263	TS-2006	31-Jul-2018

Material	Product Description	Batch	Catalog #	Expiration Date
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15611295	TS-4006	16-Oct-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15624432	TS-4006	22-Oct-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15653185	TS-4006	2-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15657827	TS-4006	5-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15675111	TS-4006	12-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15640386	TS-4006	13-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15679748	TS-4006	13-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15684667	TS-4006	15-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15698742	TS-4006	21-Nov-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15716330	TS-4006	4-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15736494	TS-4006	7-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15736689	TS-4006	7-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15746782	TS-4006	12-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15711284	TS-4006	14-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15753133	TS-4006	14-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15757057	TS-4006	17-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15760114	TS-4006	18-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15764019	TS-4006	19-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15767450	TS-4006	21-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15768507	TS-4006	26-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15770114	TS-4006	28-Dec-2015
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	15772280	TS-4006	2-Jan-2016
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	16635566	TS-4006	5-Jan-2017
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	16691146	TS-4006	26-Jan-2017
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	17148680	TS-4006	23-Jul-2017
M635TS40060	FG WATCHMAN Access Sys Anterior Curv,OUS	17628922	TS-4006	31-Jan-2018