



Fannin Ltd.  
Fannin House, South County Business Park,  
Leopardstown, Dublin 18.  
Tel: + 353 (01) 2907000  
Fax: + 353 (01) 2907118  
www.fannin.ie



Fannin Ltd.  
I.D.A. Industrial Estate,  
Ballybrit Upper, Galway.  
Tel: + 353 (01) 2907204  
Fax: + 353 (091) 771187  
www.lipdiagnostic.com

### **Urgent Field Safety Notice**

Date: 24/08/2018

W11326 Staph/Strep Selective Agar Modified

1282892533 Exp. 24/08/2018

1283904725 Exp. 01/09/2018

1284918133 Exp. 09/09/2018

W11320 Fastidious Anaerobic Agar + Nav

0361900928 Exp. 03/09/2018

An internal investigation has identified that the lots listed above may not provide the expected microbiological performance and the growth of non-target organisms could mask the recovery of target organisms.

Internal QC testing was performed on the lots listed above prior to release and results were within specification for productivity ratio and colony morphology.

Requirement for review of reported test results should be determined by the appropriate technical expert. We request that any remaining stock of the lots listed should be destroyed.

This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please maintain awareness of this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

**The HPRA and MHRA have been notified**

A handwritten signature in blue ink, appearing to read "Irene Slevin".

Irene Slevin QA Manager  
Irene.slevin@fannin.eu