



Fannin Ltd.  
Fannin House, South County Business Park,  
Leopardstown, Dublin 18.  
Tel: + 353 (01) 2907000  
Fax: + 353 (01) 2907118  
www.fannin.ie



c/o Fannin Ltd.  
Fannin House, South County Business Park,  
Leopardstown, Dublin 18.  
Tel: + 353 (01) 2907000  
Fax: + 353 (01) 2907118  
www.oxid.com



Fannin Ltd.  
I.D.A. Ballybrit Upper Ind Est,  
Ballybrit, Galway.  
Tel: + 353 (01) 2907204  
Fax: + 353 (091) 771187  
www.lipdiagnostic.com

## Urgent Field Safety Notice – Product Recall Notification

08/10/15

**W11024 BURKHOLDERIA CEPACIA SPECIFIC AGAR**

**Lot no 0962079323**

Customers are advised that upon review of the batch production record a query has arisen in relation to the addition of the selective supplement Polymixin B.

There is a risk that the supplement was prepared incorrectly resulting in reduced concentration of the selective agent in the prepared media plates. Polymixin B is included in the formulation to suppress Gram negative back-round flora.

A root cause investigation has commenced under our non-conformance procedure.

A replacement batch is now in production and will be QC tested over the coming days.

All end-users should be made aware of this notification. This lot number, delivered by courier on the morning of 08/10/15 should be immediately disposed of to prevent use.

The HPRA have been advised of this notification.

Please also refer to FaxBack Form supplied – replacements or credit will be issued as indicated by completion and return of the form.

Please accept our sincere apologies any inconvenience cause, should you require any further information please do not hesitate to contact me directly.

Best Regards,

---

Eva Lyons – QA Manager (LIP)  
eva.lyons@fannin.eu