

Martillac, October 2nd 2019

Subject : Security information n°R1916718 ANSM

Dear Madam, Dear Sir,

We hereby confirm that the Giemsa Solution lot I88021 and lot I97523 has a staining defect which could appear depending on the method used and result in stainings that may be either too dark or too pale.

These lots of Giemsa have been sold in 1 Litre bottles under the reference 75030SX1000 and are being recalled.

Therefore, all bottles that you have in stock must be withdrawn and destroyed.


You will find below a destruction certificate to be filled in and sent to your distributor. The products will be replaced as soon as possible.

The ANSM (French National Authority) has been informed of this communication.

Be assured of our total commitment in solving this problem and in the satisfaction of our clients.

Julien VEYSSY

COO



Date :

Made by :

Society / department :

I declare, to have destroyed the following products :

PRODUCTS	REFERENCES	QUANTITIES	LOTS

Signature and company stamp.

Document to be returned to your Sysmex Distributor