

«Hospital_Name»
 «Users_Name»- «Department»
 «Customer_Address»
 «Zip_Code» «City» - «Country_name»

Reference: 91107242-FA

27 November 2015

Urgent Field Safety Notice - Medical Device Removal Chariot™ Guiding Sheath

Dear «Users_Name»,

Boston Scientific is voluntarily recalling its Chariot™ Guiding Sheath. To date, Boston Scientific has received fourteen complaints for shaft separation, four for distal shaft separation, which occurred during device preparation or use. The most severe outcome of this failure is embolism of device fragments, which could lead to obstruction of blood flow, or additional intervention to remove a device fragment.

Our records indicate that your facility received some of the concerned product. **Table below provides a complete list of all affected products**, including Product Description, Material Number (UPN) and Lot/Batch numbers and expiry date. Please note that **only the material listed below is affected. No other Boston Scientific product is involved by this Field Safety Notice.**

Further distribution or use of any remaining product affected by this action should cease immediately.

Product Description	Material Number (UPN)	Batch	Batch Expiration Date Range
Chariot™ Guiding Sheath	H74939277645110	See Attached Affected Product Listing	31 March 2016 - 30 November 2018
	H74939277690210		
	H74939277690220		
	H74939277845210		
	H74939277745210		
	H74939277645210		
	H74939277745110		
	H74939277665110		
	H74939277690110		
	H74939277765110		
	H74939277545110		
	H74939277690120		
	H74939277845110		
	H74939277790110		
	H74939277865110		
	H74939277545210		
	H74939277790220		
	H74939277665120		
	H74939277645220		
	H74939277645120		
	H74939277745120		
	H74939277790210		
	H74939277790120		
	H74939277890110		
H74939277765120			
H74939277865120			
H74939277890120			
H74939277845120			
H74939277745220			

INSTRUCTIONS:

- 1- **Please immediately discontinue use of the Boston Scientific product reported in the list and remove all of the affected units from your inventory**, irregardless of where these units are stored in your facility. Segregate the units in a secure place, pending return to Boston Scientific.
- 2- **Please complete the attached Verification Form even if you do not have any product to return.**
- 3- **When completed, please return the Verification Form to your local Boston Scientific office** to the attention of Cas Sherrard - Fax: 01442 411816 - Email: UK-Quality@bsci.com; on or before **11 December 2015**.
- 4- **If you have products to return**, please package them in an appropriate shipping box and **contact** Cas Sherrard - Tel: 01442 411731 **of your local Boston Scientific office**, to arrange return.
- 5- Please pass this notice to any health professional of your organization that need to be aware and to any organization where the potentially affected devices have been transferred (If appropriate). Please provide Boston Scientific with details of any affected devices that have been transferred to other organizations (if appropriate).

Your Competent Authority is being notified of this Field Safety Notice.

We regret any inconvenience that this action may cause, and we appreciate your understanding as we take action to ensure patient safety and customer satisfaction.

If you have any questions or would like assistance with this Field Safety Notice, please contact your local Sales Representative.

Yours sincerely,



Marie Pierre Barlanga
Quality Department
Boston Scientific International S.A.

Attachments: - Verification Form
- Full affected Product listing

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797683	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797684	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797685	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797686	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797687	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797688	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797689	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797690	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797691	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797692	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797693	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797694	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797695	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797696	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17797697	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798268	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798269	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798270	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798271	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798272	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17798273	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17806530	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17829972	08714729863861	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17829976	08714729863861	31-Mar-16
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	17830219	08714729863946	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17839346	08714729863861	31-Mar-16
H74939277690220	Guiding Sheath 6F, 90 cm, MP, TB	17839347	08714729863953	31-Mar-16
H74939277845210	Guiding Sheath 8F, 45 cm, MP, CC	17839350	08714729864080	31-Mar-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17866538	08714729863861	30-Apr-16
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	17866881	08714729863984	30-Apr-16
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	17866882	08714729863984	30-Apr-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17872044	08714729863861	30-Apr-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17872045	08714729863861	30-Apr-16
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	17915368	08714729863885	30-Apr-16
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	17915369	08714729863984	30-Apr-16
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17940916	08714729863861	30-Apr-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17941111	08714729863861	30-Apr-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17943786	08714729863960	30-Apr-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17945064	08714729863960	30-Apr-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	17945224	08714729863908	30-Apr-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	17946351	08714729863922	30-Apr-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17949277	08714729863861	30-Apr-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17950404	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17951572	08714729863960	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17952749	08714729863960	31-May-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	17953534	08714729864004	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17958093	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17958489	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17961473	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17962740	08714729863861	31-May-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	17966591	08714729863847	31-May-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	17971123	08714729863847	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17972092	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17972761	08714729863861	31-May-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	17976335	08714729863939	31-May-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	17977575	08714729863908	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17978575	08714729863960	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17981064	08714729863960	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17983006	08714729863861	31-May-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	17983027	08714729864066	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17984693	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17985812	08714729863960	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17986420	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17987381	08714729863861	31-May-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	17989008	08714729864028	31-May-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	17990047	08714729863922	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17991170	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	17991321	08714729863861	31-May-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	17994529	08714729864004	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17995087	08714729863960	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17996706	08714729863960	31-May-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	17997093	08714729863885	31-May-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	17997095	08714729863885	31-May-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	17997242	08714729863908	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	17999337	08714729863960	31-May-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	17999936	08714729863922	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18002533	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18003369	08714729863960	31-May-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18003765	08714729863984	31-May-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18003862	08714729863847	31-May-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18003942	08714729864066	31-May-18
H74939277865110	Guiding Sheath 8F, 65 cm, ST, CC	18005013	08714729864097	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18007740	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18007952	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18008205	08714729863960	31-May-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18008324	08714729863984	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18010176	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18012336	08714729863861	31-May-18
H74939277545210	Guiding Sheath 5F, 45 cm, MP, CC	18012757	08714729863854	31-May-18
H74939277845210	Guiding Sheath 8F, 45 cm, MP, CC	18013142	08714729864080	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18015270	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18015790	08714729863861	31-May-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18016778	08714729863908	31-May-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18017840	08714729863939	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18020685	08714729863861	31-May-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18029579	08714729863885	31-May-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18029913	08714729863885	31-May-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18034289	08714729863885	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18034502	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18034664	08714729863960	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18036489	08714729863960	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18038122	08714729863861	31-May-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18044434	08714729863847	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18044488	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18044489	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18044490	08714729863960	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18044491	08714729863960	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18046224	08714729863861	31-May-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18046901	08714729863861	31-May-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18047347	08714729863960	31-May-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18047887	08714729863847	31-May-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18048032	08714729864066	31-May-18
H74939277790220	Guiding Sheath 7F, 90 cm, MP, TB	18049469	08714729864059	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18051103	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18051829	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18052052	08714729863960	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18052409	08714729863908	30-Jun-18
H74939277665120	Guiding Sheath 6F, 65 cm, ST, TB	18053044	08714729863915	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18056432	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18057517	08714729863861	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18058857	08714729863908	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18059057	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18059058	08714729863960	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18062333	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18062608	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18065837	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18066670	08714729863861	30-Jun-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18067077	08714729863984	30-Jun-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18067146	08714729863847	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18070015	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18072456	08714729863861	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18073066	08714729863922	30-Jun-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18073470	08714729863984	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18074999	08714729863922	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18076430	08714729863885	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18077317	08714729863885	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18078017	08714729863960	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18078082	08714729863922	30-Jun-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18078165	08714729863939	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18079822	08714729863960	30-Jun-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18082142	08714729863939	30-Jun-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18082207	08714729863984	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18082531	08714729863861	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18082921	08714729863885	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18085825	08714729863861	30-Jun-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18086686	08714729864028	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18086882	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18086932	08714729863960	30-Jun-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18092005	08714729863984	30-Jun-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18092006	08714729864004	30-Jun-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18093994	08714729863847	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18094930	08714729863861	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18096358	08714729863885	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18097378	08714729863861	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18098350	08714729863908	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18098815	08714729863908	30-Jun-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18099638	08714729864004	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18101449	08714729863922	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18102355	08714729863922	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18102910	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18103367	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18104025	08714729863960	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18104630	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18105040	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18107059	08714729863960	30-Jun-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18107083	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18107524	08714729863960	30-Jun-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18107757	08714729864028	30-Jun-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18108782	08714729864066	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18110405	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18110948	08714729863861	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18111971	08714729863908	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18112662	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18113138	08714729863960	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18115318	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18115964	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18116708	08714729863960	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18116885	08714729863922	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18116989	08714729863908	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18118304	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18118931	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18122227	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18122237	08714729863861	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18122635	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18122959	08714729863960	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18126319	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18126915	08714729863861	30-Jun-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18127339	08714729864004	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18127733	08714729863960	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18127962	08714729863960	30-Jun-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18128758	08714729863922	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18129548	08714729863861	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18130475	08714729863885	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18131052	08714729863885	30-Jun-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18131652	08714729863885	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18131911	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18132320	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18132945	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18133938	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18134577	08714729863861	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18135679	08714729863908	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18135986	08714729863908	30-Jun-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18135987	08714729864004	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18136779	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18137430	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18138447	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18139476	08714729863861	30-Jun-18
H74939277645220	Guiding Sheath 6F, 45 cm, MP, TB	18139925	08714729863892	30-Jun-18
H74939277645120	Guiding Sheath 6F, 45 cm, ST, TB	18140726	08714729863878	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18141215	08714729863861	30-Jun-18
H74939277745120	Guiding Sheath 7F, 45 cm, ST, TB	18141391	08714729863977	30-Jun-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18141405	08714729864066	30-Jun-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18141915	08714729864066	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18142142	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18142173	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18143369	08714729863861	30-Jun-18
H74939277790210	Guiding Sheath 7F, 90 cm, MP, CC	18144790	08714729864042	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18145669	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18145832	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18145833	08714729863861	30-Jun-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18146076	08714729864066	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18147635	08714729863960	30-Jun-18
H74939277790120	Guiding Sheath 7F, 90 cm, ST, TB	18147874	08714729864035	30-Jun-18
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	18149313	08714729863946	30-Jun-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18150210	08714729864066	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18150467	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18150468	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18150469	08714729863861	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18150470	08714729863861	30-Jun-18
H74939277690220	Guiding Sheath 6F, 90 cm, MP, TB	18151082	08714729863953	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18151355	08714729863960	30-Jun-18
H74939277890110	Guiding Sheath 8F, 90 cm, ST, CC	18151492	08714729864110	30-Jun-18
H74939277765120	Guiding Sheath 7F, 65 cm, ST, TB	18152060	08714729864011	30-Jun-18
H74939277845210	Guiding Sheath 8F, 45 cm, MP, CC	18155063	08714729864080	30-Jun-18
H74939277865120	Guiding Sheath 8F, 65 cm, ST, TB	18155965	08714729864103	30-Jun-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18156068	08714729863960	30-Jun-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18156069	08714729863908	30-Jun-18
H74939277890120	Guiding Sheath 8F, 90 cm, ST, TB	18156127	08714729864127	30-Jun-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18164527	08714729863861	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18166595	08714729863922	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18166935	08714729863922	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18175318	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18175368	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18175369	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18175449	08714729863861	31-Jul-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18175629	08714729863984	31-Jul-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18175633	08714729863984	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18176217	08714729863861	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18176854	08714729863922	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18177736	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18183385	08714729863861	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18184155	08714729863960	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18184219	08714729863861	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18184332	08714729863922	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18184430	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18186104	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18186317	08714729863960	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18187801	08714729863885	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18188328	08714729863885	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18189087	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18189088	08714729863861	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18189435	08714729863885	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18189886	08714729863885	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18190437	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18190498	08714729863960	31-Jul-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18195307	08714729863908	31-Jul-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18195524	08714729863908	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18196366	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18197790	08714729863960	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18198875	08714729863885	31-Jul-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18199716	08714729863984	31-Jul-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18199826	08714729863847	31-Jul-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18199853	08714729863908	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18199856	08714729863960	31-Jul-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18203679	08714729863847	31-Jul-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18210498	08714729863908	31-Jul-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18211954	08714729864028	31-Jul-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18213187	08714729863847	31-Jul-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18213425	08714729864066	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18213426	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18213429	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18213430	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18213431	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18214135	08714729863960	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221960	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221961	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221962	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221963	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221964	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221965	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221966	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221967	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221971	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221972	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221973	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18221974	08714729863861	31-Jul-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18223105	08714729864004	31-Jul-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18223233	08714729864004	31-Jul-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18226247	08714729864004	31-Jul-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18227249	08714729863908	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18227827	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18228071	08714729863960	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18230234	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18230892	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18231027	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18231034	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18231361	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18231908	08714729863861	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18232150	08714729863960	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18232319	08714729863922	31-Jul-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18232485	08714729863922	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18233549	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18234386	08714729863861	31-Jul-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18235819	08714729863861	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18236889	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18237097	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18237634	08714729863960	31-Jul-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18237659	08714729863960	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18238199	08714729863885	31-Jul-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18239959	08714729863885	31-Jul-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18241213	08714729863939	31-Jul-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18241447	08714729863847	31-Jul-18
H74939277745120	Guiding Sheath 7F, 45 cm, ST, TB	18242208	08714729863977	31-Jul-18
H74939277645220	Guiding Sheath 6F, 45 cm, MP, TB	18243084	08714729863892	31-Jul-18
H74939277790210	Guiding Sheath 7F, 90 cm, MP, CC	18266289	08714729864042	31-Aug-18
H74939277790210	Guiding Sheath 7F, 90 cm, MP, CC	18266831	08714729864042	31-Aug-18
H74939277845120	Guiding Sheath 8F, 45 cm, ST, TB	18267278	08714729864073	31-Aug-18
H74939277745220	Guiding Sheath 7F, 45 cm, MP, TB	18267593	08714729863991	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18267975	08714729863847	31-Aug-18
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	18271147	08714729863946	31-Aug-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277790120	Guiding Sheath 7F, 90 cm, ST, TB	18272347	08714729864035	31-Aug-18
H74939277890110	Guiding Sheath 8F, 90 cm, ST, CC	18275111	08714729864110	31-Aug-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18276555	08714729863922	31-Aug-18
H74939277865110	Guiding Sheath 8F, 65 cm, ST, CC	18276556	08714729864097	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18276596	08714729863847	31-Aug-18
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	18276597	08714729863946	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18282714	08714729863861	31-Aug-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18282966	08714729863922	31-Aug-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18285339	08714729863885	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18286589	08714729863847	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18286953	08714729863847	31-Aug-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18287102	08714729863922	31-Aug-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18290794	08714729863922	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18305021	08714729863847	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18309633	08714729864004	31-Aug-18
H74939277765120	Guiding Sheath 7F, 65 cm, ST, TB	18310068	08714729864011	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18310205	08714729864004	31-Aug-18
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	18314233	08714729863946	31-Aug-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18318148	08714729864028	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18319481	08714729864004	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18323631	08714729864004	31-Aug-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18323931	08714729863984	31-Aug-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18323954	08714729863984	31-Aug-18
H74939277790120	Guiding Sheath 7F, 90 cm, ST, TB	18324899	08714729864035	31-Aug-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18324915	08714729864028	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18325121	08714729863847	31-Aug-18
H74939277790210	Guiding Sheath 7F, 90 cm, MP, CC	18325669	08714729864042	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18325680	08714729863908	31-Aug-18
H74939277645120	Guiding Sheath 6F, 45 cm, ST, TB	18326280	08714729863878	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18329813	08714729863847	31-Aug-18
H74939277790210	Guiding Sheath 7F, 90 cm, MP, CC	18329858	08714729864042	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18330833	08714729863861	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18331239	08714729863861	31-Aug-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18333673	08714729863847	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18334834	08714729863861	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18335450	08714729863861	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18335511	08714729863861	31-Aug-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18335976	08714729863861	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18339203	08714729864004	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18339412	08714729864004	31-Aug-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18339474	08714729864004	31-Aug-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18343686	08714729864028	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18344722	08714729863908	31-Aug-18
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18345214	08714729863939	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18350415	08714729863908	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18351103	08714729863908	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18352058	08714729863960	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18354607	08714729863960	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18354703	08714729863908	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18354978	08714729863960	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18355365	08714729863960	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18355470	08714729863960	31-Aug-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18358197	08714729863908	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18358876	08714729863960	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18359047	08714729863960	31-Aug-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18359609	08714729863960	31-Aug-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18361656	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18363042	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18363407	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18364153	08714729863960	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18364154	08714729863908	30-Sep-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18366116	08714729863847	30-Sep-18
H74939277665120	Guiding Sheath 6F, 65 cm, ST, TB	18368060	08714729863915	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18368264	08714729863908	30-Sep-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18369606	08714729864028	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18370835	08714729863960	30-Sep-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18372718	08714729863847	30-Sep-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18372913	08714729864028	30-Sep-18
H74939277545210	Guiding Sheath 5F, 45 cm, MP, CC	18376452	08714729863854	30-Sep-18
H74939277790120	Guiding Sheath 7F, 90 cm, ST, TB	18376474	08714729864035	30-Sep-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18378781	08714729863847	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18378782	08714729863861	30-Sep-18
H74939277645120	Guiding Sheath 6F, 45 cm, ST, TB	18378785	08714729863878	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18378786	08714729863960	30-Sep-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18378787	08714729864066	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18379137	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18380109	08714729863908	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18380645	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18381111	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18382497	08714729863960	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18384084	08714729863908	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18385102	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18385156	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18385953	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18386014	08714729863960	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18386109	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18387930	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18388909	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18389826	08714729863960	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18390056	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18390426	08714729863960	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18390877	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18391074	08714729863908	30-Sep-18
H74939277845110	Guiding Sheath 8F, 45 cm, ST, CC	18391573	08714729864066	30-Sep-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18393119	08714729863984	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18396701	08714729864004	30-Sep-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18396724	08714729863984	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18396749	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18396750	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18396751	08714729863960	30-Sep-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18396783	08714729863984	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18396822	08714729863960	30-Sep-18
H74939277690220	Guiding Sheath 6F, 90 cm, MP, TB	18397703	08714729863953	30-Sep-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18398900	08714729863847	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18399821	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18400036	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18401277	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18401638	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18402267	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18402268	08714729863861	30-Sep-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277690120	Guiding Sheath 6F, 90 cm, ST, TB	18402766	08714729863939	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18404144	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18405356	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18406140	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18406155	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18406487	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18406702	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18408996	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18409056	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18410045	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18410358	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18411125	08714729863861	30-Sep-18
H74939277845120	Guiding Sheath 8F, 45 cm, ST, TB	18411541	08714729864073	30-Sep-18
H74939277865120	Guiding Sheath 8F, 65 cm, ST, TB	18412395	08714729864103	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18413522	08714729863861	30-Sep-18
H74939277890120	Guiding Sheath 8F, 90 cm, ST, TB	18414133	08714729864127	30-Sep-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18414743	08714729863847	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18416005	08714729863861	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18417037	08714729863960	30-Sep-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18417524	08714729863984	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18419071	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18419406	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18420635	08714729864004	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18421331	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18421862	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18421876	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18421913	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18422051	08714729863861	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18422346	08714729863960	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18422702	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18423073	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18423158	08714729863960	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18423159	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18423200	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18423242	08714729863960	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18423243	08714729863960	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18423244	08714729863861	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18425083	08714729863960	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18425256	08714729863861	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18425451	08714729863960	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18427324	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18427931	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18428880	08714729863960	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18428956	08714729863861	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18430109	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18430956	08714729864004	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18431440	08714729864004	30-Sep-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18432504	08714729863960	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18432606	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18433289	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18434514	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18435645	08714729863861	30-Sep-18
H74939277765110	Guiding Sheath 7F, 65 cm, ST, CC	18436440	08714729864004	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18436799	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18437551	08714729863885	30-Sep-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18439565	08714729863908	30-Sep-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18440410	08714729863922	30-Sep-18
H74939277645120	Guiding Sheath 6F, 45 cm, ST, TB	18441193	08714729863878	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18441831	08714729863861	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18442146	08714729863908	30-Sep-18
H74939277745120	Guiding Sheath 7F, 45 cm, ST, TB	18442149	08714729863977	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18442863	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18443685	08714729863861	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18448954	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18448955	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18448956	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18449140	08714729863908	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18449627	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18449713	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18449722	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18449984	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18450853	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18452131	08714729863861	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18454031	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18454223	08714729863908	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18454238	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18454418	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18454657	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18456169	08714729863861	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18456678	08714729863861	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18457133	08714729863908	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18458737	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18459199	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18459514	08714729863885	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18459820	08714729863908	30-Sep-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18460428	08714729863908	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18461523	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18463199	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18464087	08714729863885	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18464740	08714729863885	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18464753	08714729863861	30-Sep-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18464841	08714729863885	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18466537	08714729863861	30-Sep-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18467575	08714729863922	30-Sep-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18468976	08714729863922	30-Sep-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18470080	08714729863861	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18474687	08714729863885	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18475195	08714729863922	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18475394	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18475764	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18475951	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18477455	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18478618	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18480307	08714729863922	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18480349	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18480878	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18480888	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18481889	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18482629	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18483178	08714729863861	31-Oct-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18484150	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18485041	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18485145	08714729863861	31-Oct-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18485446	08714729863984	31-Oct-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18485458	08714729864028	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18485622	08714729863861	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18485980	08714729863885	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18485983	08714729863922	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18487613	08714729863885	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18488132	08714729863885	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18488580	08714729863885	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18489651	08714729863885	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18490347	08714729863885	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18490767	08714729863922	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18492409	08714729863922	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18492899	08714729863885	31-Oct-18
H74939277790220	Guiding Sheath 7F, 90 cm, MP, TB	18494839	08714729864059	31-Oct-18
H74939277865110	Guiding Sheath 8F, 65 cm, ST, CC	18494949	08714729864097	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18495727	08714729863922	31-Oct-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18496126	08714729864028	31-Oct-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18497504	08714729864028	31-Oct-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18497935	08714729864028	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18498993	08714729863847	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18499491	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18499534	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18499822	08714729863922	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18500071	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18500122	08714729863861	31-Oct-18
H74939277790110	Guiding Sheath 7F, 90 cm, ST, CC	18500214	08714729864028	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18500622	08714729863861	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18501090	08714729863847	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18501105	08714729863922	31-Oct-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18501107	08714729863984	31-Oct-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18501431	08714729863984	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18501860	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18501916	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18502768	08714729863922	31-Oct-18
H74939277745210	Guiding Sheath 7F, 45 cm, MP, CC	18502978	08714729863984	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18504335	08714729863847	31-Oct-18
H74939277865110	Guiding Sheath 8F, 65 cm, ST, CC	18505221	08714729864097	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18505426	08714729863908	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18505776	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18505777	08714729863861	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18507287	08714729863908	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18509458	08714729863908	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18509673	08714729863847	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18510382	08714729863908	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18510497	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18511022	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18512359	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18513291	08714729863861	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18514352	08714729863847	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18514888	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18515500	08714729863861	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18515987	08714729863885	31-Oct-18

91107242FA - Chariot™ Guiding Sheath Removal

Material	Product Description	Batch	GTIN	Expiration Date
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18516682	08714729863885	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18525714	08714729863847	31-Oct-18
H74939277645210	Guiding Sheath 6F, 45 cm, MP, CC	18525715	08714729863885	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18528141	08714729863908	31-Oct-18
H74939277665110	Guiding Sheath 6F, 65 cm, ST, CC	18530529	08714729863908	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18531726	08714729863960	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18532034	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18532301	08714729863861	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18533645	08714729863960	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18534161	08714729863861	31-Oct-18
H74939277690210	Guiding Sheath 6F, 90 cm, MP, CC	18534901	08714729863946	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18536914	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18537114	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18537129	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18537337	08714729863861	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18537617	08714729863960	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18539387	08714729863960	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18541194	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18543779	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18543881	08714729863922	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18545803	08714729863960	31-Oct-18
H74939277745110	Guiding Sheath 7F, 45 cm, ST, CC	18546263	08714729863960	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18546975	08714729863861	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18547188	08714729863861	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18548330	08714729863922	31-Oct-18
H74939277690110	Guiding Sheath 6F, 90 cm, ST, CC	18550813	08714729863922	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18551985	08714729863861	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18556150	08714729863847	31-Oct-18
H74939277545110	Guiding Sheath 5F, 45 cm, ST, CC	18556155	08714729863847	31-Oct-18
H74939277890110	Guiding Sheath 8F, 90 cm, ST, CC	18556164	08714729864110	31-Oct-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18585792	08714729863861	30-Nov-18
H74939277645110	Guiding Sheath 6F, 45 cm, ST, CC	18585793	08714729863861	30-Nov-18