

URGENT: FIELD SAFETY NOTICE
LIGAMAX – 5mm Endoscopic Multiple Clip Applier
Product Codes: EL5ML
(Multiple Lot Numbers) – Voluntary Product Recall (Removal)

07/11/2019

Dear Valued Customer,

Ethicon has initiated a voluntary medical device recall (removal) of specific lots of LIGAMAX – 5mm Endoscopic Multiple Clip Applier distributed in Austria, Belgium, Denmark, Estonia, Finland, France, Germany, Great Britain, Ireland, Netherlands, Slovenia, Sweden, Switzerland, Qatar, Saudi Arabia Kuwait and United Arab Emirates.

Ethicon identified that a potential exists that certain LIGAMAX devices within the impacted lots may have small holes in the Tyvek lidding which could result in a breach of sterility. In the event that the sterility barrier has been breached, there is a chance that a pathogen may be introduced and if unrecognized or untreated may cause life threatening infection. As the population at risk is most likely to receive prophylactic antibiotics the probability of harm is extremely rare. We have identified the root cause and have implemented corrective actions to address the issue and prevent reoccurrence. As of the date of this communication, no adverse events or complaints have been reported for this issue.

Our records indicate that you may have ordered or received product subject to this recall. **PLEASE DISTRIBUTE THIS INFORMATION TO ALL STAFF WITHIN YOUR FACILITY WHO USE LIGAMAX – 5mm ENDOSCOPIC MULTIPLE CLIP APPLIER.** The dates of distribution for affected products were from July 31, 2019 thru October 4, 2019.

EFFECTIVE IMMEDIATELY – DO NOT USE OR DISTRIBUTE THE FOLLOWING PRODUCT CODE/ LOT:

Table 1 – Product Subject to this Recall

PRODUCT NAME	PRODUCT CODE	PRODUCT LOT	DESCRIPTION / SIZE
LIGAMAX	EL5ML	T93X4R	5mm Endoscopic Multiple Clip Applier
		T93X9D	
		T93Y96	
		T94118	
		T9416M	
		T94G99	

Health care practitioners who have treated patients using this product should follow those patients post-operatively in the usual manner with no additional action required.

This voluntary medical device recall has been communicated to all impacted Health Authorities in EEA, Switzerland and Turkey.

URGENT: FIELD SAFETY NOTICE
LIGAMAX – 5mm Endoscopic Multiple Clip Applier
Product Codes: EL5ML
(Multiple Lot Numbers) – Voluntary Product Recall (Removal)

IDENTIFICATION OF PRODUCT SUBJECT TO THIS RECALL (Removal):

Product subject to the recall in your inventory can be identified by product code described in above table. All unused LIGAMAX – 5mm Endoscopic Multiple Clip Appliers subject to this recall are required to be returned. Please utilize attachment 1 for assistance in identifying the product lots subject to this recall.

ACTION REQUIRED:

1. Examine your inventory immediately to determine if you have product subject to this recall on hand and quarantine such product(s).
2. Remove the product subject to this recall and communicate the issue to relevant operating room or materials management personnel, or anyone else in your facility who needs to be informed.
3. If any product subject to this recall has been forwarded to another facility, contact that facility to arrange return. Please consider including a copy of this recall letter when communicating.
4. Complete the Business Reply Form (BRF) (Attachment 2) confirming receipt of this notice and return to MDFieldactionsUKIRL@its.jnj.com or fax to 0113 307 1808 within three (3) business days. **Please return the BRF even if you do not have product subject to this recall.**
5. Follow instructions in the letter and immediately return any inventory of LIGAMAX – 5mm Endoscopic Multiple Clip Appliers subject to this recall.
6. Keep this notice visibly posted for awareness until all product subject to this recall has been returned. While processing your returns, please maintain a copy of this notice with the product subject to this recall and keep a copy for your records.
7. Customers are required to return unused **impacted LIGAMAX – 5mm Endoscopic Multiple Clip Appliers** subject to this recall that are in their inventory immediately. To receive credit reimbursement, customers must return product subject to this recall.
8. To return product subject to this recall, please complete the attached BRF and return to MDFieldActionsUKIRL@its.jnj.com or fax to 0113 307 1808, our returns team will then contact you separately to arrange the return.

If you require any assistance with returning product lots subject to this recall, please contact MDFieldActionsUKIRL@its.jnj.com OR your local sales representative.

At Ethicon, our first priority is to our customers and their patients, and that includes the safe and effective use of our products. We recognize the recall of this product may be disruptive to your facility and we apologize for any inconvenience this may cause.

If you have additional questions regarding this recall, please contact MDFieldActionsUKIRL@its.jnj.com.

URGENT: FIELD SAFETY NOTICE
LIGAMAX – 5mm Endoscopic Multiple Clip Applier
Product Codes: EL5ML

(Multiple Lot Numbers) – Voluntary Product Recall (Removal)

As with any medical device, adverse reactions or quality problems experienced with the use of this product should be reported to your Sales Representative, directly to Ethicon, or your National Health Authority. If you have any further questions related to this notice or if you need any additional communications, please contact your local Sales Representative.


Attachments:

Attachment 1: Product Identification Tool
Attachment 2: Business Reply Form (BRF)


Attachment 1: Product Identification Tool for LIGAMAX– 5mm Endoscopic Multiple Clip Applier

Please refer to below in order to identify location of product code, expiration date, and lot number for LIGAMAX– 5mm Endoscopic Multiple Clip Applier subject to this recall by using the packaging labels.

Front of Device Carton:



Sealed Device Tyvek Package:



URGENT: FIELD SAFETY NOTICE
LIGAMAX – 5mm Endoscopic Multiple Clip Applier
Product Codes: EL5ML
(Multiple Lot Numbers) – Voluntary Product Recall (Removal)

Attachment 2: Business Reply Form

Business Reply Form (BRF)

Your timely response to this customer notification is requested. Please complete and fax or email this form to 0113 307 1808 or MDFieldActionsUKIRL@its.jnj.com **within 3 business days, even if you do not have product subject to this recall (removal) to return.**

If you have product subject to this recall to return, please make a photocopy of your completed Business Reply Form and enclose with your return. Thank you for your cooperation.

Print Name of Person Completing Business Reply Form:	Telephone Number:
Account Name & Address (including collection department):	Date:
Email Address:	
Signed*: <small><i>*Your signature provides confirmation that you have received and understood this notification</i></small>	
<i>Your comments are welcome.</i>	

Product Inventory – please check one

- We have **NO** inventory of product subject to this recall (removal).
 We have product subject to this recall (removal) and are returning the following products:

PRODUCT NAME/CODE	LOT NUMBER	EXPIRATION DATE	RETURNING QUANTITY (Eaches)