

---

## **Urgent Field Safety Notice**

Date: 28<sup>th</sup> January 2009

### **Products affected:**

Leckey Mobile Advance Seat, Size 1 and 2  
Leckey Contoured Advance Seat, Size 1 and 2  
Leckey Squiggles Saddle Seat  
Leckey Squiggles Early Sitting System

**FSCA-identifier – MHRA -012 28-01-09**

**Type of action – Replacement of Hip belts**

---

### **Attention:**

All users, carers, parents and health professional staff

### **Details on affected devices:**

Potential defect on plastic buckle used on hip belt on the above products supplied between September 2007 and December 2008

### **Description of the problem:**

It has come to our attention that plastic buckles on the Hip belts supplied with the above products may not be manufactured from the correct material and over time the buckle clip may not lock securely.

To ensure ongoing safe use of the product we will be replacing all Hip belts supplied on the above products between September 2007 and December 2008.

### **Advise on action to be taken by the user:**

To expedite the replacement of hip belts we are requesting customers / end users who have received one of the above products in the designated timeframe to contact our customer services department to arrange for a replacement hip belt to be sent out.

### **UK Service Centre UK 0800 318265 or ROI 1800 626020**

Until you receive your replacement hip belt we are advising customers to be vigilant and check Hip belts daily to ensure the user is secure in their seating system.

### **Transmission of this Field Safety Notice:**

This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please transfer this notice to other organisations on which this action has an impact.

Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

### **Contact reference person:**

Noel McQuaid  
Technical Director  
James Leckey Design Ltd.  
Kilwee Business Park, Dunmurry, Belfast, BT17 0HD.  
+44 (0)28 9060 22 77

The undersigned confirms that this notice has been forwarded to the appropriate Regulatory Agency

Signature

