



Roma Medical Aids Ltd  
York Road  
Bridgend Ind. Estate  
Bridgend  
CF31 3TB  
T: 01656 674488  
F: 01656 674499

FSN 2011.02  
Ref: 2010/010/013/401/005

02/03/2011

## Urgent Field Safety Notice Corrective Action Romachair Model 4730 / 4740

Attention: Ward Managers where Romachairs are in use

Roma Medical has been made aware of an incident where an arm became detached from the chair causing the user to fall. It is possible that the arm became detached due to being incorrectly fitted. The arm detaching during use could lead to injury to the user.

To overcome this problem Roma Medical has updated the user manuals to clarify the fitting procedure and developed armrests which will prevent incorrect fitting.

Please contact Roma Medical Aids to arrange replacement armrests and updated user manuals.

Distribution:

This Field Safety Notice should be distributed to all those involved with Infection Control / Ward Managers.

The undersigned confirms that this notice has been transmitted to the appropriate Regulatory Agency (M.H.R.A).

A handwritten signature in black ink, appearing to read "P. Richards".

Paul Richards  
Technical Manager  
Roma Medical Aids Ltd