

Circulation

 Feeling faint or dizzy Loss of consciousness Feeling confused Pale and clammy skin Rapid or weak pulse

C

Vomiting, abdominal pain, diarrhoea or a rash may also be present.



adrenaline (epinephrine)

Read the patient information leaflet supplied with your medicine before use.

Side effects may be reported via HPRA Pharmacovigilance, website; www.hpra.ie

Anapen[®] instructional video can be viewed at You Tube www.tiny.cc/anapen

01 686 9342 (Freephone number) in Republic of Ireland

anapen.ie

Anapen[®] adrenaline (epinephrine)

anapen.ie

The training device and the active device should not be routinely carried together to minimise the risk of confusion between the two devices in an emergency situation.

• Always carry your two Anapen® devices with you to ensure the recommended dose can be given

 Make sure you and those closest to you are familiar with these instructions, so you can avoid uncertainty in an emergency

01 686 9342 (Freephone number) in Republic of Ireland

ANA/IE/15/22 | July 2022

bioprojet

How to use **Anapen**[®] adrenaline (epinephrine) Remove the black needle Remove the arev Hold the open end shield by pulling hard, in safety cap from the (needle end) of Anapen® the direction of the arrow. red firing button. against the outer part of This also removes a grey your thigh. protective needle shield.

Top tips: Anaper clothir

 Anapen[®] can be given through clothing, including denim

 Use Anapen[®] without delay at the first onset of symptoms If possible lie down with your legs elevated following administration of Anapen[®]



You should always get medical help immediately after using Anapen. Call 999 or 112; ask for an ambulance and state 'anaphylaxis' even if symptoms are improving.

Press the red firing button so that it clicks. Keep holding the Anapen® Auto-Injector against your thigh for 10 seconds. Slowly remove Anapen® from your thigh. Then gently massage the injection area. The injection indicator will have turned red. This shows that the injection is complete. If the injection indicator is not red, you must repeat the injection with a new Anapen®.

5

After the injection, the needle sticks out. To cover it, click the wide end of the black needle shield back on the open end (needle end) of Anapen[®] Auto-Injector.

