

THE AURUM RANGE

Ready to use in a single twist

Connector compatibility.

The new Aurum Range in PRTC (plastic rigid tip cap) format is Luer Lock compatible to help aid drug administration.

The Aurum Range of pre-filled syringes are now compatible with:

- Baxter Healthcare: Clearlink and V-Link needlefree Luer activated devices
- B.Braun Medical - Ultrasite/Safeflow Connectors
- Becton, Dickinson U.K. Ltd - BD Q-Syte™ Connectors
- BMS Critical Care - EasiFlush Needlefree Connector
- CareFusion™ Maxplus Displacement Connectors
- CareFusion™ SmartSite® Valve
- CODAN Ltd - SWAN-LOCK® Connectors
- Fannin - MicroCLAVE and MicroCLAVE Clear: with the use of the CS-25 Clear Adaptor*

*CS-25 Clear Adaptor and emergency poster are available free of charge to your Pharmacy or Supplies Department. Contact Fannin, Wellingborough 0800 212 827, Belfast 02890 735 588 or email enquiries@fanninuk.com.

- Fannin - Flowart Needlefree
- ICU Medical CLC 2000 Connector
- ICU Medical Tego D1000 Connector
- Sendal UK Ltd - TPS Connectors
- Vygon Bionector and Bionector 'S': with ATP431 adaptor*

*ATP431 adaptor and information poster available free of charge from Vygon (UK) Ltd. Contact: 01793 748800 or e-mail: vygon@vygon.co.uk

Other connectors are not currently recommended. Information correct as of January 2013.

Prescribing information available on request.

Date of preparation: January 2013.

Job no. CRIT/12/2012/099.

For further information, please contact our Customer Services team on **0800 028 7933** or Medical Information on **01277 266 600**, alternatively, email us at **medinfo@martindalepharma.co.uk**.

MARTINDALE PHARMA®
Making lives better