



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

Baxter

FA-2022-016

Date

Customer Name

FSCA-identifier: FA-2022-016

Type of action: Field Safety Notice

Subject: AD271B Siderails disengagement risk

Commercial name of the impacted product:
Hill-Rom 900, Ampera and Centuris Pro

Impacted devices:

See Appendix A for impacted Serial Numbers

Name	Model
Hill-Rom 900	LI900B2
Ampera	LI901A0
Centuris Pro	LI900B1

All products manufactured between 20 April 2020 to 10 March 2022 may be impacted.

To: Facility Risk Manager/Facility Administrator, Chief Executive Facility Administrator, Facility Engineer, Vigilance Manager, Biomedical Engineering, Medical Device Liaison Officer

Description of the problem:

Hillrom, a Baxter company, has become aware of a potential issue where the siderail of the LI900B1, LI900B2 and LI901A0 Bed Assembly may not properly engage. If the siderail is not properly engaged, it could result in the siderail falling and would therefore no longer protect the patient, possibly leading to patient fall.

Background:

Hillrom produces bed assemblies (product models LI900B1, LI900B2 and LI901A0) that are equipped with AD271B ¾ siderails (patient safety barriers). In February 2022, Hillrom detected a siderail that was not properly engaged into the bed frame due to a burr in the plastic connector. See Figure 1 below.



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

Baxter

FA-2022-016



Figure 1: Siderail connector with burr


Potential Risk:

If the siderail is not properly engaged, it could result in the siderail falling and would therefore no longer protect the patient, possibly leading to patient fall.

There have been no complaints of this issue occurring in the field.

Actions to be taken by Customer:

- Please identify if you have the impacted product, per Appendix A, in your facility.
Serial number and manufacturing date can be found on frame identification plate – see example below.



- Please follow IFU instructions on how to disengage and re-engage your siderail. Please ensure there is an audible “click” to confirm it is properly engaged. See images below



Hillrom™


Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice


Baxter

FA-2022-016

Removing the long siderails



Installing the long siderails



- Return the enclosed response form within two weeks, confirming inspection complete.

Action to be taken by the Distributor:

Please share this Field Safety Notice with your end users and complete the attached response form and return to hillromPLU016@Sedgwick.com within two weeks.

Please support your end users if inspection by technician is requested.

Action to be taken by Hillrom:

Hillrom have addressed this issue to prevent burrs from occurring in future productions.

Contact Reference Person:

Should you have any questions regarding this Field Safety Notice, please contact Hillrom Technical Support using email or number below.

Market /Region/ Country	Phone Number	Technical Support Email
Austria	(+ 43) 2 243 285 50	service-administration.dach@hillrom.com
Germany	(+ 49) 2 0149869500	service.dach@hillrom.com
Switzerland	(+ 41) 8 48 811530	service-administration.dach@hillrom.com
Netherlands	(+ 31) 347 323 532	service.nl@hillrom.com
Spain	(+ 34) 9 36856000	asistencia@hillrom.com
Italy	(+ 39) 02 950541	assistenza.tecnica@hillrom.com
France	(+ 33) 0 820 012345	sav@hillrom.com

 Hillrom™ Hill-Rom S.A.S, ZI du Talhouet 56330 Pluvigner, France	Field Safety Notice	 FA-2022-016
--	----------------------------	---

Sweden	(+ 46) 20-781030	ordernordic@hillrom.com
UK/Ireland	(+ 44) 1530 562176	UKCustomerCare@hillrom.com
Eastern Europe Countries	Contact your Local Hillrom Distributor	

Transmission of this Field Safety Notice:

Please ensure this notice is circulated to all appropriate personnel. This may include, but is not limited to:

A&E departments	In-house maintenance staff
Adult intensive care units	IV nurse specialists
All wards & Clinics	Medical directors
Biomedical engineering staff	Nursing executive directors
Clinical governance leads	Oncology units
Day case theatres	Pediatric intensive care units
EBME departments	Risk managers
Equipment stores & Libraries	Supplies managers
Health and safety managers	Theatres

The undersigned confirms that this notice has been communicated with the appropriate Regulatory Agency.

Sincerely



Maurice Herve
Director, Quality Assurance



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

Baxter

FA-2022-016

Response Form / Receipt

Subject: AD271B Siderail disengagement risk (FA-2022-016)

It is important that you return this form/receipt as acknowledgement of your receipt and provide us with the necessary information.

Please complete the following with the correct information and **return this Response Form** within two weeks.

Thank you!

Hillrom account number (if known): _____

Name of the facility: _____

Address of the facility: _____

City: _____ Zip: _____ Country: _____

Facility Contact Person Name: (print)

Signature: _____ Date: ____/____/____

Title: _____ Phone: _____

Email: _____

Check actions taken:

We have reviewed and understand the attached Field Safety Notice

Results from the inspection of our product inventory show:

We do not have any potentially affected products per Appendix A.

We confirm we have inspected and re-engaged the siderail with an audible "click" as per IFU instructions. Please list serial numbers below.

We confirm we have inspected and re-engaged the siderail however, we require a Hillrom technician to confirm it is conforming per IFU instructions. Please list serial numbers below.

Please identify your impacted serial numbers in the table below.

Serial Number	Serial Number	Serial Number

Response form shall be returned to hillromPLU016@Sedgwick.com within two weeks



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

Appendix A – Impacted Serial Numbers

Serial Number	Serial Number	Serial Number	Serial Number	Serial Number	Serial Number	Serial Number
HRP004138921	HRP004252551	HRP004177852	HRP004173428	HRP004182265	HRP004181045	HRP004128753
HRP004144524	HRP004252552	HRP004177853	HRP004173429	HRP004182266	HRP004181046	HRP004128754
HRP004144525	HRP004252553	HRP004177854	HRP004173430	HRP004182267	HRP004181049	HRP004153452
HRP004210429	HRP004252559	HRP004177855	HRP004208832	HRP004153313	HRP004181050	HRP004153453
HRP004210430	HRP004252560	HRP004177856	HRP004208833	HRP004153314	HRP004181051	HRP004153454
HRP004210445	HRP004252723	HRP004177857	HRP004208834	HRP004183572	HRP004181052	HRP004153455
HRP004210446	HRP004252967	HRP004177858	HRP004208835	HRP004277612	HRP004181053	HRP004153456
HRP004164745	HRP004252968	HRP004177859	HRP004208836	HRP004277613	HRP004160677	HRP004153481
HRP004186146	HRP004252969	HRP004177860	HRP004208837	HRP004277614	HRP004160678	HRP004153482
HRP004255189	HRP004252970	HRP004177861	HRP004208842	HRP004277615	HRP004160679	HRP004153483
HRP004255190	HRP004254473	HRP004177611	HRP004208843	HRP004277616	HRP004160680	HRP004153592
HRP004255194	HRP004254474	HRP004177453	HRP004208844	HRP004277617	HRP004160681	HRP004153593
HRP004255215	HRP004254475	HRP004177847	HRP004208845	HRP004277618	HRP004160682	HRP004176861
HRP004255216	HRP004254476	HRP004177848	HRP004208846	HRP004277619	HRP004160683	HRP004176862
HRP004255356	HRP004254477	HRP004177849	HRP004209162	HRP004277620	HRP004160684	HRP004176863
HRP004255357	HRP004199037	HRP004177851	HRP004209163	HRP004277621	HRP004176499	HRP004176864
HRP004255358	HRP004199038	HRP004258002	HRP004209164	HRP004111253	HRP004176500	HRP004176865
HRP004255359	HRP004199039	HRP004258003	HRP004209165	HRP004111254	HRP004176501	HRP004176866
HRP004255605	HRP004199040	HRP004258223	HRP004177421	HRP004113599	HRP004176503	HRP004176867
HRP004255606	HRP004199041	HRP004258224	HRP004177425	HRP004113608	HRP004176889	HRP004176868
HRP004255607	HRP004111237	HRP004113940	HRP004177426	HRP004113609	HRP004176890	HRP004176869
HRP004255610	HRP004161304	HRP004113941	HRP004177427	HRP004113610	HRP004176891	HRP004176870
HRP004255622	HRP004161305	HRP004113942	HRP004176754	HRP004113611	HRP004176892	HRP004176871
HRP004255791	HRP004161306	HRP004113943	HRP004176755	HRP004113612	HRP004176893	HRP004176872
HRP004255792	HRP004161307	HRP004113948	HRP004176756	HRP004113669	HRP004176894	HRP004176873
HRP004255793	HRP004161308	HRP004113949	HRP004176757	HRP004263879	HRP004176895	HRP004176874
HRP004255794	HRP004161314	HRP004113950	HRP004176758	HRP004263880	HRP004176902	HRP004176875
HRP004255795	HRP004161315	HRP004113951	HRP004176759	HRP004263881	HRP004176903	HRP004176876
HRP004256057	HRP004161331	HRP004113952	HRP004176760	HRP004278290	HRP004176904	HRP004176881
HRP004256058	HRP004161332	HRP004113958	HRP004176761	HRP004278291	HRP004176905	HRP004176882
HRP004256059	HRP004161333	HRP004113959	HRP004176762	HRP004278292	HRP004176906	HRP004176883
HRP004256060	HRP004161340	HRP004113960	HRP004176763	HRP004278293	HRP004196014	HRP004176884
HRP004256061	HRP004161341	HRP004113961	HRP004176764	HRP004145990	HRP004196015	HRP004176885
HRP004256206	HRP004161342	HRP004113962	HRP004176765	HRP004145991	HRP004196020	HRP004176886
HRP004256207	HRP004161343	HRP004113963	HRP004176766	HRP004145992	HRP004196021	HRP004128746
HRP004256208	HRP004161344	HRP004113964	HRP004176767	HRP004145993	HRP004196022	HRP004128747
HRP004256209	HRP004161345	HRP004175624	HRP004176768	HRP004145994	HRP004113600	HRP004128748
HRP004256217	HRP004161385	HRP004175630	HRP004176769	HRP004146546	HRP004113607	HRP004128749
HRP004256218	HRP004161386	HRP004175631	HRP004211691	HRP004146547	HRP004115232	HRP004128750
HRP004261800	HRP004161392	HRP004160774	HRP004211692	HRP004146548	HRP004115235	HRP004128751
HRP004161992	HRP004161494	HRP004160775	HRP004211693	HRP004146549	HRP004115236	HRP004128752
HRP004161993	HRP004161495	HRP004239905	HRP004211801	HRP004146550	HRP004115237	HRP004128753
HRP004161994	HRP004161496	HRP004239907	HRP004211802	HRP004037054	HRP004115238	HRP004128754
HRP004161995	HRP004161497	HRP004278559	HRP004211803	HRP004190989	HRP004115243	HRP004129102
HRP004161996	HRP004161498	HRP004278580	HRP004211804	HRP004190990	HRP004115244	HRP004176505
HRP004161997	HRP004161499	HRP004278581	HRP004211805	HRP004190991	HRP004115245	HRP004176508
HRP004161998	HRP004161500	HRP004278582	HRP004212061	HRP004190992	HRP004115246	HRP004176509
HRP004161999	HRP004161501	HRP004278583	HRP004212062	HRP004190993	HRP004115247	HRP004176510
HRP004162000	HRP004161502	HRP004278584	HRP004212063	HRP004190994	HRP004115248	HRP004176511
HRP004162001	HRP004161503	HRP004278585	HRP004212064	HRP004190995	HRP004115249	HRP004176512
HRP004114560	HRP004161504	HRP004278586	HRP004212065	HRP004187424	HRP004115250	HRP004176513
HRP004114561	HRP004161229	HRP004278587	HRP004212066	HRP004187425	HRP004110573	HRP004176517



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

FA-2022-016

HRP004114562	HRP004161231	HRP004278588	HRP004212071	HRP004187426	HRP004110588	HRP004176528
HRP004114563	HRP004161232	HRP004278589	HRP004212072	HRP004187431	HRP004110590	HRP004176530
HRP004114564	HRP004161233	HRP004278590	HRP004212110	HRP004264271	HRP004110571	HRP004176531
HRP004114565	HRP004161334	HRP004278591	HRP004212111	HRP004118267	HRP004110572	HRP004176532
HRP004114566	HRP004161505	HRP004278592	HRP004212112	HRP004118268	HRP004207958	HRP004176533
HRP004114567	HRP004161506	HRP004278593	HRP004212113	HRP004118269	HRP004207959	HRP004176534
HRP004114568	HRP004161521	HRP004278594	HRP004209166	HRP004118270	HRP004207960	HRP004176535
HRP004114581	HRP004161523	HRP004278595	HRP004209167	HRP004118271	HRP004207961	HRP004176536
HRP004114582	HRP004161524	HRP004278596	HRP004209208	HRP004118272	HRP004207962	HRP004176537
HRP004114583	HRP004161534	HRP004278604	HRP004211298	HRP004118273	HRP004147893	HRP004176538
HRP004114584	HRP004161535	HRP004278605	HRP004211299	HRP004177980	HRP004147894	HRP004176539
HRP004114585	HRP004161536	HRP004278606	HRP004209168	HRP004177981	HRP004147895	HRP004176540
HRP004114586	HRP004161537	HRP004278607	HRP004209445	HRP004177982	HRP004147896	HRP004176541
HRP004114587	HRP004161555	HRP004278608	HRP004209759	HRP004177983	HRP004147947	HRP004176542
HRP004114588	HRP004163988	HRP004278609	HRP004209760	HRP004177984	HRP004148755	HRP004176543
HRP004114589	HRP004163989	HRP004278610	HRP004209761	HRP004177985	HRP004148975	HRP004176544
HRP004114595	HRP004163990	HRP004207295	HRP004209762	HRP004177986	HRP004148976	HRP004176545
HRP004114596	HRP004163991	HRP004207296	HRP004209763	HRP004177987	HRP004148978	HRP004176546
HRP004114558	HRP004163992	HRP004207297	HRP004209764	HRP004177988	HRP004147948	HRP004176857
HRP004114559	HRP004163993	HRP004207298	HRP004209765	HRP004177989	HRP004148756	HRP004176858
HRP004161039	HRP004163994	HRP004207299	HRP004209766	HRP004177990	HRP004148757	HRP004176859
HRP004161040	HRP004155245	HRP004207300	HRP004209996	HRP004177991	HRP004148758	HRP004176860
HRP004161041	HRP004155246	HRP004207301	HRP004209997	HRP004177992	HRP004148977	HRP004181375
HRP004161042	HRP004155264	HRP004207302	HRP004210002	HRP004177993	HRP004148979	HRP004181376
HRP004161050	HRP004155265	HRP004207303	HRP004210003	HRP004177994	HRP004148980	HRP004181377
HRP004161051	HRP004155266	HRP004207304	HRP004210004	HRP004223609	HRP004148981	HRP004181378
HRP004161052	HRP004155273	HRP004242379	HRP004210005	HRP004228107	HRP004148982	HRP004181379
HRP004161053	HRP004155274	HRP004242380	HRP004210006	HRP004256597	HRP004273257	HRP004181380
HRP004161057	HRP004155275	HRP004241810	HRP004210008	HRP004268635	HRP004273258	HRP004181381
HRP004161058	HRP004155276	HRP004241811	HRP004210010	HRP004268636	HRP004273259	HRP004181382
HRP004161063	HRP004155277	HRP004241812	HRP004210011	HRP004268637	HRP004273266	HRP004181383
HRP004161096	HRP004262787	HRP004241813	HRP004210016	HRP004268638	HRP004273277	HRP004181384
HRP004161097	HRP004262788	HRP004241814	HRP004210017	HRP004268639	HRP004273278	HRP004181385
HRP004161098	HRP004262789	HRP004241815	HRP004210018	HRP004268640	HRP004273279	HRP004181386
HRP004161099	HRP004262995	HRP004241816	HRP004210019	HRP004268641	HRP004273280	HRP004181387
HRP004161100	HRP004262996	HRP004241817	HRP004210020	HRP004268642	HRP004270082	HRP004181388
HRP004161101	HRP004263314	HRP004241807	HRP004136048	HRP004268643	HRP004270083	HRP004181373
HRP004161102	HRP004263315	HRP004241808	HRP004136049	HRP004268644	HRP004270084	HRP004181374
HRP004161103	HRP004263316	HRP004241809	HRP004007979.	HRP004268645	HRP004270085	HRP004181667
HRP004161104	HRP004263527	HRP004243125	HRP004269060	HRP004268646	HRP004270086	HRP004181668
HRP004161105	HRP004263528	HRP004243126	HRP004160744	HRP004268647	HRP004270087	HRP004181669
HRP004161106	HRP004263882	HRP004243127	HRP004160745	HRP004268648	HRP004270088	HRP004181670
HRP004161107	HRP004263883	HRP004243128	HRP004160746	HRP004268649	HRP004270089	HRP004181671
HRP004161108	HRP004263884	HRP004243129	HRP004160747	HRP004268650	HRP004270090	HRP004181684
HRP004161119	HRP004264040	HRP004243130	HRP004160748	HRP004268651	HRP004270091	HRP004181685
HRP004161120	HRP004264041	HRP004243271	HRP004160749	HRP004268652	HRP004270092	HRP004181686
HRP004190376	HRP004264357	HRP004243272	HRP004160750	HRP004268653	HRP004270093	HRP004181687
HRP004190386	HRP004264358	HRP004243273	HRP004160751	HRP004268654	HRP004270094	HRP004181688
HRP004190387	HRP004264359	HRP004243274	HRP004160752	HRP004158706	HRP004160780	HRP004181689
HRP004190389	HRP004264360	HRP004243275	HRP004160753	HRP004158707	HRP004160781	HRP004181690
HRP004190565	HRP004160924	HRP004243276	HRP004151967	HRP004278227	HRP004160782	HRP004181691
HRP004190576	HRP004165216	HRP004243277	HRP004151968	HRP004284251	HRP004160783	HRP004181692
HRP004190752	HRP004165211	HRP004243278	HRP004151969	HRP004284255	HRP004160784	HRP004181693
HRP004190753	HRP004165212	HRP004269058	HRP004151970	HRP004232232	HRP004160785	HRP004181694



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004190375	HRP004165213	HRP004269059	HRP004151971	HRP004234551	HRP004160786	HRP004181695
HRP004190377	HRP004165214	HRP004169686	HRP004152172	HRP004234552	HRP004160787	HRP004127906
HRP004190381	HRP004165215	HRP004169687	HRP004152173	HRP004234553	HRP004160788	HRP004127907
HRP004190391	HRP004165217	HRP004169688	HRP004152174	HRP004234554	HRP004160789	HRP004127908
HRP004190577	HRP004170978	HRP004169689	HRP004152175	HRP004232198	HRP004160790	HRP004127909
HRP004190578	HRP004170979	HRP004169690	HRP004152176	HRP004232199	HRP004160801	HRP004127910
HRP004190579	HRP004170980	HRP004169831	HRP004149190	HRP004232205	HRP004160802	HRP004127911
HRP004190580	HRP004170981	HRP004169832	HRP004149191	HRP004232206	HRP004160803	HRP004127912
HRP004190581	HRP004170982	HRP004169833	HRP004149192	HRP004232233	HRP004160804	HRP004127913
HRP004190582	HRP004170983	HRP004169834	HRP004149193	HRP004273567	HRP004160814	HRP004127914
HRP004190587	HRP004170984	HRP004169835	HRP004149194	HRP004273568	HRP004160815	HRP004127915
HRP004190588	HRP004170985	HRP004170437	HRP004152493	HRP004281250	HRP004160816	HRP004127916
HRP004190744	HRP004171454	HRP004170438	HRP004153132	HRP004281257	HRP004160817	HRP004127921
HRP004190745	HRP004182458	HRP004170439	HRP004153133	HRP004281258	HRP004160838	HRP004127922
HRP004190748	HRP004158701	HRP004170440	HRP004153134	HRP004281259	HRP004160839	HRP004127923
HRP004190749	HRP004158702	HRP004281991	HRP004153135	HRP004281260	HRP004160840	HRP004127924
HRP004190750	HRP004237944	HRP004281992	HRP004153306	HRP004281261	HRP004160841	HRP004127925
HRP004190751	HRP004237945	HRP004161144	HRP004153307	HRP004281262	HRP004160862	HRP004127926
HRP004190754	HRP004237946	HRP004161145	HRP004153308	HRP004281263	HRP004160863	HRP004127927
HRP004190755	HRP004237947	HRP004161150	HRP004153309	HRP004282499	HRP004160865	HRP004127928
HRP004273558	HRP004237948	HRP004161151	HRP004153310	HRP004282501	HRP004160878	HRP004127929
HRP004190374	HRP004141034	HRP004161152	HRP004153311	HRP004282502	HRP004160879	HRP004127935
HRP004265065	HRP004266003	HRP004161153	HRP004153312	HRP004282523	HRP004160913	HRP004127937
HRP004265066	HRP004266004	HRP004161154	HRP004157681	HRP004282524	HRP004160914	HRP004127943
HRP004265067	HRP004266005	HRP004161155	HRP004157682	HRP004282525	HRP004160915	HRP004127944
HRP004265068	HRP004266006	HRP004161156	HRP004157683	HRP004282526	HRP004160916	HRP004127945
HRP004265069	HRP004266007	HRP004161157	HRP004157684	HRP004282527	HRP004160918	HRP004127946
HRP004141496	HRP004139644	HRP004161158	HRP004157685	HRP004282528	HRP004160919	HRP004127947
HRP004141497	HRP004139645	HRP004161159	HRP004157686	HRP004282529	HRP004160920	HRP004127948
HRP004141498	HRP004139646	HRP004161160	HRP004157687	HRP004282530	HRP004160813	HRP004127953
HRP004141499	HRP004139647	HRP004161161	HRP004157688	HRP004282531	HRP004160921	HRP004128096
HRP004173169	HRP004139648	HRP004161162	HRP004157689	HRP004282533	HRP004160922	HRP004128097
HRP004209153	HRP004145368	HRP004161163	HRP004287269	HRP004282534	HRP004160923	HRP004128098
HRP004209154	HRP004145370	HRP004161164	HRP004287270	HRP004282535	HRP004114336	HRP004128099
HRP004234544	HRP004145371	HRP004161165	HRP004287271	HRP004282536	HRP004114355	HRP004128100
HRP004234547	HRP004145372	HRP004161166	HRP004287272	HRP004232860	HRP004235104	HRP004128101
HRP004234548	HRP004145373	HRP004161167	HRP004287273	HRP004274133	HRP004235105	HRP004128102
HRP004236146	HRP004197332	HRP004161168	HRP004287275	HRP004164867	HRP004235106	HRP004128103
HRP004236147	HRP004197333	HRP004161169	HRP004287276	HRP004164868	HRP004235132	HRP004128104
HRP004236148	HRP004197334	HRP004161174	HRP004200344	HRP004164869	HRP004235133	HRP004128105
HRP004236149	HRP004197335	HRP004161175	HRP004238900	HRP004164870	HRP004235134	HRP004128106
HRP004239909	HRP004282387	HRP004161176	HRP004155452	HRP004164871	HRP004235135	HRP004128107
HRP004239910	HRP004282388	HRP004161177	HRP004155453	HRP004187432	HRP004235136	HRP004128112
HRP004243061	HRP004282389	HRP004270618	HRP004155456	HRP004227263	HRP004235137	HRP004128113
HRP004243068	HRP004282390	HRP004270619	HRP004155457	HRP004227264	HRP004235138	HRP004128114
HRP004243069	HRP004282392	HRP004270620	HRP004155624	HRP004227265	HRP004273125	HRP004128115
HRP004243090	HRP004282397	HRP004270621	HRP004155630	HRP004228105	HRP004273126	HRP004128116
HRP004243091	HRP004282398	HRP004270622	HRP004155631	HRP004228106	HRP004273127	HRP004128117
HRP004243131	HRP004282399	HRP004270623	HRP004155637	HRP004278093	HRP004273129	HRP004128118
HRP004243132	HRP004282400	HRP004277293	HRP004155638	HRP004282356	HRP004273130	HRP004128119
HRP004243148	HRP004282401	HRP004277294	HRP004155639	HRP004167198	HRP004273131	HRP004128120
HRP004243149	HRP004126014	HRP004277295	HRP004155284	HRP004167199	HRP004273132	HRP004128232
HRP004243150	HRP004126015	HRP004277296	HRP004155285	HRP004167200	HRP004273133	HRP004128233
HRP004243151	HRP004126016	HRP004277347	HRP004155286	HRP004167201	HRP004274108	HRP004128284



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

FA-2022-016

HRP004243152	HRP004126017	HRP004277348	HRP004155287	HRP004167202	HRP004274122	HRP004128285
HRP004243154	HRP004126018	HRP004235917	HRP004155288	HRP004167181	HRP004140146	HRP004128317
HRP004241727	HRP004168928	HRP004235918	HRP004155458	HRP004240210	HRP004140147	HRP004128318
HRP004180852	HRP004168933	HRP004235919	HRP004155622	HRP004240211	HRP004140148	HRP004128319
HRP004141035	HRP004270281	HRP004235920	HRP004155623	HRP004240212	HRP004140149	HRP004128320
HRP004153907	HRP004270307	HRP004235921	HRP004155278	HRP004269065	HRP004140150	HRP004128321
HRP004180853	HRP004270628	HRP004235644	HRP004155279	HRP004269067	HRP004140151	HRP004128708
HRP004214065	HRP004270631	HRP004235645	HRP004241418	HRP004269206	HRP004140406	HRP004128709
HRP004214116	HRP004270636	HRP004235646	HRP004241419	HRP004269207	HRP004140407	HRP004128710
HRP004214520	HRP004270637	HRP004235647	HRP004241420	HRP004269208	HRP004140408	HRP004128711
HRP004214528	HRP004270638	HRP004235660	HRP004241421	HRP004269210	HRP004140417	HRP004128712
HRP004214531	HRP004270639	HRP004167974	HRP004241422	HRP004269211	HRP004140418	HRP004128713
HRP004214532	HRP004270640	HRP004167975	HRP004241423	HRP004269212	HRP004140419	HRP004128718
HRP004214994	HRP004270641	HRP004167976	HRP004239921	HRP004269213	HRP004178357	HRP004128719
HRP004214995	HRP004270642	HRP004191095	HRP004241804	HRP004269214	HRP004178358	HRP004128720
HRP004214996	HRP004270643	HRP004191096	HRP004241805	HRP004281668	HRP004178359	HRP004128721
HRP004214997	HRP004271741	HRP004191097	HRP004242639	HRP004281670	HRP004178360	HRP004128722
HRP004272597	HRP004271772	HRP004191098	HRP004242659	HRP004281671	HRP004178366	HRP004128723
HRP004272598	HRP004271773	HRP004191099	HRP004155448	HRP004281672	HRP004178367	HRP004128724
HRP004272599	HRP004271774	HRP004191105	HRP004155449	HRP004281846	HRP004181289	HRP004128725
HRP004272600	HRP004271775	HRP004191091	HRP004155450	HRP004131430	HRP004181290	HRP004128726
HRP004272731	HRP004271790	HRP004191092	HRP004155451	HRP004131431	HRP004181291	HRP004128727
HRP004272732	HRP004271791	HRP004191093	HRP004155454	HRP004131432	HRP004181292	HRP004128728
HRP004272733	HRP004271792	HRP004191094	HRP004155455	HRP004131433	HRP004181293	HRP004128729
HRP004272734	HRP004271793	HRP004242381	HRP004240213	HRP004131434	HRP004181294	HRP004128730
HRP004272735	HRP004271794	HRP004242382	HRP004240214	HRP004131435	HRP004167646	HRP004128731
HRP004272736	HRP004271946	HRP004242383	HRP004240215	HRP004131436	HRP004167647	HRP004128732
HRP004272737	HRP004271947	HRP004242384	HRP004240216	HRP004131437	HRP004167648	HRP004128733
HRP004272738	HRP004271948	HRP004242385	HRP004240217	HRP004131438	HRP004167649	HRP004128734
HRP004272844	HRP004271949	HRP004242660	HRP004240218	HRP004131439	HRP004167650	HRP004128735
HRP004272845	HRP004271950	HRP004242661	HRP004155280	HRP004231474	HRP004188905	HRP004128736
HRP004272846	HRP004271951	HRP004242662	HRP004155281	HRP004231475	HRP004188907	HRP004128737
HRP004272847	HRP004271952	HRP004214233	HRP004243701	HRP004231476	HRP004188908	HRP004128741
HRP004272848	HRP004271953	HRP004143989	HRP004243702	HRP004231477	HRP004188924	HRP004128742
HRP004272849	HRP004272346	HRP004143990	HRP004243703	HRP004231478	HRP004188925	HRP004128743
HRP004278918	HRP004272347	HRP004143991	HRP004214064	HRP004231479	HRP004253096	HRP004128744
HRP004280257	HRP004272348	HRP004143992	HRP004214519	HRP004231480	HRP004253097	HRP004128745
HRP004239100	HRP004272349	HRP004143993	HRP004241414	HRP004231481	HRP004253098	HRP004247422
HRP004239101	HRP004272350	HRP004143994	HRP004278225	HRP004231482	HRP004253099	HRP004247423
HRP004239351	HRP004272351	HRP004143995	HRP004278226	HRP004231483	HRP004253100	HRP004247424
HRP004239352	HRP004272352	HRP004143996	HRP004279657	HRP004231484	HRP004007980.	HRP004247425
HRP004239353	HRP004272449	HRP004143997	HRP004279658	HRP004231485	HRP004118166	HRP004247426
HRP004239495	HRP004272450	HRP004143998	HRP004279659	HRP004232183	HRP004118167	HRP004247427
HRP004239496	HRP004272451	HRP004171513	HRP004279660	HRP004232184	HRP004118183	HRP004247428
HRP004239517	HRP004272452	HRP004171518	HRP004279661	HRP004232185	HRP004118185	HRP004249397
HRP004239518	HRP004272453	HRP004171519	HRP004279663	HRP004232186	HRP004118186	HRP004249398
HRP004239519	HRP004272454	HRP004256050	HRP004279664	HRP004232187	HRP004118187	HRP004249399
HRP004239520	HRP004272455	HRP004256051	HRP004280792	HRP004187433	HRP004118203	HRP004267045
HRP004239717	HRP004270629	HRP004256052	HRP004281245	HRP004188579	HRP004118204	HRP004267046
HRP004239718	HRP004271438	HRP004256053	HRP004281246	HRP004188580	HRP004118205	HRP004267047
HRP004239719	HRP004271439	HRP004256054	HRP004281247	HRP004188581	HRP004118206	HRP004267048
HRP004239720	HRP004269061	HRP004173170	HRP004281248	HRP004188592	HRP004118207	HRP004267049
HRP004239721	HRP004269062	HRP004177665	HRP004278359	HRP004188593	HRP004118208	HRP004267050
HRP004239911	HRP004269063	HRP004177666	HRP004278360	HRP004189734	HRP004118210	HRP004267051



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004239913	HRP004269064	HRP004177667	HRP004278361	HRP004189735	HRP004118211	HRP004267053
HRP004239914	HRP004269066	HRP004177668	HRP004278362	HRP004189736	HRP004118212	HRP004267054
HRP004241806	HRP004269633	HRP004177678	HRP004278363	HRP004189737	HRP004118213	HRP004267056
HRP004257773	HRP004269634	HRP004177679	HRP004278364	HRP004190343	HRP004118214	HRP004264907
HRP004258437	HRP004269635	HRP004176101	HRP004278365	HRP004190413	HRP004118215	HRP004264915
HRP004258438	HRP004269636	HRP004176102	HRP004278366	HRP004190415	HRP004118216	HRP004264916
HRP004258664	HRP004269637	HRP004176103	HRP004278367	HRP004190416	HRP004118217	HRP004264917
HRP004258665	HRP004269638	HRP004176104	HRP004278368	HRP004190417	HRP004118218	HRP004264918
HRP004258840	HRP004269639	HRP004176105	HRP004278369	HRP004190003	HRP004118219	HRP004264919
HRP004258841	HRP004269640	HRP004234886	HRP004279662	HRP004190009	HRP004118184	HRP004264920
HRP004177527	HRP004269829	HRP004234887	HRP004280759	HRP004190020	HRP004118224	HRP004264921
HRP004177528	HRP004269830	HRP004234890	HRP004280760	HRP004190021	HRP004118225	HRP004046377
HRP004177598	HRP004269831	HRP004234891	HRP004280761	HRP004190022	HRP004118226	HRP004183337
HRP004177599	HRP004269832	HRP004234892	HRP004280762	HRP004186967	HRP004118227	HRP004183344
HRP004177600	HRP004269834	HRP004273742	HRP004280763	HRP004186968	HRP004118228	HRP004183411
HRP004254901	HRP004269839	HRP004273743	HRP004280764	HRP004186969	HRP004118229	HRP004183656
HRP004254902	HRP004269840	HRP004273744	HRP004280765	HRP004186970	HRP004118231	HRP004183659
HRP004254903	HRP004269841	HRP004273771	HRP004280766	HRP004187085	HRP004118232	HRP004183660
HRP004254904	HRP004270115	HRP004273773	HRP004280767	HRP004187086	HRP004118233	HRP004183662
HRP004254905	HRP004270116	HRP004273769	HRP004280791	HRP004187087	HRP004118234	HRP004183666
HRP004254906	HRP004270117	HRP004273770	HRP004280793	HRP004187088	HRP004118235	HRP004183667
HRP004252243	HRP004270118	HRP004273772	HRP004280797	HRP004190339	HRP004118236	HRP004183674
HRP004252244	HRP004270119	HRP004274026	HRP004280798	HRP004190340	HRP004118237	HRP004183691
HRP004252246	HRP004270120	HRP004274574	HRP004278224	HRP004190341	HRP004118238	HRP004184045
HRP004252299	HRP004270121	HRP004160685	HRP004278234	HRP004190342	HRP004118239	HRP004237196
HRP004252300	HRP004270122	HRP004160686	HRP004278235	HRP004190344	HRP004118240	HRP004237201
HRP004252301	HRP004270280	HRP004160687	HRP004278236	HRP004190412	HRP004118241	HRP004237203
HRP004252302	HRP004270302	HRP004160689	HRP004278237	HRP004190414	HRP004118242	HRP004222312
HRP004254928	HRP004270303	HRP004160688	HRP004278238	HRP004163322	HRP004118243	HRP004222313
HRP004254979	HRP004270304	HRP004160690	HRP004278239	HRP004163323	HRP004118263	HRP004222314
HRP004254981	HRP004270305	HRP004160691	HRP004278240	HRP004163324	HRP004118264	HRP004222315
HRP004254982	HRP004270306	HRP004160696	HRP004278241	HRP004163326	HRP004135468	HRP004222316
HRP004254984	HRP004270624	HRP004160697	HRP004278242	HRP004163327	HRP004135469	HRP004222321
HRP004254986	HRP004270625	HRP004160698	HRP004278243	HRP004163328	HRP004135470	HRP004222322
HRP004254987	HRP004270626	HRP004160699	HRP004278244	HRP004163329	HRP004135471	HRP004222323
HRP004254989	HRP004270627	HRP004231290	HRP004278245	HRP004163330	HRP004189134	HRP004222324
HRP004272850	HRP004270630	HRP004231291	HRP004278246	HRP004163331	HRP004189135	HRP004222325
HRP004264922	HRP004272353	HRP004231292	HRP004278248	HRP004163332	HRP004189136	HRP004222326
HRP004284245	HRP004177872	HRP004231293	HRP004278252	HRP004269858	HRP004189137	HRP004223334
HRP004284246	HRP004177873	HRP004231294	HRP004278352	HRP004269859	HRP004189138	HRP004223335
HRP004283784	HRP004258211	HRP004231295	HRP004278353	HRP004269860	HRP004189139	HRP004223336
HRP004283792	HRP004258821	HRP004231296	HRP004278355	HRP004269861	HRP004189140	HRP004223337
HRP004283793	HRP004258822	HRP004231693	HRP004278356	HRP004269862	HRP004189141	HRP004223338
HRP004283794	HRP004258823	HRP004231297	HRP004278357	HRP004269863	HRP004189142	HRP004223339
HRP004283795	HRP004256055	HRP004231298	HRP004278358	HRP004269864	HRP004189143	HRP004223340
HRP004283796	HRP004258212	HRP004231694	HRP004278916	HRP004269865	HRP004214053	HRP004223341
HRP004283797	HRP004258213	HRP004231695	HRP004282067	HRP004269866	HRP004214054	HRP004223342
HRP004283798	HRP004258214	HRP004232079	HRP004282078	HRP004269867	HRP004214055	HRP004223343
HRP004283799	HRP004258215	HRP004232080	HRP004118293	HRP004269889	HRP004214056	HRP004223344
HRP004283800	HRP004258216	HRP004232081	HRP004118294	HRP004269890	HRP004214057	HRP004223345
HRP004283801	HRP004258217	HRP004232082	HRP004118295	HRP004269891	HRP004214058	HRP004223346
HRP004283802	HRP004258218	HRP004232083	HRP004118296	HRP004269900	HRP004214059	HRP004223347
HRP004283803	HRP004258219	HRP004232084	HRP004118297	HRP004269902	HRP004214060	HRP004223348
HRP004283804	HRP004258220	HRP004232085	HRP004118298	HRP004216442	HRP004214061	HRP004223349



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004283805	HRP004258221	HRP004232086	HRP004281864	HRP004216443	HRP004214062	HRP004223350
HRP004265581	HRP004187089	HRP004231933	HRP004281865	HRP004216444	HRP004214063	HRP004223351
HRP004265582	HRP004187090	HRP004231934	HRP004281866	HRP004216445	HRP004116699	HRP004223352
HRP004265583	HRP004187091	HRP004231935	HRP004281867	HRP004216446	HRP004116705	HRP004223353
HRP004265584	HRP004187092	HRP004231936	HRP004280200	HRP004216447	HRP004116706	HRP004223967
HRP004265586	HRP004187093	HRP004231937	HRP004280201	HRP004216448	HRP004116712	HRP004223968
HRP004265587	HRP004187094	HRP004231938	HRP004282391	HRP004216449	HRP004116713	HRP004223969
HRP004265588	HRP004187095	HRP004231939	HRP004282899	HRP004216450	HRP004116986	HRP004223970
HRP004265589	HRP004187096	HRP004231940	HRP004284239	HRP004216451	HRP004116987	HRP004223971
HRP004265590	HRP004187097	HRP004231941	HRP004142624	HRP004227656	HRP004116988	HRP004223972
HRP004265591	HRP004187098	HRP004231942	HRP004142625	HRP004228087	HRP004116989	HRP004223973
HRP004265592	HRP004187099	HRP004156917	HRP004142626	HRP004228088	HRP004116990	HRP004267057
HRP004265593	HRP004187100	HRP004156918	HRP004142627	HRP004228089	HRP004239099	HRP004187084
HRP004265594	HRP004187101	HRP004156920	HRP004142628	HRP004228090	HRP004239102	HRP004240988
HRP004265595	HRP004187102	HRP004156921	HRP004230076	HRP004228091	HRP004239103	HRP004240989
HRP004265596	HRP004187103	HRP004156922	HRP004230077	HRP004227575	HRP004239104	HRP004280166
HRP004265597	HRP004187104	HRP004156923	HRP004230078	HRP004227576	HRP004239349	HRP004280167
HRP004265598	HRP004187105	HRP004156924	HRP004230079	HRP004227669	HRP004149273	HRP004146440
HRP004265599	HRP004187106	HRP004156925	HRP004230080	HRP004148983	HRP004149274	HRP004146441
HRP004265600	HRP004283971	HRP004156926	HRP004232188	HRP004148984	HRP004149275	HRP004146442
HRP004265601	HRP004283972	HRP004156927	HRP004232195	HRP004148985	HRP004149276	HRP004146443
HRP004265608	HRP004283973	HRP004247049	HRP004232196	HRP004149011	HRP004201655	HRP004146444
HRP004265609	HRP004283974	HRP004247050	HRP004232197	HRP004149012	HRP004201656	HRP004146445
HRP004265610	HRP004283975	HRP004247051	HRP004234555	HRP004149013	HRP004201678	HRP004146446
HRP004265611	HRP004283976	HRP004247052	HRP004244698	HRP004149014	HRP004201679	HRP004146447
HRP004154076	HRP004283977	HRP004247053	HRP004244699	HRP004149021	HRP004201680	HRP004146448
HRP004154089	HRP004283978	HRP004247054	HRP004244700	HRP004197315	HRP004201681	HRP004146449
HRP004154090	HRP004283979	HRP004247055	HRP004244701	HRP004197316	HRP004281898	HRP004146450
HRP004154091	HRP004283980	HRP004247056	HRP004244702	HRP004197319	HRP004281899	HRP004146451
HRP004154093	HRP004283981	HRP004247685	HRP004244703	HRP004197320	HRP004281900	HRP004146452
HRP004154094	HRP004283982	HRP004247686	HRP004281922	HRP004197321	HRP004281902	HRP004146453
HRP004154108	HRP004283983	HRP004145394	HRP004281933	HRP004197322	HRP004281906	HRP004146454
HRP004154109	HRP004283984	HRP004145395	HRP004281934	HRP004197323	HRP004281907	HRP004146455
HRP004154110	HRP004283985	HRP004145396	HRP004281935	HRP004197779	HRP004288181	HRP004146456
HRP004154111	HRP004283986	HRP004145501	HRP004281936	HRP004197780	HRP004288182	HRP004146457
HRP004154112	HRP004283987	HRP004145502	HRP004281937	HRP004197781	HRP004288203	HRP004146458
HRP004154113	HRP004283988	HRP004145503	HRP004281938	HRP004197782	HRP004288204	HRP004146459
HRP004154114	HRP004283989	HRP004287268	HRP004281939	HRP004197783	HRP004256056	HRP004146464
HRP004154115	HRP004283990	HRP004287274	HRP004281940	HRP004197784	HRP004217756	HRP004146465
HRP004154116	HRP004178802	HRP004144773	HRP004281941	HRP004197785	HRP004217757	HRP004146466
HRP004154124	HRP004263291	HRP004144774	HRP004160754	HRP004197786	HRP004217758	HRP004146467
HRP004154125	HRP004263292	HRP004145987	HRP004160755	HRP004197877	HRP004217759	HRP004146468
HRP004132658	HRP004263308	HRP004163574	HRP004160756	HRP004197879	HRP004217760	HRP004146469
HRP004132659	HRP004263309	HRP004163575	HRP004160757	HRP004197880	HRP004187812	HRP004146470
HRP004132660	HRP004263310	HRP004163576	HRP004160758	HRP004197901	HRP004188378	HRP004146471
HRP004132661	HRP004253690	HRP004163577	HRP004160759	HRP004197902	HRP004188379	HRP004146472
HRP004132662	HRP004253691	HRP004271429	HRP004160760	HRP004197903	HRP004188380	HRP004146473
HRP004145389	HRP004253692	HRP004271430	HRP004160761	HRP004197904	HRP004188381	HRP004272544
HRP004145390	HRP004253693	HRP004271431	HRP004160762	HRP004198658	HRP004188382	HRP004281628
HRP004145391	HRP004253694	HRP004271432	HRP004160763	HRP004198665	HRP004188383	HRP004193844
HRP004145392	HRP004253695	HRP004271433	HRP004160764	HRP004197317	HRP004188384	HRP004193845
HRP004145393	HRP004254315	HRP004271728	HRP004160765	HRP004197318	HRP004188577	HRP004193846
HRP004275319	HRP004254316	HRP004271729	HRP004160766	HRP004197584	HRP004188578	HRP004193847
HRP004275320	HRP004254317	HRP004271730	HRP004160767	HRP004197585	HRP004187434	HRP004193848



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004275354	HRP004254318	HRP004271739	HRP004160768	HRP004197586	HRP004187435	HRP004193849
HRP004275355	HRP004177881	HRP004271740	HRP004160769	HRP004197587	HRP004187436	HRP004193850
HRP004275358	HRP004177882	HRP004153315	HRP004160770	HRP004197588	HRP004187437	HRP004193857
HRP004275359	HRP004177883	HRP004153316	HRP004160771	HRP004197878	HRP004187603	HRP004193858
HRP004278354	HRP004177884	HRP004153317	HRP004160772	HRP004198236	HRP004187604	HRP004193859
HRP004278370	HRP004177885	HRP004152485	HRP004160773	HRP004198237	HRP004187605	HRP004193860
HRP004278371	HRP004177886	HRP004152486	HRP004202543	HRP004198238	HRP004187606	HRP004193861
HRP004278372	HRP004210789	HRP004152487	HRP004202544	HRP004198239	HRP004187808	HRP004193959
HRP004278373	HRP004210790	HRP004152488	HRP004202545	HRP004198240	HRP004187809	HRP004193960
HRP004278915	HRP004210791	HRP004152489	HRP004202546	HRP004198241	HRP004187810	HRP004193961
HRP004216452	HRP004210793	HRP004152490	HRP004202547	HRP004198242	HRP004187811	HRP004193962
HRP004216453	HRP004210794	HRP004152491	HRP004202548	HRP004198243	HRP004187985	HRP004193963
HRP004216454	HRP004210795	HRP004152492	HRP004202549	HRP004198399	HRP004187986	HRP004193964
HRP004216455	HRP004212393	HRP004153131	HRP004202550	HRP004198400	HRP004187987	HRP004193965
HRP004216456	HRP004212394	HRP004158275	HRP004202551	HRP004198401	HRP004187988	HRP004193966
HRP004216955	HRP004212395	HRP004158276	HRP004202552	HRP004198402	HRP004188210	HRP004193967
HRP004216956	HRP004212396	HRP004158277	HRP004203132	HRP004198403	HRP004188211	HRP004193968
HRP004216957	HRP004272561	HRP004158278	HRP004203133	HRP004198404	HRP004188212	HRP004193969
HRP004216958	HRP004272562	HRP004158279	HRP004203134	HRP004198405	HRP004188213	HRP004193970
HRP004216959	HRP004272563	HRP004158330	HRP004203135	HRP004198406	HRP004158708	HRP004193971
HRP004272571	HRP004272564	HRP004158331	HRP004203136	HRP004197583	HRP004158709	HRP004182453
HRP004272572	HRP004272565	HRP004158332	HRP004203137	HRP004197607	HRP004202682	HRP004182454
HRP004272574	HRP004272566	HRP004158333	HRP004203138	HRP004197608	HRP004214657	HRP004182455
HRP004272575	HRP004272567	HRP004158334	HRP004203139	HRP004198536	HRP004178561	HRP004200946
HRP004272596	HRP004272568	HRP004158335	HRP004203140	HRP004198537	HRP004178562	HRP004200947
HRP004138540	HRP004272569	HRP004158336	HRP004203141	HRP004198538	HRP004178563	HRP004200948
HRP004138541	HRP004272570	HRP004158337	HRP004203142	HRP004198539	HRP004178768	HRP004200949
HRP004138542	HRP004095490	HRP004158338	HRP004203143	HRP004198666	HRP004228098	HRP004200954
HRP004138543	HRP004095491	HRP004182426	HRP004203144	HRP004198667	HRP004228099	HRP004200955
HRP004138544	HRP004095492	HRP004182427	HRP004203145	HRP004198668	HRP004144946	HRP004200956
HRP004138546	HRP004095493	HRP004182428	HRP004203146	HRP004198669	HRP004144947	HRP004200957
HRP004138547	HRP004095494	HRP004177779	HRP004203147	HRP004198670	HRP004144948	HRP004200958
HRP004138548	HRP004095519	HRP004177780	HRP004203148	HRP004198671	HRP004144968	HRP004200959
HRP004138549	HRP004181058	HRP004177781	HRP004203149	HRP004198703	HRP004144969	HRP004201020
HRP004138550	HRP004181059	HRP004177782	HRP004203150	HRP004198704	HRP004144970	HRP004201021
HRP004164883	HRP004181065	HRP004177783	HRP004203151	HRP004198705	HRP004144971	HRP004201022
HRP004164884	HRP004181067	HRP004177808	HRP004163981	HRP004198706	HRP004144972	HRP004201023
HRP004164885	HRP004181068	HRP004177809	HRP004163982	HRP004198711	HRP004144973	HRP004201024
HRP004164886	HRP004181069	HRP004177810	HRP004207898	HRP004198712	HRP004144974	HRP004201025
HRP004164887	HRP004179391	HRP004177811	HRP004207903	HRP004198713	HRP004193994	HRP004201026
HRP004164895	HRP004179392	HRP004177812	HRP004207904	HRP004198714	HRP004193995	HRP004201027
HRP004164896	HRP004179393	HRP004177813	HRP004207905	HRP004198715	HRP004193996	HRP004201028
HRP004164897	HRP004179394	HRP004177814	HRP004207906	HRP002469397	HRP004193997	HRP004201029
HRP004164898	HRP004179395	HRP004177815	HRP004207907	HRP004210778	HRP004193998	HRP004272547
HRP004164981	HRP004179396	HRP004177816	HRP004207908	HRP004210779	HRP004105147	HRP004272552
HRP004164982	HRP004179399	HRP004178368	HRP004207909	HRP004211158	HRP004105157	HRP004272553
HRP004164983	HRP004179400	HRP004178369	HRP004207910	HRP004211159	HRP005105846	HRP004272554
HRP004164984	HRP004180013	HRP004178370	HRP004207911	HRP004211160	HRP004105845	HRP004272555
HRP004164985	HRP004180019	HRP004178371	HRP004207993	HRP004211161	HRP004105843	HRP004272556
HRP004164986	HRP004180020	HRP004178372	HRP004207994	HRP004211162	HRP004153634	HRP004272557
HRP004206974	HRP004180021	HRP004178373	HRP004207995	HRP004211163	HRP004153632	HRP004272558
HRP004206975	HRP004180634	HRP004177419	HRP004207996	HRP004211164	HRP004169160	HRP004272559
HRP004206976	HRP004180643	HRP004177420	HRP004207997	HRP004211165	HRP004169161	HRP004272560
HRP004206977	HRP004180644	HRP004177422	HRP004207998	HRP004177680	HRP004169162	HRP004278120



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004206982	HRP004180645	HRP004177423	HRP004270999	HRP004177681	HRP004169164	HRP004278121
HRP004206983	HRP004180946	HRP004177424	HRP004271000	HRP004177682	HRP004169165	HRP004244267
HRP004206984	HRP004180947	HRP004177956	HRP004271001	HRP004177683	HRP004169166	HRP004244268
HRP004206985	HRP004180948	HRP004177957	HRP004271002	HRP004177684	HRP004169167	HRP004244269
HRP004206986	HRP004180949	HRP004177958	HRP004271003	HRP004177685	HRP004169168	HRP004244270
HRP004206997	HRP004181260	HRP004177959	HRP004271004	HRP004177686	HRP004169169	HRP004244459
HRP004169530	HRP004181281	HRP004177960	HRP004271005	HRP004177687	HRP004169303	HRP004244460
HRP004169536	HRP004181282	HRP004177961	HRP004271006	HRP004177688	HRP004169304	HRP004244461
HRP004169571	HRP004181283	HRP004177962	HRP004271940	HRP004177689	HRP004169305	HRP004244462
HRP004169572	HRP004208496	HRP004177963	HRP004271941	HRP004114090	HRP004169306	HRP004164589
HRP004169573	HRP004208497	HRP004177964	HRP004271942	HRP004114091	HRP004169307	HRP004164590
HRP004169574	HRP004208498	HRP004177965	HRP004271943	HRP004114092	HRP004169308	HRP004164591
HRP004169575	HRP004208499	HRP004177966	HRP004271944	HRP004114093	HRP004169309	HRP004164592
HRP004169586	HRP004208500	HRP004177967	HRP004271945	HRP004114094	HRP004169310	HRP004164593
HRP004169587	HRP004208501	HRP004177968	HRP004285073	HRP004168237	HRP004169311	HRP004164594
HRP004169638	HRP004208502	HRP004177969	HRP004285074	HRP004168238	HRP004234227	HRP004164595
HRP004169639	HRP004208503	HRP004178374	HRP004285078	HRP004201687	HRP004234228	HRP004164732
HRP004169640	HRP004208847	HRP004178390	HRP004285099	HRP004201688	HRP004234229	HRP004164733
HRP004169670	HRP004208851	HRP004178804	HRP004285100	HRP004201690	HRP004234230	HRP004164734
HRP004169671	HRP004208862	HRP004178805	HRP004285101	HRP004201691	HRP004234231	HRP004164735
HRP004169672	HRP004208863	HRP004179116	HRP004285375	HRP004201839	HRP004234232	HRP004164736
HRP004169673	HRP004208864	HRP004179117	HRP004285376	HRP004201840	HRP004234233	HRP004164737
HRP004169674	HRP004208865	HRP004179118	HRP004285377	HRP004201841	HRP004234234	HRP004164738
HRP004169675	HRP004208866	HRP004179119	HRP004285378	HRP004201842	HRP004104686	HRP004164998
HRP004169676	HRP004208932	HRP004177414	HRP004237469	HRP004201843	HRP004104685	HRP004164999
HRP004169678	HRP004208933	HRP004177415	HRP004237490	HRP004201844	HRP004104684	HRP004165000
HRP004169679	HRP004208934	HRP004177416	HRP004237491	HRP004210782	HRP004104683	HRP004165001
HRP004169680	HRP004209156	HRP004177417	HRP004237492	HRP004210787	HRP004105148	HRP004165002
HRP004169681	HRP004209157	HRP004177418	HRP004237493	HRP004210788	HRP004105146	HRP004165003
HRP004169682	HRP004209227	HRP004275341	HRP004237494	HRP004163999	HRP004105844	HRP004165004
HRP004169683	HRP004209228	HRP004275342	HRP004237495	HRP004164000	HRP004105847	HRP004165171
HRP004169684	HRP004209229	HRP004275343	HRP004237496	HRP004164001	HRP004105842	HRP004165172
HRP004169685	HRP004209247	HRP004275344	HRP004237497	HRP004164002	HRP004105841	HRP004165173
HRP004169772	HRP004209248	HRP004275345	HRP004237498	HRP004164003	HRP004105840	HRP004165174
HRP004169773	HRP004209250	HRP004275346	HRP004237499	HRP004157212	HRP004178564	HRP004165175
HRP004168152	HRP004209449	HRP004275352	HRP004237500	HRP004157263	HRP004178773	HRP004165176
HRP004168153	HRP004209514	HRP004278917	HRP004237501	HRP004157264	HRP004178774	HRP004165177
HRP004169491	HRP004209525	HRP004278919	HRP004237532	HRP004157265	HRP004178775	HRP004165464
HRP004169492	HRP004209527	HRP004278920	HRP004237533	HRP004210415	HRP004178786	HRP004165485
HRP004169493	HRP004209532	HRP004278921	HRP004237534	HRP004210416	HRP004222645	HRP004165486
HRP004169494	HRP004209533	HRP004279151	HRP004237538	HRP004210417	HRP004222646	HRP004165487
HRP004169495	HRP004209552	HRP004279152	HRP004237539	HRP004210418	HRP004222647	HRP004165488
HRP004169520	HRP004209553	HRP004279153	HRP004237540	HRP004211243	HRP004222648	HRP004165489
HRP004169521	HRP004209554	HRP004279154	HRP004237605	HRP004211244	HRP004222649	HRP004165490
HRP004169529	HRP004209574	HRP004279155	HRP004237606	HRP004211245	HRP004222808	HRP004165613
HRP002612544.	HRP004209577	HRP004279156	HRP004237607	HRP004211247	HRP004222809	HRP004165614
HRP004132492	HRP004209578	HRP004279157	HRP004237608	HRP004211261	HRP004222810	HRP004165615
HRP004132493	HRP004209585	HRP004279158	HRP004237609	HRP004211262	HRP004222811	HRP004165616
HRP004132494	HRP004209587	HRP004272843	HRP004237610	HRP004189702	HRP004222812	HRP004224168
HRP004132495	HRP004209534	HRP004269907	HRP004237611	HRP004200301	HRP004222813	HRP004281987
HRP004132496	HRP004209540	HRP004269908	HRP004237612	HRP004200302	HRP004210423	HRP004281988
HRP004190478	HRP004209541	HRP004269909	HRP004237613	HRP004191271	HRP004210424	HRP004281990
HRP004190479	HRP004209575	HRP004288626	HRP004226908	HRP004191272	HRP004210425	HRP004281989
HRP004190480	HRP004209576	HRP004288627	HRP004182773	HRP004192359	HRP004210426	HRP002717123



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

FA-2022-016

HRP004190481	HRP004213096	HRP004288628	HRP004182928	HRP004192360	HRP004210427	HRP004171460
HRP004190482	HRP004213104	HRP004201685	HRP004182929	HRP004192436	HRP004210428	HRP004171461
HRP004190488	HRP004213105	HRP004201845	HRP004182930	HRP004192438	HRP004210780	HRP004171462
HRP004190489	HRP004213106	HRP004201846	HRP004182951	HRP004192439	HRP004210781	HRP004171463
HRP004190490	HRP004213107	HRP004201847	HRP004182952	HRP004192440	HRP004211685	HRP004171464
HRP004190491	HRP004213108	HRP004201848	HRP004182953	HRP004192441	HRP004211686	HRP004171465
HRP004190492	HRP004213109	HRP004201682	HRP004182954	HRP004192442	HRP004213474	HRP004171466
HRP004190493	HRP004213110	HRP004201683	HRP004182955	HRP004192443	HRP004213475	HRP004171467
HRP004190494	HRP004213111	HRP004201684	HRP004182960	HRP004192444	HRP004213476	HRP004171468
HRP004190495	HRP004213112	HRP004201686	HRP004182961	HRP004192445	HRP004213477	HRP004171469
HRP004190496	HRP004213440	HRP004201689	HRP004182962	HRP004192446	HRP004213478	HRP004171470
HRP004190640	HRP004213445	HRP004190978	HRP004182963	HRP004192447	HRP004138060	HRP004171471
HRP004190662	HRP004213446	HRP004190979	HRP004182964	HRP004192451	HRP004138061	HRP004171472
HRP004190663	HRP004213447	HRP004190980	HRP004183114	HRP004192452	HRP004138062	HRP004171473
HRP004190664	HRP004213468	HRP004190981	HRP004183115	HRP004192353	HRP004138063	HRP004171474
HRP004190665	HRP004213469	HRP004190982	HRP004183116	HRP004192355	HRP004138064	HRP004171475
HRP004190666	HRP004213470	HRP004190984	HRP004227204	HRP004192356	HRP004138065	HRP004171498
HRP004190667	HRP004213471	HRP004190985	HRP004227205	HRP004192357	HRP004138066	HRP004171499
HRP004190668	HRP004213472	HRP004190986	HRP004227206	HRP004192358	HRP004138067	HRP004171500
HRP004190669	HRP004213473	HRP004190987	HRP004227207	HRP004192374	HRP004138068	HRP004171501
HRP004190670	HRP004213711	HRP004190988	HRP004227212	HRP004192375	HRP004138069	HRP004171502
HRP004247691	HRP004213712	HRP004191090	HRP004227657	HRP004192376	HRP004138070	HRP004171504
HRP004231701	HRP004213713	HRP004191104	HRP004227658	HRP004192377	HRP004138071	HRP004171506
HRP004231702	HRP004213714	HRP004190977	HRP004227659	HRP004192378	HRP004138072	HRP004171507
HRP004231703	HRP004213715	HRP004207897	HRP004227660	HRP004192379	HRP004138073	HRP004171509
HRP004231704	HRP004213716	HRP004243058	HRP004227661	HRP004192380	HRP004138074	HRP004171511
HRP004231705	HRP004213766	HRP004240436	HRP004227662	HRP004192391	HRP004187081	HRP004200103
HRP004231706	HRP004213767	HRP004240437	HRP004227663	HRP004192437	HRP004187082	HRP004200109
HRP004232189	HRP004213768	HRP004240438	HRP004227664	HRP004192354	HRP004187083	HRP004183186
HRP004232228	HRP004213769	HRP004240439	HRP004227665	HRP004192450	HRP004284247	HRP004183187
HRP004232562	HRP004213770	HRP004263526	HRP004227666	HRP004192453	HRP004284248	HRP004183203
HRP004232563	HRP004213771	HRP004186719	HRP004227667	HRP004192458	HRP004284249	HRP004183204
HRP004174629	HRP004213772	HRP004186720	HRP004181764	HRP004192459	HRP004284250	HRP004183215
HRP004270607	HRP004213773	HRP004180379	HRP004277349	HRP004192460	HRP004202788	HRP004183216
HRP004270608	HRP004213774	HRP004180380	HRP004277350	HRP004192461	HRP004202789	HRP004183231
HRP004270609	HRP004208818	HRP004180381	HRP004239350	HRP004192462	HRP004203696	HRP004183232
HRP004270610	HRP004208819	HRP004180382	HRP004239922	HRP004192463	HRP004183510	HRP004191040
HRP004270611	HRP004208824	HRP004180383	HRP004253101	HRP004192464	HRP004183511	HRP004191041
HRP004270612	HRP004208825	HRP004180384	HRP004253102	HRP004192465	HRP004183512	HRP004191061
HRP004270613	HRP004209211	HRP004180389	HRP004253426	HRP004192466	HRP004183513	HRP004191062
HRP004270614	HRP004209212	HRP004180390	HRP004253427	HRP004192467	HRP004183545	HRP004191063
HRP004270615	HRP004209213	HRP004180396	HRP004253428	HRP004192468	HRP004183546	HRP004191064
HRP004270616	HRP004209218	HRP004180407	HRP004253429	HRP004192469	HRP004183549	HRP004191065
HRP004270617	HRP004209219	HRP004150046	HRP004253430	HRP004192470	HRP004183553	HRP004191066
HRP004164987	HRP004209220	HRP004150639	HRP004253431	HRP004192471	HRP004183935	HRP004191067
HRP004164988	HRP004209221	HRP004150641	HRP004253959	HRP004192472	HRP004183936	HRP004191068
HRP004164991	HRP004209226	HRP004150638	HRP004253960	HRP004194353	HRP004183938	HRP004191069
HRP004164992	HRP004209230	HRP004150640	HRP004253961	HRP004194354	HRP004183939	HRP004191070
HRP004164993	HRP004209249	HRP004150642	HRP004253962	HRP004194355	HRP004183940	HRP004135472
HRP004265602	HRP004209251	HRP004265057	HRP004253963	HRP004194356	HRP004183941	HRP004233261
HRP004265603	HRP004209440	HRP004265058	HRP004253964	HRP004194357	HRP004183942	HRP004273741
HRP004265604	HRP004209441	HRP004265060	HRP004254313	HRP004194358	HRP004183943	HRP004291417
HRP004265605	HRP004209442	HRP004265061	HRP004254314	HRP004194474	HRP004183944	HRP004186726
HRP004265606	HRP004209443	HRP004265062	HRP004138910	HRP004194475	HRP004183945	HRP004110611



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004265607	HRP004209444	HRP004265063	HRP004138911	HRP004194476	HRP004183946	HRP004110612
HRP004236157	HRP004209446	HRP004265064	HRP004138912	HRP004194497	HRP004183947	HRP004110613
HRP004236158	HRP004209447	HRP004264266	HRP004138913	HRP004194498	HRP004183948	HRP004110614
HRP004236159	HRP004209448	HRP004264267	HRP004138914	HRP004194499	HRP004183949	HRP004110615
HRP004236160	HRP004209512	HRP004264268	HRP004138915	HRP004111511	HRP004183950	HRP004110616
HRP004236161	HRP004209513	HRP004264269	HRP004138916	HRP004111512	HRP004183955	HRP004110617
HRP004236189	HRP004209528	HRP004264270	HRP004138917	HRP004111513	HRP004184373	HRP004110644
HRP004236190	HRP004209531	HRP004264272	HRP004138919	HRP004111514	HRP004183542	HRP004110645
HRP004236191	HRP004209535	HRP004265049	HRP004145364	HRP004111515	HRP004183543	HRP004110646
HRP004236192	HRP004209542	HRP004265050	HRP004145365	HRP004281632	HRP004183544	HRP004111066
HRP004236193	HRP004209543	HRP004265051	HRP004145366	HRP004281637	HRP004183547	HRP004111067
HRP004236194	HRP004209544	HRP004265052	HRP004145367	HRP004281640	HRP004183548	HRP004111068
HRP004236195	HRP004209545	HRP004265053	HRP004145369	HRP004281641	HRP004183550	HRP004111069
HRP004236196	HRP004209546	HRP004265054	HRP004280204	HRP004222329	HRP004183551	HRP004111070
HRP004236213	HRP004209547	HRP004265055	HRP004280205	HRP004222330	HRP004183552	HRP004111071
HRP004280237	HRP004209771	HRP004265056	HRP004280216	HRP004222331	HRP004183571	HRP004111072
HRP004280238	HRP004227577	HRP004265059	HRP004280217	HRP004222332	HRP004183937	HRP004111073
HRP004214652	HRP004227578	HRP004199266	HRP004280480	HRP004222333	HRP004184347	HRP004111074
HRP004214653	HRP004227579	HRP004199267	HRP004280481	HRP004222334	HRP004184348	HRP004111075
HRP004214654	HRP004227580	HRP004199268	HRP004280482	HRP004222335	HRP004184349	HRP004171848
HRP004214655	HRP004227581	HRP004199269	HRP004280483	HRP004222336	HRP004184370	HRP004171849
HRP004214656	HRP004227582	HRP004199270	HRP004280484	HRP004222327	HRP004184371	HRP004171850
HRP004241303	HRP004227583	HRP004199271	HRP004280485	HRP004222328	HRP004184372	HRP004171851
HRP004241304	HRP004227603	HRP004199272	HRP004280486	HRP004232861	HRP004184374	HRP004171852
HRP004241305	HRP004227604	HRP004199273	HRP004280487	HRP004232862	HRP004184375	HRP004171853
HRP004244873	HRP004227605	HRP004189556	HRP004280488	HRP004232863	HRP004184376	HRP004171854
HRP004244874	HRP004227670	HRP004189557	HRP004280489	HRP004232864	HRP004184377	HRP004171855
HRP004245123	HRP004227671	HRP004236150	HRP004280490	HRP004232865	HRP004184378	HRP004171856
HRP004245124	HRP004227672	HRP004236151	HRP004280491	HRP004157273	HRP004184379	HRP004172098
HRP004245125	HRP004227673	HRP004236152	HRP004280492	HRP004157274	HRP004184380	HRP004172099
HRP004245126	HRP004227674	HRP004132953	HRP004280493	HRP004157275	HRP004184381	HRP004172100
HRP004245127	HRP004227675	HRP004132954	HRP004280494	HRP004157276	HRP004184382	HRP004172101
HRP004245128	HRP004227141	HRP004132956	HRP004280495	HRP004157277	HRP004277355	HRP004172107
HRP004245129	HRP004227142	HRP004132957	HRP004280496	HRP004157280	HRP004277356	HRP004172108
HRP004245130	HRP004227143	HRP004132958	HRP004280497	HRP004157281	HRP004277357	HRP004172112
HRP004245131	HRP004227144	HRP004132959	HRP004280498	HRP004157287	HRP004277358	HRP004172113
HRP004245132	HRP004227588	HRP004132960	HRP004280499	HRP004157288	HRP004277359	HRP004172114
HRP004245315	HRP004227592	HRP004132961	HRP004280500	HRP004157289	HRP004277360	HRP004172262
HRP004245316	HRP004227593	HRP004132962	HRP004280501	HRP004202533	HRP004277361	HRP004172263
HRP004245317	HRP004227594	HRP004132963	HRP004280502	HRP004202534	HRP004277362	HRP004172264
HRP004245318	HRP004227595	HRP004133219	HRP004280503	HRP004202535	HRP004277363	HRP004172265
HRP004245319	HRP004227596	HRP004133220	HRP004282402	HRP004202536	HRP004277364	HRP004172266
HRP004245320	HRP004227597	HRP004133221	HRP004282403	HRP004202537	HRP004183332	HRP004172267
HRP004245321	HRP004227598	HRP004133222	HRP004282433	HRP004202538	HRP004183333	HRP004172268
HRP004245322	HRP004227606	HRP004133274	HRP004282434	HRP004202539	HRP004183382	HRP004172269
HRP004245323	HRP004227607	HRP004185555	HRP004282435	HRP004202540	HRP004183384	HRP004172559
HRP004286567	HRP004227608	HRP004185556	HRP004282436	HRP004202541	HRP004183657	HRP004172560
HRP004230081	HRP004227609	HRP004185557	HRP004282437	HRP004202542	HRP004183658	HRP004172561
HRP004230082	HRP004227610	HRP004185558	HRP004282438	HRP004133344	HRP004184110	HRP004172562
HRP004230083	HRP004241415	HRP004185559	HRP004282439	HRP004133345	HRP004184111	HRP004172563
HRP004230084	HRP004241416	HRP004185560	HRP004282450	HRP004133346	HRP004183294	HRP004172564
HRP004230085	HRP004241417	HRP004185561	HRP004282451	HRP004133347	HRP004183309	HRP004172565
HRP004230086	HRP004131809	HRP004185562	HRP004282456	HRP004173355	HRP004183334	HRP004172567
HRP004230087	HRP004131810	HRP004185563	HRP004282457	HRP004173356	HRP004183335	HRP004138678



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004230088	HRP004131811	HRP004185564	HRP004282458	HRP004173357	HRP004183412	HRP004138699
HRP004230089	HRP004131813	HRP004185565	HRP004282459	HRP004173358	HRP004184100	HRP004138702
HRP004181464	HRP004131820	HRP004185566	HRP004282460	HRP004173359	HRP004184101	HRP004138703
HRP004181465	HRP004131821	HRP004185567	HRP004282461	HRP004173360	HRP004184102	HRP004138704
HRP004181466	HRP004175001	HRP004185568	HRP004282462	HRP004173361	HRP004184103	HRP004138705
HRP004181467	HRP004175002	HRP004185569	HRP004282497	HRP004173400	HRP004184104	HRP004138706
HRP004181488	HRP004175003	HRP004185570	HRP004282498	HRP004173401	HRP004184105	HRP004138707
HRP004181892	HRP004220334	HRP004185571	HRP004283597	HRP004173402	HRP004184106	HRP004138708
HRP004181893	HRP004227266	HRP004185572	HRP004283598	HRP004173413	HRP004184107	HRP004138709
HRP004181894	HRP004227267	HRP004185573	HRP004283599	HRP004173414	HRP004184108	HRP004139025
HRP004181895	HRP004272545	HRP004185574	HRP004283600	HRP004226903	HRP004184109	HRP004139026
HRP004181896	HRP004272546	HRP004185575	HRP004283601	HRP004226904	HRP004237127	HRP004139027
HRP004182768	HRP004177971	HRP004214209	HRP004283602	HRP004226905	HRP004237129	HRP004139028
HRP004182769	HRP004177972	HRP004214218	HRP004283603	HRP004226906	HRP004237131	HRP004139029
HRP004182770	HRP004177973	HRP004214219	HRP004283604	HRP004226907	HRP004237132	HRP004139030
HRP004182771	HRP004177974	HRP004214220	HRP004283605	HRP004209169	HRP004237140	HRP004139031
HRP004182772	HRP004177975	HRP004214221	HRP004283606	HRP004266000	HRP004237141	HRP004139032
HRP004224171	HRP004207501	HRP004214222	HRP004283607	HRP004266001	HRP004237195	HRP004139033
HRP004224172	HRP004207502	HRP004214223	HRP004283608	HRP004266002	HRP004237197	HRP004139034
HRP004224173	HRP004207503	HRP004214224	HRP004283609	HRP004171439	HRP004237198	HRP004139328
HRP004224174	HRP004207530	HRP004214225	HRP004283610	HRP004171440	HRP004237199	HRP004139329
HRP004224175	HRP004207531	HRP004214226	HRP004283611	HRP004171441	HRP004237200	HRP004139330
HRP004224176	HRP004145374	HRP004214227	HRP004283612	HRP004171442	HRP004237202	HRP004139331
HRP004224571	HRP004145375	HRP004214228	HRP004283613	HRP004171443	HRP004188714	HRP004139332
HRP004224572	HRP004145376	HRP004214234	HRP004283614	HRP004171525	HRP004188715	HRP004139333
HRP004224573	HRP004145377	HRP004214235	HRP004283615	HRP004171562	HRP004188716	HRP004139334
HRP004224574	HRP004145378	HRP004214236	HRP004283616	HRP004173167	HRP004188717	HRP004139335
HRP004224575	HRP004145379	HRP004214237	HRP004214212	HRP004173168	HRP004188718	HRP004139336
HRP004224581	HRP004145380	HRP004214238	HRP004214217	HRP004173651	HRP004188719	HRP004139337
HRP004224612	HRP004145381	HRP004214239	HRP004171523	HRP004173652	HRP004188720	HRP004139537
HRP004224613	HRP004145382	HRP004214240	HRP004171524	HRP004173653	HRP004188721	HRP004139538
HRP004224614	HRP004145383	HRP004214241	HRP004280240	HRP004233288	HRP004188722	HRP004139539
HRP004224618	HRP004280218	HRP004214242	HRP004280241	HRP004233289	HRP004188723	HRP004139540
HRP004224619	HRP004280219	HRP004214210	HRP004280242	HRP004233290	HRP004188724	HRP004139541
HRP004224620	HRP004280220	HRP004214211	HRP004281629	HRP004233291	HRP004188725	HRP004139542
HRP004224621	HRP004280231	HRP004217491	HRP004281630	HRP004233292	HRP004188726	HRP004139543
HRP004224622	HRP004280232	HRP004217492	HRP004281631	HRP004157541	HRP004188727	HRP004139544
HRP004225698	HRP004169774	HRP004264273	HRP004283375	HRP004157542	HRP004188728	HRP004139545
HRP004225700	HRP004169775	HRP004264294	HRP004283391	HRP004157543	HRP004188729	HRP004139546
HRP004225701	HRP004203210	HRP004264295	HRP004283392	HRP004157564	HRP004188734	HRP004139787
HRP004225702	HRP004203211	HRP004264296	HRP004283393	HRP004157569	HRP004188735	HRP004139788
HRP004225703	HRP004145988	HRP004264297	HRP004283394	HRP004157570	HRP004188736	HRP004139789
HRP004225704	HRP004145989	HRP004264298	HRP004283395	HRP004157571	HRP004188738	HRP004139790
HRP004225705	HRP004157266	HRP004264299	HRP004283396	HRP004157572	HRP004188739	HRP004139791
HRP004225706	HRP004157272	HRP004264304	HRP004283397	HRP004157574	HRP004188740	HRP004139792
HRP004225707	HRP004236141	HRP004264305	HRP004283398	HRP004157575	HRP004188741	HRP004139813
HRP004225718	HRP004236142	HRP004264306	HRP004283399	HRP004158092	HRP004188742	HRP004139814
HRP004226513	HRP004236143	HRP004264307	HRP004283400	HRP004158093	HRP004188743	HRP004139815
HRP004226514	HRP004236144	HRP004264308	HRP004283401	HRP004158094	HRP004188744	HRP004139816
HRP004226515	HRP004236145	HRP004264309	HRP004283402	HRP004158095	HRP004188745	HRP004139817
HRP004226516	HRP004191419	HRP004264310	HRP004283403	HRP004158096	HRP004188746	HRP004139818
HRP004226517	HRP004191420	HRP004264311	HRP004283404	HRP004158097	HRP004188767	HRP004139819
HRP004226900	HRP004191421	HRP004264312	HRP004283405	HRP004158098	HRP004188768	HRP004139820
HRP004226901	HRP004191422	HRP004264313	HRP004283406	HRP004158099	HRP004188769	HRP004176486



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

FA-2022-016

HRP004226902	HRP004191423	HRP004264314	HRP004283736	HRP004158100	HRP004188770	HRP004268295
HRP004113936	HRP004209429	HRP004264315	HRP004283766	HRP004158101	HRP004235638	HRP004205579
HRP004113937	HRP004209434	HRP004264316	HRP004283767	HRP004161230	HRP004235639	HRP004132497
HRP004113938	HRP004209435	HRP004264317	HRP004283769	HRP004161556	HRP004235640	HRP004153130
HRP004113939	HRP004209583	HRP004264318	HRP004283770	HRP004161557	HRP004235641	HRP004135581
HRP004114081	HRP004209584	HRP004264319	HRP004283771	HRP004161558	HRP004235642	HRP004135582
HRP004114082	HRP004216950	HRP004264322	HRP004283772	HRP004162026	HRP004235643	HRP004135583
HRP004114083	HRP004216951	HRP004264323	HRP004274019	HRP004162027	HRP004211241	HRP004135584
HRP004114084	HRP004216952	HRP004264324	HRP004274020	HRP004162028	HRP004183692	HRP004135585
HRP004114085	HRP004216953	HRP004264325	HRP004274021	HRP004162029	HRP004183715	HRP004136157
HRP004114086	HRP004216954	HRP004264326	HRP004274022	HRP004162030	HRP004184047	HRP004136158
HRP004114087	HRP004249621	HRP004264327	HRP004274023	HRP004162031	HRP004137007	HRP004136159
HRP004114088	HRP004249622	HRP004264354	HRP004274024	HRP004162032	HRP004137301	HRP004136160
HRP004114089	HRP004249623	HRP004264355	HRP004274025	HRP004162033	HRP004137302	HRP004136161
HRP004242137	HRP004243057	HRP004264356	HRP004274154	HRP004162034	HRP004137341	HRP004136162
HRP004242378	HRP004243059	HRP004114235	HRP004274155	HRP004162035	HRP004137342	HRP004136723
HRP004242832	HRP004243060	HRP004114236	HRP004274156	HRP004162036	HRP004137343	HRP004136724
HRP004242833	HRP004244695	HRP004114237	HRP004274157	HRP004162037	HRP004137344	HRP004136725
HRP004242834	HRP004281661	HRP004114238	HRP004274158	HRP004162038	HRP004137345	HRP004136726
HRP004242835	HRP004281844	HRP004114239	HRP004274159	HRP004162039	HRP004137346	HRP004136727
HRP004242836	HRP004281845	HRP004114245	HRP004274160	HRP004162040	HRP004137347	HRP004286568
HRP004243124	HRP004281847	HRP004114246	HRP004274161	HRP004162041	HRP004137348	HRP004286569
HRP004242063	HRP004281848	HRP004114247	HRP004274569	HRP004263312	HRP004137349	HRP004113424
HRP004242065	HRP004281853	HRP004114248	HRP004274570	HRP004263313	HRP004137350	HRP004113425
HRP004242116	HRP004281854	HRP004114249	HRP004274571	HRP004263311	HRP004137351	HRP004113427
HRP004242151	HRP004281856	HRP004114250	HRP004274572	HRP004272447	HRP004137352	HRP004113428
HRP004242152	HRP004281857	HRP004114255	HRP004274573	HRP004272448	HRP004137353	HRP004113429
HRP004242153	HRP004281858	HRP004114256	HRP004281264	HRP004153640	HRP004183336	HRP004113430
HRP004242154	HRP004281849	HRP004114257	HRP004281265	HRP004153645	HRP004183343	HRP004113426
HRP004243123	HRP004281850	HRP004114258	HRP004281266	HRP004153646	HRP004183345	HRP004113431
HRP004207491	HRP004281851	HRP004114259	HRP004281633	HRP004209133	HRP004183383	HRP004113432
HRP004207492	HRP004281852	HRP004114261	HRP004281636	HRP004209134	HRP004183385	HRP004113433
HRP004207493	HRP004281855	HRP004114262	HRP004281639	HRP004209150	HRP004183661	HRP004113434
HRP004207494	HRP004158710	HRP004114263	HRP004281642	HRP004209151	HRP004183673	HRP004113435
HRP004207495	HRP004158711	HRP004114264	HRP004281643	HRP004209152	HRP004183675	HRP004113436
HRP004207496	HRP004158712	HRP004114267	HRP004281644	HRP004209155	HRP004183676	HRP004202790
HRP004207497	HRP004287263	HRP004114268	HRP004281645	HRP004209424	HRP004183687	HRP004164605
HRP004207498	HRP004287264	HRP004114269	HRP004282500	HRP004209425	HRP004184044	HRP004164615
HRP004207499	HRP004287265	HRP004114270	HRP004282532	HRP004209426	HRP004184046	HRP004164617
HRP004207500	HRP004287266	HRP004114271	HRP004283737	HRP004209427	HRP004115477	HRP004164641
HRP004249624	HRP004177326	HRP004114272	HRP004283738	HRP004209428	HRP004115478	HRP004164642
HRP004249625	HRP004177327	HRP004114273	HRP004283740	HRP004203643	HRP004115479	HRP004164643
HRP004249626	HRP004177328	HRP004114274	HRP004283741	HRP004203644	HRP004115480	HRP004164644
HRP004249627	HRP004177329	HRP004114275	HRP004283742	HRP004203569	HRP004115481	HRP004164645
HRP004194214	HRP004177330	HRP004114276	HRP004283744	HRP004203620	HRP004115482	HRP004164646
HRP004194215	HRP004177331	HRP004114277	HRP004283765	HRP004203621	HRP004115483	HRP004164649
HRP004194216	HRP004177332	HRP004114278	HRP004283773	HRP004203642	HRP004115484	HRP004164650
HRP004194217	HRP004177333	HRP004179050	HRP004240601	HRP004180667	HRP004115485	HRP004164651
HRP004194218	HRP004177334	HRP004179051	HRP004240602	HRP004180668	HRP004115486	HRP004164652
HRP004194219	HRP004177335	HRP004179052	HRP004240603	HRP004180669	HRP004127906	HRP004164653
HRP004194220	HRP004177336	HRP004179053	HRP004240604	HRP004180670	HRP004127907	HRP004164665
HRP004194221	HRP004177337	HRP004179082	HRP004240605	HRP004180674	HRP004127908	HRP004164666
HRP004194222	HRP004177338	HRP004179084	HRP004240606	HRP004180675	HRP004127909	HRP004164782
HRP004194223	HRP004177339	HRP004179085	HRP004240607	HRP004180676	HRP004128096	HRP004164783



Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice



FA-2022-016

HRP004194224	HRP004177340	HRP004179106	HRP004240608	HRP004180838	HRP004128097	HRP004164784
HRP004194225	HRP004177341	HRP004179107	HRP004240609	HRP004235862	HRP004128098	HRP004164785
HRP004194226	HRP004166075	HRP004179108	HRP004240610	HRP004235865	HRP004128099	HRP004164786
HRP004194227	HRP004166076	HRP004179114	HRP004240611	HRP004235866	HRP004128100	HRP004164787
HRP004194228	HRP004166077	HRP004179115	HRP004166196	HRP004227654	HRP004128101	HRP004164788
HRP004194229	HRP004166078	HRP004181295	HRP004166197	HRP004227655	HRP004128102	HRP004164789
HRP004194230	HRP004166079	HRP004181296	HRP004166198	HRP004228092	HRP004128103	HRP004164790
HRP004194231	HRP004194369	HRP004181297	HRP004166199	HRP004144776	HRP004128104	HRP004164791
HRP004194232	HRP004194370	HRP004181298	HRP004166200	HRP004144775	HRP004128105	HRP004164792
HRP004194233	HRP004167182	HRP004282258	HRP004166201	HRP004153881	HRP004128106	HRP004164793
HRP004194234	HRP004167183	HRP004282259	HRP004166202	HRP004153882	HRP004128107	HRP004164794
HRP004194235	HRP004271938	HRP004282260	HRP004166203	HRP004153898	HRP004128112	HRP004164795
HRP004218513	HRP004271939	HRP004282261	HRP004166698	HRP004153899	HRP004128113	HRP004164796
HRP004218514	HRP004289860	HRP004282262	HRP004166699	HRP004228072	HRP004128114	HRP004164797
HRP004218515	HRP004275353	HRP004282263	HRP004227128	HRP004228073	HRP004128115	HRP004193827
HRP004218516	HRP004275356	HRP004282264	HRP004227129	HRP004228074	HRP004128116	HRP004193828
HRP004218517	HRP004275357	HRP004282266	HRP004227130	HRP004228075	HRP004128117	HRP004193829
HRP004218518	HRP004278228	HRP004282267	HRP004227131	HRP004177970	HRP004128118	HRP004193830
HRP004218519	HRP004281244	HRP004282268	HRP004227136	HRP004216948	HRP004128119	HRP004193831
HRP004218520	HRP004270832	HRP004282269	HRP004227137	HRP004216949	HRP004128120	HRP004193842
HRP004218521	HRP004270833	HRP004282270	HRP004227138	HRP004113966	HRP004128232	HRP004193843
HRP004218522	HRP004270834	HRP004282271	HRP004227139	HRP004113967	HRP004128233	HRP004238603
HRP004218527	HRP004270835	HRP004282272	HRP004227140	HRP004113968	HRP004128284	HRP004238604
HRP004218528	HRP004270836	HRP004282273	HRP004227121	HRP004113969	HRP004128285	HRP004238605
HRP004218529	HRP004270837	HRP004282274	HRP004290408	HRP004113970	HRP004128317	HRP004238625
HRP004218530	HRP004153594	HRP004282275	HRP004290454	HRP004113971	HRP004128318	HRP004238631
HRP004218531	HRP004153598	HRP004282276	HRP004290457	HRP004113972	HRP004128319	HRP004238632
HRP004218532	HRP004153599	HRP004282281	HRP004290458	HRP004113973	HRP004128320	HRP004238633
HRP004218533	HRP004153600	HRP004282282	HRP004290459	HRP004114095	HRP004128321	HRP004238634
HRP004218635	HRP004153602	HRP004282283	HRP004290460	HRP004114096	HRP004128708	HRP004238635
HRP004218636	HRP004153603	HRP004282284	HRP004290461	HRP004207483	HRP004128709	HRP004238636
HRP004218637	HRP004153604	HRP004282285	HRP004290462	HRP004240579	HRP004128710	HRP004238638
HRP004189129	HRP004153605	HRP004282286	HRP004290463	HRP004240580	HRP004128711	HRP004238639
HRP004189130	HRP004153626	HRP004282287	HRP004290464	HRP004240581	HRP004128712	HRP004238652
HRP004189131	HRP004153627	HRP004282314	HRP004139821	HRP004240582	HRP004128713	HRP004238654
HRP004189132	HRP004171444	HRP004282318	HRP004226512	HRP004240583	HRP004128718	HRP004238655
HRP004189133	HRP004171445	HRP004282319	HRP004169066	HRP004240584	HRP004128719	HRP004238656
HRP004189349	HRP004171446	HRP004283147	HRP004169067	HRP004240585	HRP004128720	HRP004238657
HRP004189350	HRP004171447	HRP004283148	HRP004169068	HRP004240586	HRP004128725	HRP004237949
HRP004189351	HRP004171448	HRP004282277	HRP004192062	HRP004277288	HRP004128726	HRP004237950
HRP004189353	HRP004171449	HRP004282278	HRP004192075	HRP004277289	HRP004128727	HRP004237951
HRP004189354	HRP004171450	HRP004282279	HRP004192076	HRP004277290	HRP004128728	HRP004237952
HRP004189355	HRP004171451	HRP004282288	HRP004192077	HRP004277291	HRP004128729	HRP004237964
HRP004189531	HRP004171452	HRP004282289	HRP004192078	HRP004277292	HRP004128730	HRP004238015
HRP004189532	HRP004171453	HRP004282290	HRP004192079	HRP004281942	HRP004129102	HRP004238016
HRP004189533	HRP004115233	HRP004282291	HRP004192080	HRP004281943	HRP004127910	HRP004238416
HRP004189534	HRP004115234	HRP004282292	HRP004192081	HRP004281964	HRP004127911	HRP004238417
HRP004189535	HRP004155640	HRP004282293	HRP004192082	HRP004281973	HRP004127912	HRP004238418
HRP004189536	HRP004155641	HRP004282294	HRP004192083	HRP004281974	HRP004127913	HRP004238419
HRP004189566	HRP004155642	HRP004282295	HRP004192084	HRP004281975	HRP004127914	HRP004238420
HRP004189567	HRP004155643	HRP004282296	HRP004192085	HRP004281976	HRP004127915	HRP004238421
HRP004189568	HRP004155670	HRP004282297	HRP004194347	HRP004281977	HRP004127916	HRP004238422
HRP004189569	HRP004155671	HRP004282298	HRP004194348	HRP004281978	HRP004127921	HRP004238423
HRP004189570	HRP004155672	HRP004282299	HRP004194349	HRP004281979	HRP004127922	HRP004238424



Hillrom™

Hill-Rom S.A.S, ZI du Talhouet
56330 Pluvigner, France

Field Safety Notice

FA-2022-016

HRP004189571	HRP004155673	HRP004282300	HRP004194350	HRP004118086	HRP004127923	HRP004238425
HRP004189886	HRP004155674	HRP004282301	HRP004194352	HRP004118087	HRP004127924	HRP004238426
HRP004189887	HRP004155675	HRP004282302	HRP004163983	HRP004118088	HRP004127925	HRP004238427
HRP004189888	HRP004155681	HRP004282303	HRP004163984	HRP004118089	HRP004127926	HRP004238428
HRP004189889	HRP004155682	HRP004282304	HRP004163985	HRP004118090	HRP004127927	HRP004238429
HRP004189890	HRP004155683	HRP004282305	HRP004163986	HRP004118091	HRP004127928	HRP004238430
HRP004189891	HRP004155684	HRP004282306	HRP004163987	HRP004118092	HRP004127929	HRP004238431
HRP004240587	HRP004155685	HRP004282307	HRP004166080	HRP004118093	HRP004127935	HRP004238432
HRP004240588	HRP004155686	HRP004282308	HRP004166081	HRP004118094	HRP004127937	HRP004238433
HRP004240589	HRP004156038	HRP004282309	HRP004174999	HRP004118095	HRP004127943	HRP004238434
HRP004240590	HRP004156039	HRP004282310	HRP004175000	HRP004118096	HRP004127944	HRP004238435
HRP004240591	HRP004156040	HRP004282311	HRP004249357	HRP004118097	HRP004127945	HRP004238450
HRP004240592	HRP004156041	HRP004282312	HRP004249364	HRP004118098	HRP004127946	HRP004238572
HRP004240593	HRP004156042	HRP004282313	HRP004249385	HRP004118099	HRP004127947	HRP004238573
HRP004240594	HRP004156043	HRP004282315	HRP004249386	HRP004118100	HRP004127948	HRP004238574
HRP004240595	HRP004156044	HRP004282316	HRP004249387	HRP004118164	HRP004127953	HRP004238575
HRP004240596	HRP004156045	HRP004282317	HRP004249628	HRP004118165	HRP004128721	HRP004238593
HRP004244713	HRP004156046	HRP004283149	HRP004249629	HRP004118230	HRP004128722	HRP004238594
HRP004244714	HRP004231680	HRP004174420	HRP004249630	HRP004118265	HRP004128723	HRP004238595
HRP004244715	HRP004231681	HRP004174421	HRP004249631	HRP004118266	HRP004128724	HRP004238596
HRP004244716	HRP004231682	HRP004174422	HRP004249632	HRP004218638	HRP004128731	HRP004238597
HRP004244717	HRP004231683	HRP004174423	HRP004258222	HRP004228100	HRP004128732	HRP004238658
HRP004244905	HRP004231684	HRP004174424	HRP004280202	HRP004228101	HRP004128733	HRP004238659
HRP004244906	HRP004231685	HRP004173354	HRP004280203	HRP004228102	HRP004128734	HRP004238660
HRP004244907	HRP004231687	HRP004173403	HRP004281859	HRP004228103	HRP004128735	HRP004238661
HRP004244908	HRP004231688	HRP004173404	HRP004281860	HRP004228104	HRP004128736	HRP004238662
HRP004244909	HRP004231689	HRP004173412	HRP004281861	HRP004181022	HRP004128737	HRP004238663
HRP004251841	HRP004231690	HRP004173415	HRP004281862	HRP004181023	HRP004128741	HRP004238664
HRP004251864	HRP004231691	HRP004173416	HRP004281863	HRP004181024	HRP004128742	HRP004238665
HRP004251965	HRP004231692	HRP004173417	HRP004282900	HRP004181035	HRP004128743	HRP004238666
HRP004251966	HRP004273559	HRP004173418	HRP004284237	HRP004181036	HRP004128744	HRP004238667
HRP004251967	HRP004273560	HRP004173419	HRP004284238	HRP004181037	HRP004128745	HRP004238668
HRP004252305	HRP004273561	HRP004173420	HRP004182258	HRP004181038	HRP004128746	HRP004238669
HRP004252403	HRP004273562	HRP004173421	HRP004182259	HRP004181039	HRP004128747	HRP004238670
HRP004252404	HRP004273563	HRP004173422	HRP004182260	HRP004181040	HRP004128748	HRP004238671
HRP004252406	HRP004177454	HRP004173423	HRP004182261	HRP004181041	HRP004128749	HRP004238672
HRP004252407	HRP004177455	HRP004173424	HRP004182262	HRP004181042	HRP004128750	HRP004238673
HRP004252408	HRP004177456	HRP004173425	HRP004182263	HRP004181043	HRP004128751	HRP004238674
HRP004252409	HRP004177457	HRP004173426	HRP004182264	HRP004181044	HRP004128752	HRP004238675
HRP004252540	HRP004177850	HRP004173427				