

Summary of Product Characteristics

1 NAME OF THE MEDICINAL PRODUCT

innohep 10,000 IU/ml solution for injection

2 QUALITATIVE AND QUANTITATIVE COMPOSITION

Tinzaparin sodium 10,000 anti-Factor Xa IU/ml

Excipients with known effect:

Benzyl alcohol (10 mg/mL) and sodium (in total < 23 mg/mL).

For the full list of excipients, see section 6.1.

3 PHARMACEUTICAL FORM

Solution for injection.

Colourless or straw-coloured sterile aqueous solution.

4 CLINICAL PARTICULARS

4.1 Therapeutic Indications

Prophylaxis of venous thromboembolism in adult patients undergoing surgery, particularly orthopaedic, general or oncological surgery.

Prophylaxis of venous thromboembolism in non-surgical adult patients immobilised due to acute medical illness including: acute heart failure, acute respiratory failure, severe infections, active cancer, as well as exacerbation of rheumatic diseases.

Prevention of clotting in extracorporeal circuits during haemodialysis and haemofiltration in adults.

4.2 Posology and method of administration

Posology

Prophylaxis of thromboembolic events in adults:

Administration is by subcutaneous injection.

Surgical patients at moderate risk of thromboembolic events:

3,500 anti-Xa IU given SC 2 hours before surgery and then once daily for as long as the patient is considered to be at risk of VTE.

Surgical patients at high risk of thromboembolic events e.g. undergoing orthopaedic or cancer surgery:

4,500 anti-Xa IU given SC 12 hours before surgery and then once daily for as long as the patient is considered to be at risk of VTE.

Non-surgical patients immobilised due to acute medical illness:

3,500 anti-Xa IU given SC once daily in patients at moderate risk of VTE, or 4,500 anti-Xa IU given SC once daily in patients at high risk of VTE. Administration should continue for as long as the patient is considered to be at risk of VTE.

Neuraxial anaesthesia

Caution is advised when performing neuraxial anaesthesia or lumbar puncture in patients receiving prophylactic doses of innohep, see section 4.4: Neuraxial anaesthesia. If neuraxial anaesthesia is planned, a minimum delay of 12 hours should be allowed between the last prophylactic dose and the needle or catheter placement. innohep should not be resumed until at least 4-6 hours after the use of spinal anaesthesia or after the catheter has been removed. Thus, the 2 hours preoperative initiation of thromboprophylaxis with innohep is not compatible with neuraxial anaesthesia.

Haemodialysis and haemofiltration in adults:

Duration of 4 hours or less:

A bolus injection of 2,000 to 2,500 anti-Xa IU at the start of dialysis.

Duration of more than 4 hours:

A bolus injection of 2,500 anti-Xa IU at the start of dialysis/filtration, followed by 750 anti-Xa IU/hour as a continuous infusion.

Dose adjustment:

If necessary, the bolus dose may be increased or decreased gradually in increments of 500 anti-Xa IU until a satisfactory response is obtained. The usual dose is within 2,000 - 4,500 anti-Xa IU.

If case of concomitant transfusion of blood or concentrated red corpuscles, an extra bolus injection of 500 - 1,000 anti-Xa IU can be administered.

Dose monitoring:

Determination of plasma anti-Xa activity can be used to monitor the innohep dose during haemodialysis/haemofiltration. The plasma anti-Xa level should be approximately 0.5 anti-Xa IU/ml one hour after administration.

Interchangeability

For interchangeability with other LMWHs, see section 4.4.

Special populations*Paediatric population*

The safety and efficacy of innohep in children below 18 years have not yet been established. Currently available data are described in section 5.2, but no recommendation on a posology can be made.

Renal impairment

If renal impairment is suspected, renal function should be assessed using a formula based on serum creatinine to estimate creatinine clearance level.

Use in patients with a creatinine clearance level < 30 ml/minute is not recommended, as dosage in this population has not been established. Available evidence demonstrates no accumulation in patients with creatinine clearance levels down to 20 ml/min. When required in these patients, innohep administration can be initiated with anti-Xa monitoring, if the benefit outweighs the risk (see section 4.4: Renal impairment).

Elderly

innohep should be used in the elderly in standard doses. Precaution is recommended in the treatment of elderly patients with renal impairment. If renal impairment is suspected, see section 4.2: Renal impairment and section 4.4: Renal impairment.

Weight

For patients with very low or very high body weight, 50 anti-Xa IU per kg body weight once daily may be considered as an alternative to fixed dosing. For surgical patients, the first dose is given SC 2 hours before surgery. The administration should continue once daily for as long as the patient is considered to be at risk of VTE.

Method of administration

Parenteral products should be inspected visually prior to administration. Do not use if cloudiness or precipitate is observed. The liquid may turn yellow during storage but is still useable.

Administration is by subcutaneous injection when given as prophylaxis of thromboembolic events in adults. This can be done in abdominal skin, the outer side of the thigh, lower back, upper leg or upper arm. Do not inject in the area around the navel, near scars or in wounds.

For abdominal injections, the patient should be in a supine position, alternating the injections between the left and right side. The air-bubble within the syringe should not be removed. During the injection, the skin should be held in a fold.

For haemodialysis, the dose of innohep should be given into the arterial side of the dialyser or intravenously. The dialyser can be primed by flushing with 500–1,000 ml isotonic sodium chloride (9 mg/ml) containing 5,000 anti-Xa IU innohep per litre.

4.3 Contraindications

Hypersensitivity to the active substance or to any of the excipients listed in section 6.1.

Current or history of immune-mediated heparin-induced thrombocytopenia (type II) (see section 4.4).

Active major haemorrhage or conditions predisposing to major haemorrhage. Major haemorrhage is defined as fulfilling any one of these three criteria:

- a) occurs in a critical area or organ (e.g. intracranial, intraspinal, intraocular, retroperitoneal, intra-articular or pericardial, intra-uterine or intramuscular with compartment syndrome),
- b) causes a fall in haemoglobin level of 20 g/L (1.24 mmol/L) or more, or
- c) leads to transfusion of two or more units of whole blood or red cells.

Septic endocarditis

This multidose vial formulation of innohep contains 10 mg/ml of the preservative benzyl alcohol. Therefore, this formulation must not be given to premature babies and neonates due to the risk of gasping syndrome.

4.4 Special warnings and precautions for use

Neuraxial anaesthesia

In patients undergoing epidural/spinal anaesthesia or spinal puncture, the prophylactic use of heparins/low molecular weight heparins may be very rarely associated with epidural/spinal haematoma resulting in prolonged or permanent paralysis. This risk is increased by the use of: epidural/spinal catheter for anaesthesia; the concomitant use of drugs affecting haemostasis, such as non-steroidal anti-inflammatory drugs (NSAIDs), platelet inhibitors or anticoagulants; traumatic or repeated neuraxial puncture.

In decision-making on the interval between the last administration of heparin/low molecular weight heparin at prophylactic doses and the placement or removal of an epidural/spinal catheter, the pharmacokinetic profile of tinzaparin sodium (see section 5.2) and the patient profile should be taken into account.

A minimum delay of 12 hours should be allowed between the last prophylactic dose of innohep and the needle or catheter placement. For continuous techniques, a similar delay should be observed before removing the catheter. Moreover, innohep should not be resumed until at least 4-6 hours after the use of spinal anaesthesia or after the catheter has been removed.

Should a physician decide to administer anticoagulation in the context of epidural/spinal anaesthesia, extreme vigilance and frequent monitoring must be exercised to detect any signs and symptoms of neurologic impairment, such as back pain, sensory and motor deficits (numbness and weakness in lower limbs), bowel and/or bladder dysfunction. Nurses should be trained to detect such signs and symptoms. Patients should be instructed to immediately inform a nurse or a clinician if they experience any of these symptoms.

If signs or symptoms of epidural/spinal haematoma are suspected, urgent diagnosis and treatment including spinal cord decompression should be initiated.

Haemorrhage

Caution is advised when administering innohep to patients at risk of haemorrhage. For patients at risk of major haemorrhage, see section 4.3. The combination with medicinal products affecting platelet function or the coagulation system should be avoided or carefully monitored (see section 4.5).

Intramuscular injections

innohep must be administered subcutaneously and not via intramuscular injection. Due to the risk of haematoma, concomitant intramuscular injections should be avoided.

Heparin-induced thrombocytopenia

Platelet count should be measured before the start of treatment and periodically thereafter because of the risk of immune-mediated heparin-induced thrombocytopenia (type II). innohep must be discontinued in patients who develop immune-mediated heparin-induced thrombocytopenia (type II) (see sections 4.3 and 4.8). Platelet counts will usually normalise within 2 to 4 weeks after withdrawal.

Hyperkalaemia

Heparin products can suppress adrenal secretion of aldosterone, leading to hyperkalaemia. Risk factors include diabetes mellitus, chronic renal failure, pre-existing metabolic acidosis, raised plasma potassium at pre-treatment, concomitant therapy with drugs that may elevate plasma potassium, and long-term use of innohep. In patients at risk, potassium levels should be measured before starting innohep and monitored regularly thereafter. Heparin-related hyperkalaemia is usually reversible upon treatment discontinuation, though other approaches may need to be considered if innohep treatment is considered lifesaving (e.g. decreasing potassium intake, discontinuing other drugs that may affect potassium balance).

Prosthetic heart valves

There have been no adequate studies to assess the safe and effective use of innohep in preventing valve thrombosis in patients with prosthetic heart valves. The use of innohep cannot be recommended for this purpose.

Renal impairment

Use in patients with a creatinine clearance level < 30 ml/minute is not recommended, as dosage in this population has not been established. Available evidence demonstrates no accumulation in patients with creatinine clearance levels down to 20 ml/minute. When required in these patients, innohep administration can be used cautiously with anti-Xa monitoring, if the benefit outweighs the risk (see section 4.2). Although anti-Xa monitoring remains a poor predictor of haemorrhage risk, it is the most appropriate measure of the pharmacodynamic effects of innohep.

Elderly

Elderly are more likely to have reduced renal function (see section 4.4: Renal impairment); therefore caution should be exercised when prescribing innohep to the elderly.

Interchangeability

Low molecular weight heparins should not be used interchangeably because of differences in pharmacokinetics and biological activities. Switching to an alternative low molecular weight heparin, especially during extended use, must be exercised with particular caution and specific dosing instructions for each proprietary product must be followed.

Excipient warnings

This formulation of innohep contains 10mg/ml of the preservative, benzyl alcohol. Benzyl alcohol may cause allergic reactions. Benzyl alcohol may cause toxic and anaphylactoid reactions in infants and children up to 3 years old. High volumes should be used with caution and only if necessary, especially in subjects with liver or kidney impairment because of the risk of accumulation and toxicity (metabolic acidosis).

This medicinal product contains less than 1 mmol sodium (23 mg) per mL, i.e. essentially "sodium free".

4.5 Interaction with other medicinal products and other forms of interactions

The anticoagulant effect of innohep may be enhanced by other drugs affecting the coagulation system, such as those inhibiting platelet function (e.g. acetylsalicylic acid and other non-steroidal anti-inflammatory drugs), thrombolytic agents, vitamin K antagonists, activated protein C, small molecule anti-Xa and IIa inhibitors. Such combinations should be avoided or carefully monitored (see section 4.4).

4.6 Fertility, pregnancy and lactation

Pregnancy

Anticoagulant treatment of pregnant women requires specialist involvement.

Animal studies do not indicate direct or indirect harmful effects with respect to reproductive toxicity.

A large amount of data on pregnant women (more than 2,200 pregnancy outcomes) indicate no malformative or feto/neonatal toxicity of tinzaparin. Tinzaparin does not cross the placenta. innohep can be used during all trimesters of pregnancy, if clinically needed.

Epidural anaesthesia

Due to the risk of spinal haematoma, treatment doses of innohep are contraindicated in patients who receive neuraxial anaesthesia. Therefore, epidural anaesthesia in pregnant women should always be delayed until at least 24 hours after administration of the last treatment dose of innohep. Prophylactic doses may be used as long as a minimum delay of 12 hours is allowed between the last administration of innohep and the needle or catheter placement (see section 4.4).

Pregnant women with prosthetic heart valves

Therapeutic failures have been reported in pregnant women with prosthetic heart valves on full anti-coagulant doses of innohep and other low molecular weight heparins. innohep cannot be recommended for use in this population.

In the absence of clear dosing, efficacy and safety information in this circumstance, any attempt to anti-coagulate such patients must only be undertaken by medical practitioners with expertise and experience in this clinical area, and only if no safer alternative is available.

Excipients

innohep vials contain benzyl alcohol. As this preservative may cross the placenta and may cause accumulation and toxicity (metabolic acidosis), innohep formulations without benzyl alcohol (syringes) should be used during pregnancy.

Breast-feeding

Animal data indicate that innohep excretion into breast milk is minimal.

It is unknown whether tinzaparin is excreted into human milk. Although oral absorption of low molecular weight heparins is unlikely, a risk to the breast-fed child cannot be excluded.

In patients at risk the incidence of venous thromboembolism is particularly high during the first six weeks after child birth.

A decision must be made whether to discontinue breast-feeding or to discontinue/abstain from innohep therapy, taking into account the benefit of breast-feeding for the child and the benefit of therapy for the woman.

Excipients:

innohep vials contain benzyl alcohol. Due to a risk of accumulation and toxicity (metabolic acidosis), innohep vial formulations without benzyl alcohol (pre-filled syringes) are the preferred choice during breast-feeding.

Fertility

There are no clinical studies with innohep regarding fertility.

4.7 Effects on ability to drive and use machines

innohep has no or negligible influence on the ability to drive and use machines.

4.8 Undesirable effects

The most frequently reported undesirable effects are haemorrhage events, anaemia secondary to haemorrhage and injection site reactions.

Haemorrhage may present in any organ and have different degrees of severity. Complications may occur particularly when high doses are administered. Although major haemorrhages are uncommon, death or permanent disability has been reported in some cases.

Immune-mediated heparin-induced thrombocytopenia (type II) largely manifests within 5 to 14 days of receiving the first dose. Furthermore, a rapid-onset form has been described in patients previously exposed to heparin. innohep must be discontinued in all cases of immune-mediated heparin-induced thrombocytopenia (see section 4.4).

In rare cases, innohep may cause hypoaldosteronism associated with hyperkalaemia and metabolic acidosis. Patients at risk include those with diabetes mellitus or renal impairment (see section 4.4).

Serious allergic reactions may sometimes occur. These include rare cases of skin necrosis, toxic skin eruption (e.g. Stevens-Johnson syndrome), angioedema and anaphylaxis. Treatment should be promptly discontinued at the slightest suspicion of such severe reactions.

Priapism has been reported rarely.

The estimation of the frequency of undesirable effects is based on a pooled analysis of data from clinical studies and from spontaneous reporting.

Undesirable effects are listed by MedDRA SOC and the individual undesirable effects are listed starting with the most frequently reported. Within each frequency grouping, adverse reactions are presented in the order of decreasing seriousness.

Very common $\geq 1/10$

Common $\geq 1/100$ and $< 1/10$

Uncommon $\geq 1/1,000$ and $< 1/100$

Rare $\geq 1/10,000$ and $< 1/1,000$

Very rare $< 1/10,000$

Blood and lymphatic system disorders	
Common $\geq 1/100$ and $< 1/10$	Anaemia (incl. haemoglobin decreased)
Uncommon $\geq 1/1,000$ and $< 1/100$	Thrombocytopenia (type I) (incl. platelet count decreased)
Rare $\geq 1/10,000$ and $< 1/1,000$	Heparin induced thrombocytopenia (type II) Thrombocytosis
Immune system disorders	
Uncommon $\geq 1/1,000$ and $< 1/100$	Hypersensitivity
Rare $\geq 1/10,000$ and $< 1/1,000$	Anaphylactic reaction
Metabolism and nutrition disorders	
Rare $\geq 1/10,000$ and $< 1/1,000$	Hypoadosteronism associated with hyperkalaemia and metabolic acidosis
Vascular disorders	
Common $\geq 1/100$ and $< 1/10$	Haemorrhage Haematoma
Uncommon $\geq 1/1,000$ and $< 1/100$	Bruising, ecchymosis and purpura
Hepatobiliary disorders	
Uncommon $\geq 1/1,000$ and $< 1/100$	Hepatic enzyme increased (incl. increased transaminases, ALT, AST and GGT)
Skin and subcutaneous tissue disorders	
Uncommon $\geq 1/1,000$ and $< 1/100$	Dermatitis (incl. allergic dermatitis and bullous dermatitis) Rash Pruritus
Rare $\geq 1/10,000$ and $< 1/1,000$	Toxic skin eruption (including Stevens-Johnson syndrome) Skin necrosis Angioedema Urticaria
Musculoskeletal and connective tissue disorders	
Rare $\geq 1/10,000$ and $< 1/1,000$	Osteoporosis (in connection with long-term treatment)
Reproductive system and breast disorders	
Rare $\geq 1/10,000$ and $< 1/1,000$	Priapism
General disorders and administration site conditions	
Common $\geq 1/100$ and $< 1/10$	Injection site reaction (incl. injection site haematoma, haemorrhage, pain, pruritus, nodule, erythema and extravasation)

Paediatric population

Limited information derived from one study and postmarketing data indicates that the pattern of adverse reactions in children and adolescents is comparable to that in adults.

Reporting of suspected adverse reactions

Reporting suspected adverse reactions after authorisation of the medicinal product is important. It allows continued monitoring of the benefit/risk balance of the medicinal product. Healthcare professionals are asked to report any suspected adverse reactions via HPRC Pharmacovigilance, Earlsfort Terrace, IRL - Dublin 2; Tel: +353 1 6764971; Fax: +353 1 6762517. Website: www.hpra.ie; E-mail: medsafety@hpra.ie.

4.9 Overdose

Haemorrhage is the main complication of overdose. Due to the relatively short pharmacokinetic half-life of innohep (see section 5.2), minor haemorrhages can be managed conservatively following treatment discontinuation. Serious haemorrhage may require the administration of the antidote protamine sulfate. Patients should be carefully monitored.

5 PHARMACOLOGICAL PROPERTIES

5.1 Pharmacodynamic properties

Pharmacotherapeutic group: Antithrombotic Agents, ATC code: B01AB10

innohep is a low molecular weight heparin produced by controlled enzymatic depolymerisation of conventional sodium heparin derived from porcine mucosa.

innohep is an anticoagulant like conventional heparin. It potentiates the inhibition of several activated coagulation factors, especially Factor Xa, its activity being mediated via Antithrombin III.

The biological activity of innohep is standardised against the current international standard for low molecular weight heparins and is expressed in anti-factor Xa international units.

5.2 Pharmacokinetic properties

The absolute bioavailability based on anti-Xa activity after subcutaneous administration is approximately 90% and the time to reach maximal activity is 4-6 hours. The terminal elimination half-life is approximately 3.7 hours. Due to the long half-life of the pharmacological effect for innohep, once daily administration is sufficient.

The pharmacokinetic activities of tinzaparin have been studied in pregnancy. Data from sequential pharmacokinetic monitoring in 55 pregnancies suggests that pharmacokinetic properties do not differ from the non-pregnant state. There was a small, but non-statistically significant, decrease in anti-Xa levels with advancing gestation.

Paediatric population

Preliminary data on the use of tinzaparin suggest that younger children including neonates and infants clear tinzaparin faster and therefore might require higher doses than older children. However, data are not sufficient to allow for dosing recommendations, see section 4.2.

5.3 Preclinical safety data

No additional data.

6 PHARMACEUTICAL PARTICULARS

6.1 List of excipients

Sodium acetate trihydrate
Benzyl alcohol
Sodium hydroxide
Water for injections

6.2 Incompatibilities

In the absence of compatibility studies, this medicinal product must not be mixed with other medicinal products.

6.3 Shelf life

2 years.

Chemical and physical in use stability has been demonstrated for 28 days at 30°C.

From a microbiological point of view, once opened, the product may be stored for a maximum of 28 days at 30°C. Other in-use storage times and conditions are the responsibility of the user.

6.4 Special precautions for storage

This medicinal product does not require any special storage conditions.
For storage conditions after first opening of the medicinal product, see section 6.3.

6.5 Nature and contents of container

Colourless type I glass vial, sealed with a bromobutyl rubber stopper.
Each vial contains 2 ml of solution containing tinzaparin sodium 10,000 anti-Factor Xa IU/ml.
Pack sizes: 10 vial.

6.6 Special precautions for disposal and other handling

Do not use if cloudiness or particles are visible in the liquid. Any unused medicinal product or waste material should be disposed of in accordance with local requirements.

7 MARKETING AUTHORISATION HOLDER

Leo Laboratories Limited
285 Cashel Road
Dublin 12
Ireland

8 MARKETING AUTHORISATION NUMBER

PA0046/060/002

9 DATE OF FIRST AUTHORISATION/RENEWAL OF THE AUTHORISATION

Date of first authorisation: 21 October 1991

Date of last renewal: 14 February 2010

10 DATE OF REVISION OF THE TEXT

January 2020